



VERSION 24



Capabilities Catalog

INORGANIC CUSTOM & STOCK
CERTIFIED REFERENCE MATERIALS

ISO 17034 | ISO 17025 | ISO 9001

Our Employees Are to Our



Taylor Akers



Brian Alexander, PhD



Joseph Burns



Micah Cornett



Aaron Craggett



Christopher Cruz



Elizabeth Day



Madeline Denny



Justin Dirico



Christopher Gaines



Paul Gaines



Courtney Gibson



Madeline Gozzi



Uyen Gravley



Lee Hawthorne



Danielle Hinkley



Muzzammil Khan



James King



Thomas Kozikowski



Patrick Macdonald



William Marble



Ashley Michael



Amanda Miller



Shalin Presgraves



Courtney Rainer



Joshua Rancourt



Ethan Real



Mollie Reid



Abigail Resner



Angela Robson



Karen Sporakowski



Adam Stevens



Bethany Stewart



Kayla Stroupe



Kelsey Stroupe



Joshua Underwood



Zachary Underwood



Abigail Walker



Jodie Wall

the Key Element Success



Valerie Ballard



Michael Booth



Caitlin Dunavant



Marshall Durrett



Christopher Estes



Anna Falls



Jessica Farris



Liv Forbes



Brenda Francis



James Holcomb



Madeleine Hughes



Grace Hurst



Jeff Itle



Alyssa Johnson



Ashley Jones



James Jones



Erik Miller



Susan Newton



Lesley Owens, PhD



Autumn Phillips



Jesse Phillips



Ragan Phillips



Nick Plymale



Hunter Rogers



Steffy Sam



Donna Senn



Tammy Shepherd



Kevin Shrestha



Kshitish Shukla



Dixon Smiley



Kristin Stroupe



Christopher Sumner



Daniel Taylor



Katie Tindall



Daniel Todd



Rebecca Weddle



Andrew Wood



Emily-Colleen
Worthington-Deeble



Justin Yalung

- ★ Chief Executive Officer
- ★ Chief Operating Officer
- ▲ Chief Technical Officer
- ◆ Director

TABLE OF CONTENTS

WHY CHOOSE INORGANIC VENTURES?

■ Quality, Customs & More	5
■ Certificate of Analysis	6
■ Technical Support	7
■ Our Guarantee	7
■ Online Tech Center	8
■ IV Ignite	10
■ Annual ICP Conference.....	11
■ Transpiration Control Technology (TCT)	12

CUSTOM STANDARDS

■ Magic Happens with a Custom Solution	14
■ Benefits of Ordering a Custom	15
■ Customs Ordering Process	16
■ Total Volume Pricing	17

CROSS-REFERENCE TABLE

■ Cross-Reference Tables	18
--------------------------------	----

ICP & ICP-MS

■ Single-Element Standards	26
■ Cyanide Standards	40
■ Multi-Element Standards	41
■ High-Purity Ionization Buffers	55
■ USP Standards	56
■ Cannabis Standards	58

EPA STANDARDS

■ ILM03.0	62
■ ILM04.0	64
■ ILM05.2 & ILM05.3	66
■ Method 200.7	69
■ Method 200.8	79
■ Method 1311	82
■ Method 6020	83

ION CHROMATOGRAPHY

■ Anion Standards	89
■ Cation Standards	91
■ Multi-Ion Standards	92
■ Eluent Concentrates	93
■ EPA Standards	94
■ Instrument Cross-Reference Table – IONS	24

ATOMIC ABSORPTION

■ Single-Element Standards	98
■ Modifiers, Buffers & Releasing Agents	101

WATER QC

■ Potable Water Standards	103
■ Wastewater Standards	105

WET CHEMISTRY

Wet Chemical Standards

■ Conductivity Standards	108
■ pH Standards and Colored pH Standards	109
■ Cyanide Standards	110
■ pH Buffers Specially Formulated for USP <791>	111

Sample Preparation

■ Dissolution Reagents & Stabilizers	112
--	-----

Certified Titrants and Reagents

■ Certified Titrants	113
■ Reagents.....	114
■ Discrete Analyzer Reagents	114

INDEX

■ Index by Subject	115
■ Index by Catalog Number	116
■ Ordering, Terms & Conditions	121

Quality

A history of accreditation. For more than 20 years, Inorganic Ventures has been accredited by A2LA to ISO 17034 (formerly ISO Guide 34) & ISO 17025. These are the core standards of the analytical testing community, and Inorganic Ventures continues to lead the way in compliance to these quality standards. This means CRMs that are engineered to be stable, compatible, SI traceable and manufactured and tested under ISO 17034 & ISO 17025 guidelines.



Customs

Custom standards are Inorganic Ventures' specialty. Our catalog reveals only a fraction of the inorganic reference materials we can prepare. More than two thirds of our business is devoted entirely to preparing custom standards. As the leading manufacturer of custom inorganic standards, we've produced tens of thousands of unique blends for laboratories worldwide. It's our area of expertise, and perhaps the most prominent way in which we refine your results and redefine your industry.




And More...


On the web. Our technical library has been expanding for over a decade. Topics include ICP operations, sample preparation, trace metals analysis and much more*. There you'll discover the best online tool for analytical chemists with our Interactive Periodic Table. It includes chemical compatibilities, preferred lines, major interferences and additional data for 70+ elements. inorganicventures.com/tech-center

Additionally, our stock SDSs and CoAs can be found on our website for current lots as well as many older ones.

*For more new features, see pages 8-9.




WHY CHOOSE INORGANIC VENTURES? CERTIFICATE OF ANALYSIS



Compliant with ISO Guide 31:2015

You'll wonder how you ever got along without such a thorough certificate.
Contact us for a sample.



Certificate of Analysis (CoA)

Nearly every CRM we manufacture includes a highly detailed Certificate of Analysis. As an ISO 17034, A2LA accredited manufacturer, we provide certificates that include extensive data to meet the quality requirements of any laboratory:


- **Traceability** – to specific NIST SRMs and lots
- **Certified Values** – based on two independent methods
- **Trace Impurities** – listed with the actual values
- **Uncertainties** – detailed information reported

ONLINE

All CoA and Safety Data Sheet (SDS) information is now available online, 24/7. Inorganic Ventures is also pleased to share that all of our products are GHS compliant and our SDSs are available in 14 different languages.

inorganicventures.com/inorganic-standards

Inorganic Ventures Label



We're here to help. We don't just manufacture inorganic CRMs, we also provide technical support when it is needed so you can do your job. Because inorganic chemistry is all we do, Inorganic Ventures has a dedicated technical support team that can assist you with hundreds of topics: sample preparation, method development, ICP-OES and ICP-MS measurement issues and much more. You'll be amazed when you talk to a real person with a technical background ready to help you.

Our technical advisors are available to assist you Monday through Friday, 8:00 a.m. to 5:00 p.m. EST.

We can assist you with...

- Sample preparation
- Spectral interferences
- Chemical compatibilities
- Various ICP-OES & ICP-MS measurement issues



Technical Questions Answered

We've posted a variety of technical questions and answers pertaining to sample preparation, chemical stability and measurement.

inorganicventures.com/tech-center

Phone

- 800.669.6799 (US & Canada)
- +1.540.585.3030 (International)

Email

- info@inorganicventures.com

Online

- inorganicventures.com/forum

OUR GUARANTEE

Unquestionable integrity.

We believe in our products, and we value our customers. That is why every order leaving our facilities includes our "Declaration of Integrity." This document guarantees your satisfaction. Simply said, if you're dissatisfied with your order for any reason and we cannot work through the problem with you, a full refund will be issued, no questions asked.



WHY CHOOSE INORGANIC VENTURES? ONLINE TECH CENTER

Education

Chemistry Resources for Lab Techs and Spectroscopists

Inorganic Ventures' online technical library has been expanding for over a decade. Topics include ICP operations, sample preparation, trace metals analysis, and much more.



Interactive Periodic Table

Entries for each element, with details like storage and handling recommendations, chemical compatibility, and stability data intended specifically for use in spectroscopy.

[Learn More](#)



Technical Questions Forum

Whether you need detailed assistance or quick troubleshooting, the team at Inorganic Ventures is here to help with this searchable database. Can't find what you need? [Send us an email!](#)

[Learn More](#)



Technical Videos

Questions about our science or our products? Our PhDs break it down for you in a series of brief videos.

[Learn More](#)



ICP Operations Guide

A must-have guide for anyone operating and preparing samples and standards for measurement using ICP-MS and ICP-OES.

[Learn More](#)



Sample Preparation Guide

Our sample preparation guide provides specific, highly-detailed information about certain elements in regard to sample preparation.

[Learn More](#)



Trace Analysis Guide

An essential resource for trace analysts at any experience level, written by Paul Gaines, Ph.D.

[Learn More](#)



Periodic Table Guide

The Periodic Table for ICP Users Guide provides essential data for 70+ elements for every ICP user, including: chemical compatibilities, preferred emission lines, as well as major interferences and detection limits for both ICP and ICP-MS.

[Learn More](#)



Instrument Cross Reference Guide

Our cross-reference guide provides a quick key to Inorganic Ventures' part numbers to match the CRMs for your specific instrument.

[Learn More](#)



Guides and Papers

We believe that sharing knowledge lifts the scientific community as a whole. We hope you find this information helpful.

[Learn More](#)



TCT Information

Inorganic Ventures is committed to ensuring the accuracy of your CRMs, regardless of shipping or storage conditions. Learn more about the value of our industry-leading packaging technique.

[Learn More](#)



Bios

Learn more about our extraordinary technical team.

[Learn More](#)



Fun for Chemists

If you're looking for a laugh, check out our chemistry jokes and riddles!

[Learn More](#)



Links for Chemists

Professional organizations and associations, as well as EPA resources and other useful links.

[Learn More](#)



Frequently Asked Questions

Answers to questions we get asked a lot.

[Learn More](#)

[inorganicventures.com/tech-center](https://www.inorganicventures.com/tech-center)


Visit us online to see all of our upgraded features.



Bench Boost by Inorganic Ventures


Join our host and Inorganic Ventures Technical Director, Mike Booth, as he sparks insightful conversations with our panel of experts on Bench Boost, the go-to podcast for analytical chemistry enthusiasts. Together, they explore the realm of ICP, sharing tips and tricks and diving deep into all things analytical chemistry.






Transpiration Control Technology

With TCT, concerns about shipping or storage conditions are eliminated, as transpiration is no longer an issue.




Technical Videos

Watch technical videos pertaining to some of the most common questions in our recorded webinars and "Ask a Chemist" video series.




Interactive Periodic Table

Discover the best online tool for analytical chemists. Includes chemical compatibilities, preferred lines, major interferences and additional data for 70+ elements.



Guides and Papers

Inorganic Ventures' online technical library has been expanding for more than a decade. Topics include ICP operations, sample preparation, trace metals analysis and much more.



Experience
 Chemistry
 in a New

light



IV Ignite is a virtual training academy, covering topics related to
ICP-MS, and ICP-OES.

EXPERIENCE ALL THE COMMUNITY HAS TO OFFER

- Extensive Course Catalog
- *Bench Boost Podcast*: New episodes every Tuesday
- Downloadable Resources: interference tables, slides, etc.
- Forums to engage with fellow learners and instructors
- Earn PDUs and certificates as you complete milestones

COURSE CATALOG INCLUDES

- ICP Instrument Maintenance
- Controlling Sample Washout
- Common Sources of Contamination in the Lab

And so much more!



Try it FREE! To learn more or visit
pages.inorganicventures.com/iv-ignite



Inorganic Ventures' Annual ICP Conference

Calling all ICP users – don't miss our annual ICP Conference held in the fall.

You will hear from Inorganic Ventures' experts on a wide range of topics, including Sample Preparation Basics for ICP, Sample and CRM Stability Considerations, Trace Metals Analysis and much more.

TCT™

TRANSPIRATION CONTROL
TECHNOLOGY

The cornerstone of the scientific community is accuracy. That's why Inorganic Ventures has always been committed to producing the industry's most exact Certified Reference Materials.

But our control...and the control of every standard manufacturer...ends shortly after a standard is calibrated and packaged. We are improving the way we deliver our quality standards.



What is transpiration?

Transpiration refers to the passage of water vapor through the walls of a container and/or evaporation from the container opening. Transpiration results in an increase in the concentration of the CRM/RM.

What is the solution?

Transpiration Control Technology extends the shelf-life of the product. Inorganic Ventures uses a specially designed aluminized bag that prevents an increase in concentration of the CRM/RM until the TCT bag is opened.

How it works.

The sealed TCT bag stops the loss of water vapor from the bottle when equilibrium is reached inside the bag.

Has the product changed?

The product has not changed, it is the same high-quality product you have come to depend on from Inorganic Ventures. TCT is an investment we are making to extend shelf life and give you more control at no extra charge. Our products and unconditional guarantee remain the same.

What this means for you.

When you order standards from Inorganic Ventures, your product will be delivered in the TCT bag. This means you will be in control of the expiration date. Upon receiving the product, do not open the TCT bag until you are ready to use. To find out how long the product can be in the TCT bag before it expires, simply check the lot expiration found on the bottom left of the front label. Your product will expire on that date or one year after opening the sealed TCT bag, whichever comes first.

For more information on TCT, visit
inorganicventures.com/tct



Refine Your Results. Redefine Your Industry.

Inorganic Ventures' capabilities are not limited to a stock catalog. In fact, manufacturing custom standards is our passion and area of expertise. Let us lead the way as you refine your results and redefine your industry with our precise customizations.




Contents

Magic Happens with a Custom Solution	13
Benefits of Ordering a Custom.....	14
Primary Certified Reference Materials.....	15
Customs Ordering Process	16
Total Volume Pricing	17

Custom Standards

- ✓ Made to your exact specifications
- ✓ Save time and money
- ✓ Traceable to NIST SRMs
- ✓ Produced under ISO 9001
- ✓ Produced under ISO 17025
- ✓ Produced under ISO 17034
- ✓ 5-year shelf life with TCT



Mixing individual, single element standards into a working solution requires time, resources and money.

It involves preparation costs like time intensive labor and equipment, administrative costs associated with stocking and measuring re-order points for each solution, and the responsibility of handling all supporting documentation.



Our expert chemists make magic happen in the lab so that you can put these burdens by the wayside!

Each custom blend is calculated, formulated, articulated and regulated.

All of our custom blends are certified, NIST-traceable and have been put through a rigorous stability and method validation process.

Each custom standard is produced under ISO 9001, ISO 17025 and ISO 17034.



Our custom solutions are manufactured according to your exact specifications and ready for use with no preparation required!

In addition, each custom is guaranteed a 5-year shelf, thanks to our Transpiration Control Technology (TCT) which puts you in control of the expiration date!

Our team of expert chemists is available to facilitate sample preparation and troubleshoot any problems that may arise during testing.

BENEFITS OF ORDERING A CUSTOM

70,000 BLENDS AND COUNTING

Join thousands of laboratories worldwide in purchasing our custom standards.
Check out all the benefits of ordering a custom.

Potential issues building from single element stock products		BENEFITS OF ORDERING A CUSTOM	
High Preparation Costs		<p>Save on labor and equipment costs. Ready for immediate use with no prep required.</p> <p>Also save with custom pricing, calculated around your volume requirements.</p>	✓
Documentation Responsibilities		<p>All documentation and associated paperwork is handled for you and immediately available if you face an audit.</p>	✓
Uncertainty and Instability		<p>Certified, NIST traceable product from experts in elemental compatibility give you peace of mind.</p>	✓
Storage and Transpiration Issues		<p>Transpiration Control Technology (TCT) provides a 5-year shelf life and allows for flexible storage by increasing allowable temperature range.</p>	✓
Contamination		<p>Take advantage of our clean bottles and starting materials.</p>	✓
Limited Support Availability		<p>Work with dedicated chemists one-on-one to discuss your unique testing requirements.</p>	✓

How do I request a custom CRM?

Custom solutions can be requested through our convenient online quotation form. If you are unable to use our website, you may submit custom quotes through email, fax or by phone. All we need from you is a previously quoted IV part number, or for a brand new solution, the desired analytes, concentrations and matrix requirements.

www.inorganicventures.com/quote/instrumentsetup/index

email: quotes@inorganicventures.com

What happens after I submit my request?

Your custom quote is put together by our experienced chemists.

First, they go through our extensive library of more than 70,000 blends that have been developed over the last 38 years. If they don't find a match, they start formulating your exact custom standard. During this process, the blend is reviewed for stability and chemical compatibility. Your quote will be processed within two business days!*

I received the quote and I'm ready to order my custom!

If you like what you see, place your order via phone, fax or email. Your standard will be manufactured, packaged, and ready to ship within 10 business days. The custom standard is protected by our Transpiration Control Technology (TCT) and retains scientific integrity for five years from the date of manufacture.**



Custom Product Development Team

Representatives are available Monday through Friday, between 8:00 a.m. and 5:00 p.m. EST.

PHONE:
800.669.6799 / +1.540.585.3030

FAX:
540.585.3012

EMAIL:
quotes@inorganicventures.com

ONLINE:
inorganicventures.com


* Quotes and orders with many solutions, complex blends, or unique requirements may take longer to process.

** Based on stability data, some blends may receive a shorter lot expiration date.

Effortless Savings. Efficient Ordering.


With Total Volume Pricing on Custom Standards

1. Consolidate into ready-made.




Consolidate your single element standards into one ready-made custom standard. Compatibility concerns? Our world-class chemists will use your specifications and do the work for you.

2. Plan ahead with extended shelf life.



Calculate your annual usage of a custom product before you submit your request. TCT guarantees product shelf life for five years.

3. Enjoy zero-risk savings.



Single bottle price
\$\$\$ /bottle



Total Volume Pricing
\$ /bottle



Potential Savings
\$\$\$ /bottle

Order in volume to secure your savings and on-hand supply. We stand by our **100% guarantee** and will work with you to make sure your order suits your exact testing needs.

CROSS-REFERENCE TABLES

Inorganic Ventures is not affiliated with the companies and brands referenced on the following pages (other than Inorganic Ventures), and their names and marks are owned by the respective company and/or brand. The names appear solely for the purpose of permitting cross-referencing and comparison of products and standards.

AV Agilent/Varian	S Spectro
D DIONEX	T Thermo Scientific
H HACH	C Common Multi-Element Standards
JY HORIBA Jobin Yvon	I Common Multi-Ion Standards
M Merck/MilliporeSigma	U USP Method <232>
N NIST	
PE Perkin Elmer	

Common Multi-Element Standards C	
Inorganic Ventures#	Page
IV-STOCK-2	p.41
IV-STOCK-3	p.41
IV-STOCK-4	p.41
IV-STOCK-14	p.43
IV-STOCK-24	p.44
IV-STOCK-31	p.45
IV-STOCK-34	p.45
IV-STOCK-36	p.46
IV-STOCK-56	p.47
IV-STOCK-57	p.47
IV-STOCK-58	p.47
THM-TS-1	p.53

USP Method <232> U	
Inorganic Ventures#	Page
IV-STOCK-38	p.56
IV-STOCK-40	p.56
IV-STOCK-41	p.56
IV-STOCK-60	p.56
IV-STOCK-65	p.56
IV-STOCK-66	p.56
IV-STOCK-67	p.57
IV-STOCK-68	p.57
IV-STOCK-69	p.57
IV-STOCK-70	p.57
IV-STOCK-78	p.57
IV-STOCK-79	p.57

Products in bold are **near identical** formulations due to small differences in matrix percentages or additional elements.

Agilent/Varian AV		
Agilent/Varian#	Inorganic Ventures#	Page
5183-4681	IV-STOCK-53 IV-26191	p.46 <i>Custom</i>
5183-4682	IV-13676	<i>Custom</i>
5183-4687	IV-47679	<i>Custom</i>
5183-4688	IV-STOCK-50 IV-10934	p.46 <i>Custom</i>
5184-3566	IV-11304	<i>Custom</i>
5185-5959	IV-STOCK-74 IV-19645	p.47 <i>Custom</i>
5188-6524	IV-STOCK-51 IV-45975	p.46 <i>Custom</i>
5188-6525	IV-STOCK-75	p.47
5188-6526	IV-17685	<i>Custom</i>
5188-6527	6020ICS-9B IV-20338	p.85 <i>Custom</i>
5188-6564	AGI-TS-1 IV-43197	p.52 <i>Custom</i>
5190-0465	IV-37576	<i>Custom</i>
5190-7001	IV-ACID-BLANK	p.112
5190-8338	MM-MG-10	p.101
5190-8593	IV-63636	<i>Custom</i>
5190-8596	2008TS IV-20474	p.68, 81 <i>Custom</i>
5190-8599	CLPP-ICS-A IV-24398	p. 63, 65, 67 <i>Custom</i>
5190-9409	IV-STOCK-34 IV-28604	p.45 <i>Custom</i>
5190-9418	IV-12022	<i>Custom</i>
5190-9423	IV-36669	<i>Custom</i>

Agilent/Varian AV		
Agilent/Varian#	Inorganic Ventures#	Page
5190-9766	IV-STOCK-65	p.56, 60
5190-9770	IV-45454	<i>Custom</i>
5191-3933	IV-68161	<i>Custom</i>
5191-3934	IV-68162	<i>Custom</i>
5191-3935	IV-68163	<i>Custom</i>
5191-3936	IV-68164	<i>Custom</i>
5191-3937	IV-68165	<i>Custom</i>
5191-3938	IV-68166-A	<i>Custom</i>
5191-3939	IV-1739	<i>Custom</i>
8500-6940	IV-48812-A or IV-STOCK-27 & MSHGN-10PPM	<i>Custom</i> 44 27, 58
8500-6942	IV-STOCK-29	p.45
8500-6944	IV-STOCK-26	p.44
8500-6948	IV-STOCK-28	p.45
6610030000	IV-STOCK-24	p.44
6610030100	VAR-CAL-7 IV-8628	p.54 <i>Custom</i>
6610030400	VAR-IS-1 IV-31096	p.55 <i>Custom</i>
6610030500	VAR-CAL-1	p.54
6610030600	VAR-CAL-2 IV-45704	p.54 <i>Custom</i>
6610030700	IV-STOCK-33 IV-3105	p.45 <i>Custom</i>
G1820-60259	IV-DI-BLANK	p.112
ICM-240A	WW-IPC-1	p.74

To request a custom quote, contact quotes@inorganicventures.com, or visit www.inorganicventures.com/quote/instrumentsetup/index/.

CROSS-REFERENCE TABLES

DIONEX D		
DIONEX#	Inorganic Ventures#	Page
037157	IV-31102	<i>Custom</i>
037158	ICF1	p.89
037159	ICCL1	p.89
037160	ICSO41	p.90
037162	CARB	p.94
037163	BICARB	p.93
039513	ELUENT1817	p.93, 94
040187	IC-SCS1	p.92
046070	IV-8214	<i>Custom</i>
053560	ELUENT3510	p.93
056497	ICNO31	p.89
056933	IC-FAS-1A	p.92
056937	IV-39245	<i>Custom</i>
057590	IV-19190	<i>Custom</i>
063965	ELUENT4514	p.93
064161	IV-50921	<i>Custom</i>
088957	IV-77702	<i>Custom</i>
088964	IV-77704	<i>Custom</i>
088969	IV-77705	<i>Custom</i>
302511	IV-77703	<i>Custom</i>
302952	IV-50922	<i>Custom</i>
303167	ICCLO21	p.89, 96, 103
303168	ICBRO31	p.89, 96, 103
303169	ICNO21	p.89
303170	ICCLO31	p.89, 96, 103
303172	ICPO41	p.90
303173	IV-77701	<i>Custom</i>

HACH H		
HACH#	Inorganic Ventures#	Page
1222149	PH-10.01	p.110, 111
1222249	PH-7	p.109
1222256	PH-7	p.109
1222349	PH-4	p.109
1261742	MSPB-100PPM	p.28, 58
1279142	CGMN1	p.32
1279249	ICNNO31	p.89
12842	MSCU-100PPM	p.28
12932	MSCU-10PPM	p.26
13949	IV-23200	<i>Custom</i>
1402442	MSCD-100PPM	p.28, 58
14049	MSFE-10PPM	p.26
141453	PT-CO	p.114
1415142	IV-41139	<i>Custom</i>
1417442	MSAL-100PPM	p.28
1417542	MSFE-100PPM	p.28
1417642	CGNI1	p.32
1417742	CGZN1	p.35
1418642	CGMO1	p.32
1418742	MSMO-10PPM	p.27
1424232	IV-31277	<i>Custom</i>
1424342	IV-6164	<i>Custom</i>
1425010	IV-54485	<i>Custom</i>
1425249	IV-24939	<i>Custom</i>
1425310	IV-48775	<i>Custom</i>
1425610	IV-56147	<i>Custom</i>
1426010	IV-25195	<i>Custom</i>
1436716	IV-19009	<i>Custom</i>
1436832	IV-38443	<i>Custom</i>
1457142	CGAS1	p.30, 59
1461342	CGAG1	p.34
1466442	CGCR(6)1	p.31
1474949	CGNA1	p.34
15349	IV-14880	<i>Custom</i>
17149	IV-19189	<i>Custom</i>

CROSS-REFERENCE TABLES

HACH H		
HACH#	Inorganic Ventures#	Page
18349	CGICCL1	p.31
189149	IV-46777	Custom
191442	CGB1	p.30
194749	IV-31177	Custom
204649	IV-13101	Custom
2150342	CGC01	p.31
2175749	ICSO41	p.90
218710	CGCA10	p.36
2240349	CGCA1	p.30
227142	CGFE1	p.32
2283449	PHRED-4	p.109
2283456	PHRED-4	p.109
2283549	PHYELLOW-7	p.110
2283556	PHYELLOW-7	p.110
2283649	PHBLUE-10	p.110
2283656	PHBLUE-10	p.110
2305442	MSCA-10PPM	p.26
2305842	MSAL-10PPM	p.26
2318153	MSNA-100PPM	p.29
2321142	CGP1	p.33
23249	IV-31134	Custom
2340249	IV-23051	Custom
2351749	MSK-100PPM	p.29
2354153	IV-6827	Custom
2370853	IV-31116	Custom
2374820	MSPB-10PPM	p.26, 58
237842	MSZN-100PPM	p.29
2406549	IV-4035	Custom
2415132	IV-12426	Custom
256949	IV-18016	Custom
257849	IV-21404	Custom
259342	CGCU1	p.31
2605620	IV-29669	Custom
2605820	MSMN-10PPM	p.26
2714349	IV-60854	Custom

HACH H		
HACH#	Inorganic Ventures#	Page
2756549	IV-76390	Custom
2797149	IV-26620	Custom
2833049	IV-53879	Custom
2833449	IV-14891	Custom
2833349	IV-14890	Custom
2833649	IV-53878	Custom
2833749	IV-53877	Custom
2835153	MSNA-10PPM	p.27
29149	IV-18579	Custom
29153	IV-18579	Custom
2974549	IV-25155	Custom
2974649	IV-7532	Custom
2974749	IV-53788	Custom
30749	IV-12425	Custom
35949	IV-33253	Custom
40502	IV-20029	Custom
40505	IV-31132	Custom
40508	IV-18131	Custom
40515	IV-10137	Custom
40520	IV-30828	Custom
4457649	MSCA-100PPM	p.28
81042H	IV-30957	Custom
89149	IV-31293	Custom
C20C250	IV-45265	Custom
C20C270	CON1413-25	p.108
C20C280	CON147-25	p.108
S11M001	PH-1.68	p.109, 111
S11M002	PH-4	p.109
S11M003	PH-6.86	p.109, 111
S11M004	PH-7	p.109
S11M006-US	PH-9.18	p.110, 111
S11M008	PH-12.47	p.1010
S51M002	IV-49342	Custom
S51M003	IV-26474	Custom
S51M004	IV-2465	Custom

CROSS-REFERENCE TABLES

HORIBA Jobin Yvon JY		
Jobin Yvon#	Inorganic Ventures#	Page
100-2	PH-2	p.109
100-4	PH-4	p.109
100-7	PH-7	p.109
100-9	PH-9.18	p.110, 111
100-10	PH-10	p.110
500-2	PH-1.68	p.109, 111
500-4	PH-4	p.109
500-7	PH-7	p.109
500-9	PH-9.18	p.110, 111
500-10	PH-10.01	p.110, 111
500-21	CON84-25	p.108
500-22	CON1413-25	p.108
500-23	IV-45265	Custom
500-24	IV-61179	Custom
500-686	PH-6.86	p.109, 111
500-CA-SH	ICCA1	p.91
500-CA-SL	MSCA-100PPM	p.28
500-CL-SH	ICCL1	p.89
500-CL-SL	IV-31116	Custom
500-F-SH	ICF1	p.89
500-F-SL	IV-31134	Custom
500-K-SH	ICK1	p.91
500-K-SL	MSK-100PPM	p.29
500-NA-SH	ICNA1	p.91
500-NA-SL	MSNA-100PPM	p.29
500-NH4-SH	ICNH41	p.91
500-NH4-SL	IV-12420	Custom
500-NO3-SH	ICNO31	p.89
500-NO3-SL	IV-31186	Custom
JYICP-MIX7	IV-7	p.71, 77
JYICP-MIX7HSI	IV-21575	Custom
JYICP-MIX21	IV-21	p.72, 78
JYICP-MIX23	IV-STOCK-4 IV-19762	p.41 Custom

HORIBA Jobin Yvon JY		
Jobin Yvon#	Inorganic Ventures#	Page
JYICP-MIXHM	IV-STOCK-4 IV-56240	p.41 Custom
JYICP-MIXMAJ	IV-STOCK-34 IV-28604	p.45 Custom
JYICP-QC1	IV-43846	Custom
JYICP-QCACT	IV-73264	Custom

Merck/MilliporeSigma M		
Merck#	Inorganic Ventures#	Page
92091	IV-19459	Custom
104551	IV-8217	Custom
108704	ICPO41	p.90
109410	IV-STOCK-23	p.44
109411	IV-STOCK-24	p.44
109480	IV-8217	Custom
109481	IV-STOCK-14	p.43
109482	IV-STOCK-55	p.47
109487	IV-21	p.72, 78
109490	IV-19770	Custom
109491	IV-73446	Custom
109492	IV-STOCK-8 IV-63792	p.42 Custom
109493	IV-STOCK-10	p.42
109494	IV-72380	Custom
109495	IV-STOCK-17	p.43
109497	IV-16603	Custom
109498	IV-STOCK-21	p.43
109500	IV-STOCK-18	p.43
110322	IV-STOCK-7	p.42, 92
110580	IV-STOCK-6	p.41
110714	IV-STOCK-5	p.41
111355	IV-STOCK-4	p.41
115474	IV-STOCK-1	Custom
115626	IV-STOCK-16	p.43

Products in bold are **near identical** formulations due to small differences in matrix percentages or additional elements.

Perkin Elmer 		
Perkin Elmer#	Inorganic Ventures#	Page
N930021	IV-35379	<i>Custom</i>
B01906345	MM-MG-10	p.101
N0582152	IV-32705	<i>Custom</i>
N0681470	IV-STOCK-14	p.43
N0691579	IV-32706	<i>Custom</i>
N0691580	IV-37123	<i>Custom</i>
N8122017	IV-57482	<i>Custom</i>
N8122038	IV-48108	<i>Custom</i>
N8125030	IV-12922	<i>Custom</i>
N8125032	IV-STOCK-22	p.44
N8125033	IV-39880	<i>Custom</i>
N8125035	IV-9442	<i>Custom</i>
N8125040	IV-10890	<i>Custom</i>
N8125041	IV-9443	<i>Custom</i>
N8145050	IV-48108	<i>Custom</i>
N8145051	IV-STOCK-77	p.47
N8145052	IV-152	<i>Custom</i>
N8145053	IV-56364	<i>Custom</i>
N8145054	IV-13610	<i>Custom</i>
N8145055	IV-39880	<i>Custom</i>
N8145056	IV-41888	<i>Custom</i>
N8145059	IV-162	<i>Custom</i>
N8145060	IV-13373	<i>Custom</i>
N8145289	IV-66502	<i>Custom</i>
N9300200	IV-56238	<i>Custom</i>
N9300205	2007ICS-3	p.70
N9300208	IV-STOCK-54	p.46
N9300211	IV-56263	<i>Custom</i>
N9300213	IV-32109	<i>Custom</i>
N9300214	IV-19038	<i>Custom</i>
N9300215	IV-41880	<i>Custom</i>
N9300216	IV-56261-A	<i>Custom</i>
N9300217	IV-56262	<i>Custom</i>
N9300218	IV-STOCK-34	p.45

Perkin Elmer 		
Perkin Elmer#	Inorganic Ventures#	Page
N9300220	IV-35380	<i>Custom</i>
N9300221	IV-35377	<i>Custom</i>
N9300223	MSHGN-100PPM	p.29, 58
N9300224	IV-13160	<i>Custom</i>
N9300225	IV-9752	<i>Custom</i>
N9300226	CLPP-ICS-A	p.63, 65, 67
N9300227	CLPP-ICS-B4	p.65, 67
N9300230	IV-22905	<i>Custom</i>
N9300231	IV-STOCK-30	p.45
N9300232	IV-STOCK-26	p.44
N9300233	IV-48	<i>Custom</i>
N9300234	IV-STOCK-28	p.45
N9300235	IV-STOCK-29	p.45
N9300241	IV-46289	<i>Custom</i>
N9300244	IV-STOCK-18	p.43
N9300253	MSHGN-10PPM	p. 27, 58
N9300280	IV-21575	<i>Custom</i>
N9300281	IV-21	p. 72, 78
N9301720	IV-STOCK-21	p.43
N9301721	IV-14208	<i>Custom</i>
N9302946	IV-STOCK-55	p.47
N9303732	IV-48108	<i>Custom</i>
N9303813 N9303814	IV-DI-BLANK	p.114
N9303816	IV-1868	<i>Custom</i>
N9303818	IV-STOCK-35	p.45
N9303821	PE-CHK-1	p.52
N9303822	IV-15473	<i>Custom</i>
N9303824	IV-9636	<i>Custom</i>
N9303825	IV-75141-A	<i>Custom</i>
N9303827	6020ICS-9A	p.84
N9303828	6020ICS-0A	p.86
N9303831	IV-30429 & IV-56290-A	<i>Custom</i>

CROSS-REFERENCE TABLES

Products in bold are **near identical** formulations due to small differences in matrix percentages or additional elements.

Perkin Elmer PE		
Perkin Elmer#	Inorganic Ventures#	Page
N9303832 N9303833	IV-STOCK-53 IV-STOCK-53	p.46
N9303834	IV-13360	<i>Custom</i>
N9303835	IV-56512	<i>Custom</i>
N9303837	IV-21	p.72, 78
N9303839	IV-9633	<i>Custom</i>
N9303843	PE-TS-1	p.53
N9303940	IV-STOCK-1	<i>Custom</i>
N9303941	IV-STOCK-4	p.41
N9303942	IV-STOCK-8	p.42
N9303945	IV-STOCK-11	<i>Custom</i>
N9303946	IV-STOCK-13	p.42
N9303948	IV-56292	<i>Custom</i>
N9303949	WW-MSCAL-1	p.79
N9304131	IV-66049-A & IV-66049-B	<i>Custom</i>
N9307113	IV-25755	<i>Custom</i>
N9307114	IV-18652	<i>Custom</i>
N9307115	IV-33350	<i>Custom</i>
N9307116	IV-18653	<i>Custom</i>
N9307738	IV-56521	<i>Custom</i>
N9307741	IV-STOCK-6	p.41
N9307805	IV-STOCK-3	p.41
N9307806	IV-44364	<i>Custom</i>
N9307807	IV-56519	<i>Custom</i>
N9307808	IV-12490-A	<i>Custom</i>
N9307809 N9308571	IV-ACID-BLANK	p.112
N9308543	IV-4444	<i>Custom</i>

Spectro S		
Spectro#	Inorganic Ventures#	Page
USA00875	CIROS-OES-TS	p.52
USA00888	GENESIS-ICAL	p.52

NIST Multi-Element Standards N		
NIST#	Inorganic Ventures#	Page
SRM1643f	IV-STOCK-1643	p.52

Thermo Scientific T		
Thermo Scientific#	Inorganic Ventures#	Page
42885	IV-51075	<i>Custom</i>
12956213	IV-41701	<i>Custom</i>
1323760	THERMO-5A	p.53
1323770	THERMO-4AREV	p.53
ZG22950	TUNE F-X-SERIES	p.54
BRE0009578	IV-STOCK-80	p.47
4301 228 21401	IV-STOCK-31	p.45
4301 228 21411	IV-25579	<i>Custom</i>

IONS

Common Multi-Ion Standards I	
Inorganic Ventures#	Page
IC-FAS-1A	p.92
IC-SCS1	p.92
IV-STOCK-7	p.42, 92
IV-STOCK-59	p.92

To request a custom quote, contact quotes@inorganicventures.com, or visit www.inorganicventures.com/quote/instrumentsetup/index/.

Whether you use ICP-OES or ICP-MS, we offer a wide selection of Certified Reference Materials. At your request, we've expanded our line with new instrument setup standards. And we'll continue to improve our selection based on your feedback.



Contents

Single-Element Standards	19
Cyanide Standards	34
Instrument Cross-Reference Table.....	34
Multi-Element Standards	38
ICP-OES & ICP-MS	38
High-Purity Ionization Buffers.....	52
USP Standards	53
Cannabis Standards	55
Need a Custom CRM?	16

- ✓ **Five-year shelf life**
- ✓ **Traceable to NIST SRMs**
- ✓ **Produced under ISO 9001**
- ✓ **Produced under ISO 17025**
- ✓ **Produced under ISO 17034**
- ✓ **Assayed by validated wet chemical procedures**
- ✓ **Assayed by validated ICP-OES procedures**
- ✓ **Trace metallic impurities determined by ICP and ICP-MS**

SINGLE-ELEMENT STANDARDS

10 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

10 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Aluminum, Al	HNO ₃	125 mL	MSAL-10PPM-125ML
Antimony, Sb	HNO ₃ / Tartaric Acid	125 mL	MSSB-10PPM-125ML
Arsenic, As	HNO ₃	125 mL	MSAS-10PPM-125ML
Barium, Ba	HNO ₃	125 mL	MSBA-10PPM-125ML
Beryllium, Be	HNO ₃	125 mL	MSBE-10PPM-125ML
Bismuth, Bi	HNO ₃	125 mL	MSBI-10PPM-125ML
Boron, B	HNO ₃	125 mL	MSB-10PPM-125ML
¹⁰ Boron, ¹⁰ B	HNO ₃	100 mL*	MS10B-10PPM-100ML
¹¹ Boron, ¹¹ B	HNO ₃	100 mL*	MS11B-10PPM-100ML
Cadmium, Cd	HNO ₃	125 mL	MSCD-10PPM-125ML
Calcium, Ca	HNO ₃	125 mL 500 mL	MSCA-10PPM-125ML MSCA-10PPM-500ML
Cerium, Ce	HNO ₃	125 mL	MSCE-10PPM-125ML
Cesium, Cs	HNO ₃	125 mL	MSCS-10PPM-125ML
Chromium ⁺³ , Cr ⁺³	HNO ₃	125 mL	MSCR(3)-10PPM-125ML
Chromium ⁺⁶ , Cr ⁺⁶	H ₂ O	125 mL	MSCR(6)-10PPM-125ML
Cobalt, Co	HNO ₃	125 mL	MSCO-10PPM-125ML
Copper, Cu	HNO ₃	125 mL	MSCU-10PPM-125ML
Germanium, Ge	HNO ₃ / HF	125 mL	MSGE-10PPM-125ML
Gold, Au	HCl	125 mL 500 mL	MSAU-10PPM-125ML MSAU-10PPM-500ML
Hafnium, Hf	HNO ₃ / HF	125 mL 500 mL	MSHF-10PPM-125ML MSHF-10PPM-500ML
Holmium, Ho	HNO ₃	125 mL	MSHO-10PPM-125ML
Indium, In	HNO ₃	125 mL	MSIN-10PPM-125ML
Iron, Fe	HNO ₃	125 mL	MSFE-10PPM-125ML
Lead, Pb	HNO ₃	125 mL	MSPB-10PPM-125ML
Lithium, Li	HNO ₃	125 mL	MSLI-10PPM-125ML
⁶ Lithium, ⁶ Li	HNO ₃	125 mL	MS6LI-10PPM-125ML
Magnesium, Mg	HNO ₃	125 mL 500 mL	MSMG-10PPM-125ML MSMG-10PPM-500ML
Manganese, Mn	HNO ₃	125 mL 500 mL	MSMN-10PPM-125ML MSMN-10PPM-500ML

*Note: Size is 100 mL not 125 mL.

SINGLE-ELEMENT STANDARDS

10 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

10 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Mercury, Hg	HCl	125 mL	MSHG-10PPM-125ML
		500 mL	MSHG-10PPM-500ML
Mercury, Hg	HNO ₃	125 mL	MSHGN-10PPM-125ML
		500 mL	MSHGN-10PPM-500ML
Molybdenum, Mo	NH ₄ OH	125 mL	MSMO-10PPM-125ML
Nickel, Ni	HNO ₃	125 mL	MSNI-10PPM-125ML
Osmium, Os	HCl	125 mL	MSOS-10PPM-125ML
Phosphorus, P	H ₂ O	125 mL	MSP-10PPM-125ML
Platinum, Pt	HCl	125 mL	MSPT-10PPM-125ML
Potassium, K	HNO ₃	125 mL	MSK-10PPM-125ML
Rhodium, Rh	HCl	125 mL	MSRH-10PPM-125ML
Rhodium, Rh	HNO ₃	125 mL	MSRHN-10PPM-125ML
Scandium, Sc	HNO ₃	125 mL	MSSC-10PPM-125ML
Selenium, Se	HNO ₃	125 mL	MSSE-10PPM-125ML
Silicon, Si	HNO ₃ / HF	125 mL	MSSI-10PPM-125ML
Silver, Ag	HNO ₃	125 mL	MSAG-10PPM-125ML
Sodium, Na	HNO ₃	125 mL	MSNA-10PPM-125ML
		500 mL	MSNA-10PPM-500ML
Strontium, Sr	HNO ₃	125 mL	MSSR-10PPM-125ML
Sulfur, S	H ₂ O	125 mL	MSS-10PPM-125ML
Tellurium, Te	HNO ₃	125 mL	MSTEN-10PPM-125ML
Terbium, Tb	HNO ₃	125 mL	MSTB-10PPM-125ML
Thallium, Tl	HNO ₃	125 mL	MSTL-10PPM-125ML
Thorium, Th	HNO ₃	125 mL	MSTH-10PPM-125ML
Tin, Sn	HNO ₃ / HF	125 mL	MSSN-10PPM-125ML
		500 mL	MSSN-10PPM-500ML
Titanium, Ti	HNO ₃ / HF	125 mL	MSTI-10PPM-125ML
Tungsten, W	HNO ₃ / HF	125 mL	MSW-10PPM-125ML
Uranium, U	HNO ₃	125 mL	MSU-10PPM-125ML
		500 mL	MSU-10PPM-500ML
Vanadium, V	HNO ₃	125 mL	MSV-10PPM-125ML
Yttrium, Y	HNO ₃	125 mL	MSY-10PPM-125ML
Zinc, Zn	HNO ₃	125 mL	MSZN-10PPM-125ML
		500 mL	MSZN-10PPM-500ML

SINGLE-ELEMENT STANDARDS

100 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

100 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Aluminum, Al	HNO ₃	125 mL	MSAL-100PPM-125ML
Antimony, Sb	HNO ₃ / Tartaric Acid	125 mL	MSSB-100PPM-125ML
Arsenic, As	HNO ₃	125 mL	MSAS-100PPM-125ML
Barium, Ba	HNO ₃	125 mL	MSBA-100PPM-125ML
Beryllium, Be	HNO ₃	125 mL	MSBE-100PPM-125ML
Bismuth, Bi	HNO ₃	125 mL	MSBI-100PPM-125ML
Boron, B	HNO ₃	125 mL	MSB-100PPM-125ML
Cadmium, Cd	HNO ₃	125 mL	MSCD-100PPM-125ML
Calcium, Ca	HNO ₃	125 mL 500 mL	MSCA-100PPM-125ML MSCA-100PPM-500ML
Cerium, Ce	HNO ₃	125 mL	MSCE-100PPM-125ML
Cesium, Cs	HNO ₃	125 mL	MSCS-100PPM-125ML
Chromium ⁺³ , Cr ⁺³	HNO ₃	125 mL	MSCR(3)-100PPM-125ML
Chromium ⁺⁶ , Cr ⁺⁶	H ₂ O	125 mL	MSCR(6)-100PPM-125ML
Cobalt, Co	HNO ₃	125 mL	MSCO-100PPM-125ML
Copper, Cu	HNO ₃	125 mL	MSCU-100PPM-125ML
Germanium, Ge	HNO ₃ / HF	125 mL	MSGE-100PPM-125ML
Gold, Au	HCl	125 mL 500 mL	MSAU-100PPM-125ML MSAU-100PPM-500ML
Hafnium, Hf	HNO ₃ / HF	125 mL 500 mL	MSHF-100PPM-125ML MSHF-100PPM-500ML
Holmium, Ho	HNO ₃	125 mL	MSHO-100PPM-125ML
Indium, In	HNO ₃	125 mL	MSIN-100PPM-125ML
Iron, Fe	HNO ₃	125 mL 500 mL	MSFE-100PPM-125ML MSFE-100PPM-500ML
Lead, Pb	HNO ₃	125 mL 500 mL	MSPB-100PPM-125ML MSPB-100PPM-500ML
Lithium, Li	HNO ₃	125 mL 500 mL	MSLI-100PPM-125ML MSLI-100PPM-500ML
⁶ Lithium, ⁶ Li	HNO ₃	125 mL	MS6LI-100PPM-125ML
Magnesium, Mg	HNO ₃	125 mL	MSMG-100PPM-125ML
Manganese, Mn	HNO ₃	125 mL	MSMN-100PPM-125ML
Mercury, Hg	HCl	125 mL	MSHG-100PPM-125ML

SINGLE-ELEMENT STANDARDS

100 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

100 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Mercury, Hg	HNO ₃	125 mL	MSHGN-100PPM-125ML
Molybdenum, Mo	NH ₄ OH	125 mL	MSMO-100PPM-125ML
Nickel, Ni	HNO ₃	125 mL	MSNI-100PPM-125ML
Osmium, Os	HCl	125 mL	MSOS-100PPM-125ML
Phosphorus, P	H ₂ O	125 mL 500 mL	MSP-100PPM-125ML MSP-100PPM-500ML
Platinum, Pt	HCl	125 mL	MSPT-100PPM-125ML
Potassium, K	HNO ₃	125 mL 500 mL	MSK-100PPM-125ML MSK-100PPM-500ML
Rhodium, Rh	HCl	125 mL	MSRH-100PPM-125ML
Rhodium, Rh	HNO ₃	125 mL	MSRHN-100PPM-125ML
Scandium, Sc	HNO ₃	125 mL 500 mL	MSSC-100PPM-125ML MSSC-100PPM-500ML
Selenium, Se	HNO ₃	125 mL	MSSE-100PPM-125ML
Silicon, Si	HNO ₃ / HF	125 mL 500 mL	MSSI-100PPM-125ML MSSI-100PPM-500ML
Silver, Ag	HNO ₃	125 mL 500 mL	MSAG-100PPM-125ML MSAG-100PPM-500ML
Sodium, Na	HNO ₃	125 mL 500 mL	MSNA-100PPM-125ML MSNA-100PPM-500ML
Strontium, Sr	HNO ₃	125 mL	MSSR-100PPM-125ML
Sulfur, S	H ₂ O	125 mL	MSS-100PPM-125ML
Tellurium, Te	HNO ₃	125 mL	MSTEN-100PPM-125ML
Terbium, Tb	HNO ₃	125 mL	MSTB-100PPM-125ML
Thallium, Tl	HNO ₃	125 mL	MSTL-100PPM-125ML
Thorium, Th	HNO ₃	125 mL	MSTH-100PPM-125ML
Tin, Sn	HNO ₃ / HF	125 mL	MSSN-100PPM-125ML
Titanium, Ti	HNO ₃ / HF	125 mL	MSTI-100PPM-125ML
Tungsten, W	HNO ₃ / HF	125 mL	MSW-100PPM-125ML
Uranium, U	HNO ₃	125 mL 500 mL	MSU-100PPM-125ML MSU-100PPM-500ML
Vanadium, V	HNO ₃	125 mL	MSV-100PPM-125ML
Yttrium, Y	HNO ₃	125 mL 500 mL	MSY-100PPM-125ML MSY-100PPM-500ML
Zinc, Zn	HNO ₃	125 mL 500 mL	MSZN-100PPM-125ML MSZN-100PPM-500ML

SINGLE-ELEMENT STANDARDS

1,000 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Aluminum, Al	HNO ₃	30 mL	CGAL1-30ML
		125 mL	CGAL1-125ML
		500 mL	CGAL1-500ML
Aluminum, Al	HCl	30 mL	CGALCL1-30ML
		125 mL	CGALCL1-125ML
		500 mL	CGALCL1-500ML
Antimony, Sb	HNO ₃ / Tartaric Acid	125 mL	CGSB1-30ML
			CGSB1-125ML
			CGSB1-500ML
Antimony, Sb	HNO ₃ / HF	125 mL	CGSBF1-125ML
Arsenic, As	HNO ₃	30 mL	CGAS1-30ML
		125 mL	CGAS1-125ML
		500 mL	CGAS1-500ML
Arsenic ⁺³ , As ⁺³	HCl / NaOH / NaHCO ₃	30 mL	CGAS(3)1-30ML
		125 mL	CGAS(3)1-125ML
		500 mL	CGAS(3)1-500ML
Arsenic ⁺⁵ , As ⁺⁵	H ₂ O	30 mL	CGAS(5)1-30ML
		125 mL	CGAS(5)1-125ML
		500 mL	CGAS(5)1-500ML
Barium, Ba	HNO ₃	30 mL	CGBA1-30ML
		125 mL	CGBA1-125ML
		500 mL	CGBA1-500ML
Beryllium, Be	HNO ₃	30 mL	CGBE1-30ML
		125 mL	CGBE1-125ML
		500 mL	CGBE1-500ML
Bismuth, Bi <small>Commonly used as an Internal Standard for ICP-MS.</small>	HNO ₃	30 mL	CGBI1-30ML
		125 mL	CGBI1-125ML
		500 mL	CGBI1-500ML
Boron, B	NH ₄ OH	30 mL	CGB1-30ML
		125 mL	CGB1-125ML
		500 mL	CGB1-500ML
Bromide, Br- <small>Suitable for analyzing Bromide by ICP-OES.</small>	H ₂ O	30 mL	CGICBR1-30ML
		125 mL	CGICBR1-125ML
		500 mL	CGICBR1-500ML
Cadmium, Cd	HNO ₃	30 mL	CGCD1-30ML
		125 mL	CGCD1-125ML
		500 mL	CGCD1-500ML
Calcium, Ca	HNO ₃	30 mL	CGCA1-30ML
		125 mL	CGCA1-125ML
		500 mL	CGCA1-500ML
Carbon, C	HNO ₃	125 mL	CGC1-125ML
		500 mL	CGC1-500ML
Carbon, C <small>Suitable for TOC applications per Standard Methods.</small>	H ₂ O	125 mL	TOCKHP1-125ML
		500 mL	TOCKHP1-500ML
Cerium, Ce	HNO ₃	30 mL	CGCE1-30ML
		125 mL	CGCE1-125ML
		500 mL	CGCE1-500ML

SINGLE-ELEMENT STANDARDS

1,000 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Cesium, Cs	HNO ₃	30 mL	CGCS1-30ML
		125 mL	CGCS1-125ML
Chloride, Cl⁻ <small>Suitable for analyzing Chloride by ICP-OES.</small>	H ₂ O	125 mL	CGICCL1-125ML
		500 mL	CGICCL1-500ML
Chromium⁺³, Cr⁺³	HNO ₃	30 mL	CGCR(3)1-30ML
		125 mL	CGCR(3)1-125ML
		500 mL	CGCR(3)1-500ML
Chromium⁺⁶, Cr⁺⁶	H ₂ O	30 mL	CGCR(6)1-30ML
		125 mL	CGCR(6)1-125ML
		500 mL	CGCR(6)1-500ML
Cobalt, Co	HNO ₃	30 mL	CGCO1-30ML
		125 mL	CGCO1-125ML
		500 mL	CGCO1-500ML
Copper, Cu	HNO ₃	30 mL	CGCU1-30ML
		125 mL	CGCU1-125ML
		500 mL	CGCU1-500ML
Dysprosium, Dy	HNO ₃	30 mL	CGDY1-30ML
		125 mL	CGDY1-125ML
		500 mL	CGDY1-500ML
Erbium, Er	HNO ₃	30 mL	CGER1-30ML
		125 mL	CGER1-125ML
		500 mL	CGER1-500ML
Europium, Eu	HNO ₃	30 mL	CGEU1-30ML
		125 mL	CGEU1-125ML
		500 mL	CGEU1-500ML
Gadolinium, Gd	HNO ₃	30 mL	CGGD1-30ML
		125 mL	CGGD1-125ML
		500 mL	CGGD1-500ML
Gallium, Ga	HNO ₃	30 mL	CGGA1-30ML
		125 mL	CGGA1-125ML
		500 mL	CGGA1-500ML
Germanium, Ge	HNO ₃ / HF	30 mL	CGGE1-30ML
		125 mL	CGGE1-125ML
		500 mL	CGGE1-500ML
Gold, Au <small>Can also be used to stabilize low-level Hg for ICP-MS analysis.</small>	HCl	30 mL	CGAU1-30ML
		125 mL	CGAU1-125ML
		500 mL	CGAU1-500ML
Gold, Au	HNO ₃	30 mL	CGAUN1-30ML
		125 mL	CGAUN1-125ML
		500 mL	CGAUN1-500ML
Hafnium, Hf	HNO ₃ / HF	30 mL	CGHF1-30ML
		125 mL	CGHF1-125ML
		500 mL	CGHF1-500ML
Holmium, Ho <small>Commonly used as an Internal Standard for ICP-MS.</small>	HNO ₃	30 mL	CGHO1-30ML
		125 mL	CGHO1-125ML
		500 mL	CGHO1-500ML

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Indium, In Commonly used as an Internal Standard for ICP-MS.	HNO ₃	30 mL	CGIN1-30ML
		125 mL	CGIN1-125ML
		500 mL	CGIN1-500ML
Iodide, I- Suitable for analyzing Iodide by ICP-OES.	H ₂ O / TEA	30 mL	CGICI1-30ML
		125 mL	CGICI1-125ML
		500 mL	CGICI1-500ML
Iridium, Ir	HCl	30 mL	CGIR1-30ML
		125 mL	CGIR1-125ML
		500 mL	CGIR1-500ML
Iron, Fe	HNO ₃	30 mL	CGFE1-30ML
		125 mL	CGFE1-125ML
		500 mL	CGFE1-500ML
Lanthanum, La	HNO ₃	30 mL	CGLA1-30ML
		125 mL	CGLA1-125ML
		500 mL	CGLA1-500ML
Lead, Pb	HNO ₃	30 mL	CGPB1-30ML
		125 mL	CGPB1-125ML
		500 mL	CGPB1-500ML
Lithium, Li	HNO ₃	30 mL	CGLI1-30ML
		125 mL	CGLI1-125ML
		500 mL	CGLI1-500ML
⁶Lithium, ⁶Li Commonly used as an Internal Standard for ICP-MS.	HNO ₃	30 mL	CG6LI1-30ML
		125 mL	CG6LI1-125ML
Lutetium, Lu	HNO ₃	30 mL	CGLU1-30ML
		125 mL	CGLU1-125ML
		500 mL	CGLU1-500ML
Magnesium, Mg	HNO ₃	30 mL	CGMG1-30ML
		125 mL	CGMG1-125ML
		500 mL	CGMG1-500ML
Manganese, Mn	HNO ₃	30 mL	CGMN1-30ML
		125 mL	CGMN1-125ML
		500 mL	CGMN1-500ML
Mercury, Hg	HNO ₃	30 mL	CGHG1-30ML
		125 mL	CGHG1-125ML
		500 mL	CGHG1-500ML
Molybdenum, Mo	NH ₄ OH	30 mL	CGMO1-30ML
		125 mL	CGMO1-125ML
		500 mL	CGMO1-500ML
Neodymium, Nd	HNO ₃	30 mL	CGND1-30ML
		125 mL	CGND1-125ML
		500 mL	CGND1-500ML
Nickel, Ni	HNO ₃	30 mL	CGNI1-30ML
		125 mL	CGNI1-125ML
		500 mL	CGNI1-500ML
Niobium, Nb	HNO ₃ / HF	30 mL	CGNB1-30ML
		125 mL	CGNB1-125ML
		500 mL	CGNB1-500ML

SINGLE-ELEMENT STANDARDS

1,000 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Niobium, Nb High purity, low Tantalum	HNO ₃ / HF	125 mL	CGNB2051-125ML
		500 mL	CGNB2051-500ML
Osmium, Os	HCl	30 mL	CGOS1-30ML
		125 mL	CGOS1-125ML
		500 mL	CGOS1-500ML
Palladium, Pd	HCl	30 mL	CGPD1-30ML
		125 mL	CGPD1-125ML
		500 mL	CGPD1-500ML
Palladium, Pd	HNO ₃	30 mL	CGPDN1-30ML
		125 mL	CGPDN1-125ML
		500 mL	CGPDN1-500ML
Phosphorus, P	H ₂ O	30 mL	CGP1-30ML
		125 mL	CGP1-125ML
		500 mL	CGP1-500ML
Platinum, Pt Contains Chloride	HNO ₃	30 mL	CGPTN1-30ML
		125 mL	CGPTN1-125ML
		500 mL	CGPTN1-500ML
Platinum, Pt	HCl	30 mL	CGPT1-30ML
		125 mL	CGPT1-125ML
		500 mL	CGPT1-500ML
Platinum, Pt Chloride Free	HNO ₃	30 mL	CGPTNO31-30ML
		125 mL	CGPTNO31-125ML
		500 mL	CGPTNO31-500ML
Potassium, K	HNO ₃	30 mL	CGK1-30ML
		125 mL	CGK1-125ML
		500 mL	CGK1-500ML
Praseodymium, Pr	HNO ₃	30 mL	CGPR1-30ML
		125 mL	CGPR1-125ML
		500 mL	CGPR1-500ML
Rhenium, Re	HNO ₃	30 mL	CGRE1-30ML
		125 mL	CGRE1-125ML
		500 mL	CGRE1-500ML
Rhodium, Rh Commonly used as an Internal Standard for ICP-MS.	HCl	30 mL	CGRH1-30ML
		125 mL	CGRH1-125ML
		500 mL	CGRH1-500ML
Rhodium, Rh Commonly used as an Internal Standard for ICP-MS.	HNO ₃	30 mL	CGRHN1-30ML
		125 mL	CGRHN1-125ML
		500 mL	CGRHN1-500ML
Rubidium, Rb	HNO ₃	30 mL	CGRB1-30ML
		125 mL	CGRB1-125ML
Ruthenium, Ru	HCl	30 mL	CGRU1-30ML
		125 mL	CGRU1-125ML
		500 mL	CGRU1-500ML
Samarium, Sm	HNO ₃	30 mL	CGSM1-30ML
		125 mL	CGSM1-125ML
		500 mL	CGSM1-500ML

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Scandium, Sc Commonly used as an Internal Standard for ICP-MS.	HNO ₃	30 mL	CGSC1-30ML
		125 mL	CGSC1-125ML
		500 mL	CGSC1-500ML
Selenium⁺⁴, Se⁺⁴	HNO ₃	30 mL	CGSE(4)1-30ML
		125 mL	CGSE(4)1-125ML
		500 mL	CGSE(4)1-500ML
Selenium⁺⁶, Se⁺⁶	H ₂ O	30 mL	CGSE(6)1-30ML
		125 mL	CGSE(6)1-125ML
Silica, SiO₂	HNO ₃ / HF	30 mL	CGSIO1-30ML
		125 mL	CGSIO1-125ML
		500 mL	CGSIO1-500ML
Silica, SiO₂	NaOH	125 mL	CGSIONA1-125ML
		500 mL	CGSIONA1-500ML
Silicon, Si	HNO ₃ / HF	30 mL	CGSI1-30ML
		125 mL	CGSI1-125ML
		500 mL	CGSI1-500ML
Silicon, Si	NaOH	125 mL	CGSINA1-125ML
		500 mL	CGSINA1-500ML
Silver, Ag	HNO ₃	30 mL	CGAG1-30ML
		125 mL	CGAG1-125ML
		500 mL	CGAG1-500ML
Sodium, Na	HNO ₃	30 mL	CGNA1-30ML
		125 mL	CGNA1-125ML
		500 mL	CGNA1-500ML
Strontium, Sr	HNO ₃	30 mL	CGSR1-30ML
		125 mL	CGSR1-125ML
		500 mL	CGSR1-500ML
Sulfur, S Compatible if mixed with Ba and Pb.	H ₂ O	125 mL	CGMSA1-125ML
		500 mL	CGMSA1-500ML
Sulfur, S	H ₂ O	30 mL	CGS1-30ML
		125 mL	CGS1-125ML
		500 mL	CGS1-500ML
Tantalum, Ta	HNO ₃ / HF	30 mL	CGTA1-30ML
		125 mL	CGTA1-125ML
		500 mL	CGTA1-500ML
Tellurium, Te	HCl	30 mL	CGTE1-30ML
		125 mL	CGTE1-125ML
		500 mL	CGTE1-500ML
Tellurium, Te	HNO ₃	30 mL	CGTEN1-30ML
		125 mL	CGTEN1-125ML
		500 mL	CGTEN1-500ML
Terbium, Tb Commonly used as an Internal Standard for ICP-MS.	HNO ₃	30 mL	CGTB1-30ML
		125 mL	CGTB1-125ML
		500 mL	CGTB1-500ML

SINGLE-ELEMENT STANDARDS

1,000 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Thallium, Tl	HNO ₃	30 mL	CGTL1-30ML
		125 mL	CGTL1-125ML
		500 mL	CGTL1-500ML
Thorium, Th	HNO ₃	30 mL	CGTH1-30ML
		125 mL	CGTH1-125ML
		500 mL	CGTH1-500ML
Thulium, Tm	HNO ₃	30 mL	CGTM1-30ML
		125 mL	CGTM1-125ML
		500 mL	CGTM1-500ML
Tin, Sn	HCl	30 mL	CGSNCL1-30ML
		125 mL	CGSNCL1-125ML
Tin, Sn	HNO ₃ / HF	30 mL	CGSN1-30ML
		125 mL	CGSN1-125ML
		500 mL	CGSN1-500ML
Titanium, Ti	HNO ₃ / HF	30 mL	CGT11-30ML
		125 mL	CGT11-125ML
		500 mL	CGT11-500ML
Tungsten, W	HNO ₃ / HF	30 mL	CGW1-30ML
		125 mL	CGW1-125ML
		500 mL	CGW1-500ML
Tungsten, W	H ₂ O	125 mL	CGWH201-125ML
Uranium, U	HNO ₃	30 mL	CGU1-30ML
		125 mL	CGU1-125ML
		500 mL	CGU1-500ML
Vanadium, V	HNO ₃	30 mL	CGV1-30ML
		125 mL	CGV1-125ML
		500 mL	CGV1-500ML
Ytterbium, Yb	HNO ₃	30 mL	CGYB1-30ML
		125 mL	CGYB1-125ML
		500 mL	CGYB1-500ML
Yttrium, Y <small>Commonly used as an Internal Standard for ICP-MS.</small>	HNO ₃	30 mL	CGY1-30ML
		125 mL	CGY1-125ML
		500 mL	CGY1-500ML
Zinc, Zn	HNO ₃	30 mL	CGZN1-30ML
		125 mL	CGZN1-125ML
		500 mL	CGZN1-500ML
Zirconium, Zr	HF	30 mL	CGZR1-30ML
		125 mL	CGZR1-125ML
		500 mL	CGZR1-500ML

See pg. 34 for our HF-free Zirconium, part number CGZRCL10-125ML or CGZRCL10-500ML.

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

10,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Aluminum, Al	HNO ₃	30 mL	CGAL10-30ML
		125 mL	CGAL10-125ML
		500 mL	CGAL10-500ML
Antimony, Sb	HNO ₃ / Tartaric Acid	125 mL	CGSB10-125ML
		500 mL	CGSB10-500ML
Arsenic, As	HNO ₃	30 mL	CGAS10-30ML
		125 mL	CGAS10-125ML
		500 mL	CGAS10-500ML
Barium, Ba	HNO ₃	125 mL	CGBA10-125ML
		500 mL	CGBA10-500ML
Beryllium, Be	HNO ₃	125 mL	CGBE10-125ML
		500 mL	CGBE10-500ML
Bismuth, Bi	HNO ₃	30 mL	CGBI10-30ML
		125 mL	CGBI10-125ML
		500 mL	CGBI10-500ML
Boron, B	NH ₄ OH	125 mL	CGB10-125ML
		500 mL	CGB10-500ML
Cadmium, Cd	HNO ₃	125 mL	CGCD10-125ML
		500 mL	CGCD10-500ML
Calcium, Ca	HNO ₃	30 mL	CGCA10-30ML
		125 mL	CGCA10-125ML
		500 mL	CGCA10-500ML
Carbon, C	HNO ₃	125 mL	CGC10-125ML
Cerium, Ce	HNO ₃	30 mL	CGCE10-30ML
		125 mL	CGCE10-125ML
		500 mL	CGCE10-500ML
Cesium, Cs	HNO ₃	125 mL	CGCS10-125ML
		500 mL	CGCS10-500ML
Chromium ⁺³ , Cr ⁺³	HNO ₃	30 mL	CGCR(3)10-30ML
		125 mL	CGCR(3)10-125ML
		500 mL	CGCR(3)10-500ML
Cobalt, Co <small>Commonly used as an Internal Standard for ICP-OES.</small>	HNO ₃	30 mL	CGCO10-30ML
		125 mL	CGCO10-125ML
		500 mL	CGCO10-500ML
Copper, Cu	HNO ₃	30 mL	CGCU10-30ML
		125 mL	CGCU10-125ML
		500 mL	CGCU10-500ML
Dysprosium, Dy	HNO ₃	30 mL	CGDY10-30ML
		125 mL	CGDY10-125ML
		500 mL	CGDY10-500ML
Erbium, Er	HNO ₃	30 mL	CGER10-30ML
		125 mL	CGER10-125ML
		500 mL	CGER10-500ML
Europium, Eu	HNO ₃	30 mL	CGEU10-30ML
		125 mL	CGEU10-125ML
		500 mL	CGEU10-500ML

SINGLE-ELEMENT STANDARDS

10,000 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

10,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Gadolinium, Gd	HNO ₃	30 mL	CGGD10-30ML
		125 mL	CGGD10-125ML
		500 mL	CGGD10-500ML
Gallium, Ga	HNO ₃	125 mL	CGGA10-125ML
		500 mL	CGGA10-500ML
Germanium, Ge	HNO ₃ / HF	125 mL	CGGE10-125ML
		500 mL	CGGE10-500ML
Gold, Au	HCl	30 mL	CGAU10-30ML
		125 mL	CGAU10-125ML
		500 mL	CGAU10-500ML
Hafnium, Hf	HNO ₃ / HF	125 mL	CGHF10-125ML
		500 mL	CGHF10-500ML
Holmium, Ho	HNO ₃	30 mL	CGHO10-30ML
		125 mL	CGHO10-125ML
		500 mL	CGHO10-500ML
Indium, In <small>Commonly used as an Internal Standard for ICP-OES.</small>	HNO ₃	125 mL	CGIN10-125ML
		500 mL	CGIN10-500ML
Iridium, Ir	HCl	30 mL	CGIR10-30ML
		125 mL	CGIR10-125ML
		500 mL	CGIR10-500ML
Iron, Fe	HNO ₃	30 mL	CGFE10-30ML
		125 mL	CGFE10-125ML
		500 mL	CGFE10-500ML
Lanthanum, La	HNO ₃	30 mL	CGLA10-30ML
		125 mL	CGLA10-125ML
		500 mL	CGLA10-500ML
Lead, Pb	HNO ₃	30 mL	CGPB10-30ML
		125 mL	CGPB10-125ML
		500 mL	CGPB10-500ML
Lithium, Li	HNO ₃	30 mL	CGLI10-30ML
		125 mL	CGLI10-125ML
		500 mL	CGLI10-500ML
Lutetium, Lu	HNO ₃	30 mL	CGLU10-30ML
		125 mL	CGLU10-125ML
		500 mL	CGLU10-500ML
Magnesium, Mg	HNO ₃	30 mL	CGMG10-30ML
		125 mL	CGMG10-125ML
		500 mL	CGMG10-500ML
Manganese, Mn	HNO ₃	30 mL	CGMN10-30ML
		125 mL	CGMN10-125ML
		500 mL	CGMN10-500ML
Mercury, Hg	HNO ₃	125 mL	CGHG10-125ML
		500 mL	CGHG10-500ML
Molybdenum, Mo	NH ₄ OH	30 mL	CGMO10-30ML
		125 mL	CGMO10-125ML
		500 mL	CGMO10-500ML

SINGLE-ELEMENT STANDARDS

10,000 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

10,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Neodymium, Nd	HNO ₃	30 mL	CGND10-30ML
		125 mL	CGND10-125ML
		500 mL	CGND10-500ML
Nickel, Ni	HNO ₃	30 mL	CGNI10-30ML
		125 mL	CGNI10-125ML
		500 mL	CGNI10-500ML
Niobium, Nb	HNO ₃ / HF	125 mL	CGNB10-125ML
Niobium, Nb High purity, low Tantalum	HNO ₃ / HF	125 mL	CGNB20510-125ML
Palladium, Pd	HCl	30 mL	CGPD10-30ML
		125 mL	CGPD10-125ML
		500 mL	CGPD10-500ML
Phosphorus, P	H ₂ O	30 mL	CGP10-30ML
		125 mL	CGP10-125ML
		500 mL	CGP10-500ML
Platinum, Pt	HCl	30 mL	CGPT10-30ML
		125 mL	CGPT10-125ML
		500 mL	CGPT10-500ML
Potassium, K	HNO ₃	30 mL	CGK10-30ML
		125 mL	CGK10-125ML
		500 mL	CGK10-500ML
Praseodymium, Pr	HNO ₃	30 mL	CGPR10-30ML
		125 mL	CGPR10-125ML
		500 mL	CGPR10-500ML
Rhenium, Re	HNO ₃	125 mL	CGRE10-125ML
		500 mL	CGRE10-500ML
Rhodium, Rh	HCl	30 mL	CGRH10-30ML
		125 mL	CGRH10-125ML
		500 mL	CGRH10-500ML
Rubidium, Rb	HNO ₃	125 mL	CGRB10-125ML
		500 mL	CGRB10-500ML
Ruthenium, Ru	HCl	30 mL	CGRU10-30ML
		125 mL	CGRU10-125ML
		500 mL	CGRU10-500ML
Samarium, Sm	HNO ₃	30 mL	CGSM10-30ML
		125 mL	CGSM10-125ML
Scandium, Sc Commonly used as an Internal Standard for ICP-OES.	HNO ₃	30 mL	CGSC10-30ML
		125 mL	CGSC10-125ML
		500 mL	CGSC10-500ML
Selenium, Se	HNO ₃	30 mL	CGSE10-30ML
		125 mL	CGSE10-125ML
		500 mL	CGSE10-500ML
Silicon, Si	HNO ₃ / HF	30 mL	CGSI10-30ML
		125 mL	CGSI10-125ML
		500 mL	CGSI10-500ML

SINGLE-ELEMENT STANDARDS

10,000 µg/mL Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

10,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Silver, Ag	HNO ₃	125 mL	CGAG10-125ML
		500 mL	CGAG10-500ML
Sodium, Na	HNO ₃	30 mL	CGNA10-30ML
		125 mL	CGNA10-125ML
		500 mL	CGNA10-500ML
Strontium, Sr	HNO ₃	125 mL	CGSR10-125ML
		500 mL	CGSR10-500ML
Sulfur, S Compatible if mixed with Ba and Pb.	H ₂ O	125 mL	CGMSA10-125ML
		500 mL	CGMSA10-500ML
Sulfur, S	H ₂ O	30 mL	CGS10-30ML
		125 mL	CGS10-125ML
		500 mL	CGS10-500ML
Tantalum, Ta	HNO ₃ / HF	125 mL	CGTA10-125ML
Tellurium, Te	HCl	125 mL	CGTE10-125ML
		500 mL	CGTE10-500ML
Terbium, Tb	HNO ₃	30 mL	CGTB10-30ML
		125 mL	CGTB10-125ML
		500 mL	CGTB10-500ML
Thallium, Tl	HNO ₃	125 mL	CGTL10-125ML
		500 mL	CGTL10-500ML
Thorium, Th	HNO ₃	125 mL	CGTH10-125ML
Thulium, Tm	HNO ₃	30 mL	CGTM10-30ML
		125 mL	CGTM10-125ML
Tin, Sn	HNO ₃ / HF	30 mL	CGSN10-30ML
		125 mL	CGSN10-125ML
		500 mL	CGSN10-500ML
Titanium, Ti	HNO ₃ / HF	30 mL	CGT110-30ML
		125 mL	CGT110-125ML
		500 mL	CGT110-500ML
Tungsten, W	HNO ₃ / HF	125 mL	CGW10-125ML
		500 mL	CGW10-500ML
Uranium, U	HNO ₃	30 mL	CGU10-30ML
		125 mL	CGU10-125ML
		500 mL	CGU10-500ML
Vanadium, V	HNO ₃	30 mL	CGV10-30ML
		125 mL	CGV10-125ML
		500 mL	CGV10-500ML
Ytterbium, Yb	HNO ₃	30 mL	CGYB10-30ML
		125 mL	CGYB10-125ML
		500 mL	CGYB10-500ML
Yttrium, Y Commonly used as an Internal Standard for ICP-OES.	HNO ₃	30 mL	CGY10-30ML
		125 mL	CGY10-125ML
		500 mL	CGY10-500ML

SINGLE-ELEMENT STANDARDS

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Need a specific starting material or matrix? Custom 10 ppm single element solutions available upon request. Need a multielement solution? We can do that too!

10,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Zinc, Zn	HNO ₃	30 mL	CGZN10-30ML
		125 mL	CGZN10-125ML
		500 mL	CGZN10-500ML
Zirconium, Zr	HF	125 mL	CGZR10-125ML
Zirconium, Zr HF free	HCl	30 mL	CGZRCL10-30ML
		125 mL	CGZRCL10-125ML
		500 mL	CGZRCL10-500ML

CYANIDE STANDARDS

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Certificate of Analysis includes lot specific trace metal impurity analysis.

Custom cyanide standards are available upon request.

1,000 µg/mL

ANALYTE	µg/mL	MATRIX	VOLUME	CATALOG #
Copper, Cu	1,000	NaCN	125 mL	AACUCN-125ML
			500 mL	AACUCN-500ML
Gold, Au	1,000	NaCN	125 mL	AAAUCN-125ML
			500 mL	AAAUCN-500ML
Silver, Ag	1,000	NaCN	125 mL	AAAGCN-125ML
			500 mL	AAAGCN-500ML
Zinc, Zn	1,000	NaCN	125 mL	AAZNCN-125ML
			500 mL	AAZNCN-500ML

MULTI-ELEMENT STANDARDS

Identical or near identical formulations | Multi-Element Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

ICP Calibration Standard			
IV-STOCK-2 C		Matrix: HNO ₃	
IV-STOCK-2-125ML		Volume: 125 mL	
IV-STOCK-2-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ca	10,000	Mg	10,000
K	10,000	Na	10,000

ICP Calibration Standard			
IV-STOCK-3 C PE		Matrix: HNO ₃	
IV-STOCK-3-125ML		Volume: 125 mL	
IV-STOCK-3-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ca	1,000	Mg	1,000
K	1,000	Na	1,000

ICP Calibration Standard			
IV-STOCK-4 C JY M PE		Matrix: HNO ₃	
IV-STOCK-4-125ML		Volume: 125 mL	
IV-STOCK-4-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	1,000	In	1,000
Al	1,000	K	1,000
B	1,000	Li	1,000
Ba	1,000	Mg	1,000
Bi	1,000	Mn	1,000
Ca	1,000	Na	1,000
Cd	1,000	Ni	1,000
Co	1,000	Pb	1,000
Cr	1,000	Sr	1,000
Cu	1,000	Tl	1,000
Fe	1,000	Zn	1,000
Ga	1,000		

Wavelength Calibration Standard			
IV-STOCK-5 M		Matrix: HCl / HF	
IV-STOCK-5-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	20	Mg	1
As	20	Mn	1
B	2	Na	20
Ba	2	Ni	5
Be	1	P	10
Ca	10	Pb	20
Cd	2	Sc	1
Cr	2	Se	20
Cu	2	Sr	1
Fe	2	Te	20
Hg	5	Ti	2
K	100	Y	1
Li	2	Zn	2

ICP Calibration Standard			
IV-STOCK-6 M PE		Matrix: HNO ₃	
IV-STOCK-6-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Li	10
Al	10	Mg	10
As	100	Mn	10
B	100	Mo	10
Ba	10	Na	10
Be	100	Ni	10
Bi	10	Pb	10
Ca	1,000	Rb	10
Cd	10	Se	100
Co	10	Sr	10
Cr	10	Te	10
Cu	10	Tl	10
Fe	100	U	10
Ga	10	V	10
K	10	Zn	100

C Common Multi-Element Standards

JY HORIBA Jobin Yvon

PE Perkin Elmer

M Merck/MilliporeSigma

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

Cation Calibration Standard			
IV-STOCK-7 		Matrix: HNO ₃	
IV-STOCK-7-125ML		Volume: 125 mL	
IV-STOCK-7-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba ²⁺	100	Mn ²⁺	100
Ca ²⁺	100	Na ⁺	100
K ⁺	100	NH ₄ ⁺	100
Li ⁺	100	Sr ²⁺	100
Mg ²⁺	100		

ICP Calibration Standard			
IV-STOCK-8 		Matrix: HNO ₃	
IV-STOCK-8-125ML		Volume: 125 mL	
IV-STOCK-8-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	100	K	100
B	100	Li	100
Ba	100	Mg	100
Be	100	Mn	100
Bi	100	Na	100
Ca	100	Ni	100
Cd	100	Pb	100
Co	100	Se	100
Cr	100	Sr	100
Cu	100	Te	100
Fe	100	Tl	100
Ga	100	Zn	100

ICP Calibration Standard – Toxic Elements			
IV-STOCK-9		Matrix: HNO ₃	
IV-STOCK-9-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	100	Pb	100
Be	100	Se	100
Cd	100	Tl	100
Ni	100		

 Merck/MilliporeSigma

 Common Multi-Ion Standards

 Perkin Elmer

ICP Calibration Standard – Surface Water			
IV-STOCK-10 		Matrix: HNO ₃	
IV-STOCK-10-125ML		Volume: 125 mL	
Analyte	µg/L*	Analyte	µg/L*
As	50	Mg	15,000
B	100	Mn	30
Ba	50	Mo	100
Be	20	Na	8,000
Bi	10	Ni	50
Ca	35,000	Pb	25
Cd	20	Se	10
Co	25	Sr	100
Cr	20	Tl	10
Cu	20	V	50
Fe	100	Zn	50
K	3,000		

*Parts per billion

ICP-MS Calibration Standard			
IV-STOCK-12		Matrix: HNO ₃	
IV-STOCK-12-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba	10	In	10
Be	10	Li	10
Bi	10	Ni	10
Ce	10	Pb	10
Co	10	U	10

ICP Calibration Standard – Trace Metals			
IV-STOCK-13 		Matrix: HNO ₃	
IV-STOCK-13-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	500	Fe	100
As	100	Mn	100
Be	100	Ni	100
Cd	25	Pb	100
Co	100	Se	25
Cr	100	V	250
Cu	100	Zn	100

MULTI-ELEMENT STANDARDS

Identical or near identical formulations | Multi-Element Standards

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

Wavelength Calibration Standard			
IV-STOCK-14 C M PE		Matrix: HCl / HNO ₃ / HF	
IV-STOCK-14-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
As	20	Na	20
K	100	Ni	20
La	20	P	100
Li	20	S	100
Mn	20	Sc	20
Mo	20		

ICP-MS Calibration Standard			
IV-STOCK-15		Matrix: HNO ₃	
IV-STOCK-15-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ca	10	Li	10
Fe	10	Na	10
K	10		

ICP Calibration Standard – Alkaline Earth Element			
IV-STOCK-16 M		Matrix: HNO ₃	
IV-STOCK-16-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba	1,000	Mg	1,000
Ca	1,000	Sr	1,000

ICP Calibration Standard – HCl Soluble Elements			
IV-STOCK-17 M		Matrix: HCl/HNO ₃ /HF	
IV-STOCK-17-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Hf	100	Ta	100
Ir	100	Ti	100
Sb	100	Zr	100
Sn	100		

GFAA Calibration Standard			
IV-STOCK-18 M PE		Matrix: HNO ₃	
IV-STOCK-18-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Cu	50
Al	100	Fe	20
As	100	Mn	20
Ba	50	Ni	50
Be	5	Pb	100
Cd	5	Sb	100
Co	50	Se	100
Cr	20	Tl	100

ICP Calibration Standard			
IV-STOCK-21 M PE		Matrix: HNO ₃	
IV-STOCK-21-125ML		Volume: 125 mL	
IV-STOCK-21-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	K	10
Al	10	Li	10
As	10	Mg	10
Ba	10	Mn	10
Be	10	Na	10
Bi	10	Ni	10
Ca	10	Pb	10
Cd	10	Rb	10
Co	10	Se	10
Cr	10	Sr	10
Cs	10	Tl	10
Cu	10	U	10
Fe	10	V	10
Ga	10	Zn	10
In	10		

C Common Multi-Element Standards

M Merck/MilliporeSigma

PE Perkin Elmer

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

ICP Calibration Standard			
IV-STOCK-22 		Matrix: HNO ₃	
IV-STOCK-22-125ML		Volume: 125 mL	
Analyte	µg/L*	Analyte	µg/L*
Cd	200	Pb	200
Cu	200	Rh	200
Mg	200		

*Parts per billion

ICP Calibration Standard			
IV-STOCK-23 		Matrix: HNO ₃	
IV-STOCK-23-500ML		Volume: 500 mL	
Analyte	µg/L*	Analyte	µg/L*
B	1	Lu	1
Ba	1	Na	1
Co	1	Rh	1
Fe	1	Sc	1
Ga	1	Tl	1
In	1	U	1
K	1	Y	1
Li	1		

*Parts per billion

Tuning Solution			
IV-STOCK-24   		Matrix: HNO ₃	
IV-STOCK-24-125ML		Volume: 125 mL	
IV-STOCK-24-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	50	Mn	50
As	50	Mo	50
Ba	50	Ni	50
Cd	50	Pb	50
Co	50	Se	50
Cr	50	Sr	50
Cu	50	Zn	50
K	500		

ICP Calibration Standard			
IV-STOCK-26  		Matrix: HNO ₃	
IV-STOCK-26-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ce	10	Pr	10
Dy	10	Sc	10
Er	10	Sm	10
Eu	10	Tb	10
Gd	10	Th	10
Ho	10	Tm	10
La	10	Y	10
Lu	10	Yb	10
Nd	10		

ICP Calibration Standard			
IV-STOCK-27 		Matrix: HNO ₃	
IV-STOCK-27-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Li	10
Al	10	Mg	10
As	10	Mn	10
Ba	10	Na	10
Be	10	Ni	10
Ca	10	Pb	10
Cd	10	Rb	10
Co	10	Se	10
Cr	10	Sr	10
Cs	10	Tl	10
Cu	10	U	10
Fe	10	V	10
Ga	10	Zn	10
K	10		

 Agilent/Varian Common Multi-Element Standards Merck/MilliporeSigma Perkin Elmer

MULTI-ELEMENT STANDARDS

Multi-Element Standards

Identical or near identical formulations

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

ICP Calibration Standard			
IV-STOCK-28 AV PE		Matrix: HCl / HNO ₃	
IV-STOCK-28-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	10	Rh	10
Hf	10	Ru	10
Ir	10	Sb	10
Pd	10	Sn	10
Pt	10	Te	10

ICP Calibration Standard			
IV-STOCK-29 AV PE		Matrix: HNO ₃ / HF	
IV-STOCK-29-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
B	10	S	10
Ge	10	Si	10
Mo	10	Ta	10
Nb	10	Ti	10
P	10	W	10
Re	10	Zr	10

ICP Calibration Standard			
IV-STOCK-30 PE		Matrix: HNO ₃	
IV-STOCK-30-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Be	10	Mg	10
Bi	10	Ni	10
Ce	10	Pb	10
Co	10	U	10
In	10		

AV Agilent/Varian**C** Common Multi-Element Standards**JY** HORIBA Jobin Yvon**PE** Perkin Elmer**T** Thermo Scientific

ICP Calibration Standard			
IV-STOCK-31 C T		Matrix: HNO ₃	
IV-STOCK-31-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	1	Mg	0.2
Ba	0.2	Mn	1
Ca	0.2	Ni	5
Cu	1	P	10
K	5	Zn	0.2

Calibration Standard – Mix Majors			
IV-STOCK-33 AV		Matrix: HNO ₃	
IV-STOCK-33-125ML		Volume: 125 mL	
IV-STOCK-33-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ca	500	Mg	500
Fe	500	Na	500
K	500		

ICP Calibration Standard			
IV-STOCK-34 AV C JY PE		Matrix: HNO ₃	
IV-STOCK-34-125ML		Volume: 125 mL	
IV-STOCK-34-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ca	5,000	Mg	5,000
K	5,000	Na	5,000

ICP Calibration Standard			
IV-STOCK-35 PE		Matrix: HNO ₃	
IV-STOCK-35-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ca	1,000	Mg	1,000
Fe	1,000	Na	1,000
K	1,000		

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

ICP Calibration Standard			
IV-STOCK-36 C		Matrix: HCl	
IV-STOCK-36-125ML		Volume: 125 mL	
IV-STOCK-36-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	100	Pt	100
Pd	100		

Environmental Calibration Standard			
IV-STOCK-50 AV		Matrix: HNO ₃ / HF	
IV-STOCK-50-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Mn	10
Al	10	Mo	10
As	10	Na	1,000
Ba	10	Ni	10
Be	10	Pb	10
Ca	1,000	Sb	10
Cd	10	Se	10
Co	10	Th	10
Cr	10	Tl	10
Cu	10	U	10
Fe	1,000	V	10
K	1,000	Zn	10
Mg	1,000		

Internal Standard			
IV-STOCK-53 AV PE		Matrix: HNO ₃ / HF	
IV-STOCK-53-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	10	Sc	10
Ge	10	Tb	10
In	10	Y	10
⁶ Li	10		

- AV** Agilent/Varian
C Common Multi-Element Standards
PE Perkin Elmer

7500 Series PA Tuning Solution 1 (commonly used with IV-Stock-52)			
IV-STOCK-51 AV		Matrix: HNO ₃	
IV-STOCK-51-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	5	Mn	5
As	20	Na	5
Ba	5	Ni	10
Be	20	Pb	10
Bi	5	Sc	5
Cd	20	Sr	5
Co	5	Th	5
Cr	5	Tl	5
Cu	5	U	5
In	5	V	5
⁶ Li	5	Y	2.5
Lu	5	Yb	2.5
Mg	10	Zn	20

7500 Series PA Tuning Solution 2 (commonly used with IV-Stock-51)			
IV-STOCK-52		Matrix: HCl	
IV-STOCK-52-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ge	10	Ru	10
Ir	5	Sb	10
Mo	10	Sn	10
Pd	10	Ti	5

Interference Check Standard			
IV-STOCK-54 PE		Matrix: HNO ₃	
IV-STOCK-54-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	1,200	Mg	3,000
Ca	6,000	Na	1,000
Fe	5,000		

MULTI-ELEMENT STANDARDS

Multi-Element Standards

Identical or near identical formulations |

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

Wavecal Standard			
IV-STOCK-55 M PE		Matrix: HNO ₃	
IV-STOCK-55-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba	1	Li	10
Ca	1	Mn	10
K	50	Na	10
La	10	Sr	10

ICP Calibration Standard			
IV-STOCK-56 C		Matrix: HNO ₃ / HF	
IV-STOCK-56-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Mo	100	Sn	100
Sb	100	Ti	100
Si	100		

ICP Calibration Standard			
IV-STOCK-57 C		Matrix: HNO ₃ / HF	
IV-STOCK-57-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Mo	10	Sn	10
Sb	10	Ti	10

ICP Calibration Standard			
IV-STOCK-58 C		Matrix: HCl	
IV-STOCK-58-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	100	Pt	100
Ir	100	Re	100
Os	100	Rh	100
Pd	100	Ru	100

- AV** Agilent/Varian
C Common Multi-Element Standards
M Merck/MilliporeSigma

- PE** Perkin Elmer
T Thermo Scientific

ICP-MS Tuning Solution			
IV-STOCK-74 AV		Matrix: HNO ₃	
IV-STOCK-74-500ML		Volume: 500 mL	
Analyte	µg/L*	Analyte	µg/L*
Ce	1	Mg	1
Co	1	Tl	1
Li	1	Y	1

*Parts per billion

ICP-MS Internal Standard			
IV-STOCK-75 AV		Matrix: HNO ₃ / HF	
IV-STOCK-75-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	100	Lu	100
Ge	100	Rh	100
In	100	Sc	100
⁶ Li	100	Tb	100

ICP-MS Tuning Solution			
IV-STOCK-77 PE		Matrix: HNO ₃	
IV-STOCK-77-500ML		Volume: 500 mL	
Analyte	µg/L*	Analyte	µg/L*
Be	1	Li	1
Ce	1	Mg	1
Fe	1	Pb	1
In	1	U	1

*Parts per billion

Multi-Element ICP-MS Tuning Solution for Thermo iCAP TQ			
IV-STOCK-80 T		Matrix: 2% v/v HNO ₃ / 0.5% v/v/HCl	
IV-STOCK-80-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba	1	Li	1
Bi	1	Mg	1
Ce	1	Ti	1
Co	1	U	1
Ho	1	Y	
In	1		

These elements are grouped for ease of use. Intended for ICP-MS and ICP-OES, they can be used individually or in any combination upon dilution into 1% HNO₃. Custom ICP-MS/OES calibration standards are available upon request.

65-Element Group

Rare Earth ICP-MS Standard			
CMS-1 C		Matrix: HNO ₃	
CMS-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ce	10	Pr	10
Dy	10	Sc	10
Er	10	Sm	10
Eu	10	Tb	10
Gd	10	Th	10
Ho	10	Tm	10
La	10	U	10
Lu	10	Y	10
Nd	10	Yb	10

For ICP analysis of all rare earth elements plus U and Th.

Precious Metals ICP-MS Standard			
CMS-2 C		Matrix: HCl	
CMS-2-125ML		Volume: 125 mL	
CMS-2-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	10	Re	10
Ir	10	Rh	10
Pd	10	Ru	10
Pt	10	Te	10

For ICP analysis of precious metals plus Re and Te.

Fluoride Soluble ICP-MS Standard			
CMS-3 C		Matrix: HNO ₃ /HF	
CMS-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ge	10	Ta	10
Hf	10	Ti	10
Mo	10	W	10
Nb	10	Zr	10
Sn	10		

For ICP analysis of elements that tolerate or require HF for stability

Hot Plasma ICP-MS Complete Standard			
CMS-4 C		Matrix: HNO ₃	
CMS-4-125ML		Volume: 125 mL	
CMS-4-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
As	10	In	10
B	10	Pb	10
Ba	10	Sb	10
Be	10	Se	10
Bi	10	Tl	10
Cd	10	V	10
Ga	10		

For direct use in ICP analysis or as stock concentrate.

Cool Plasma ICP-MS Complete Standard			
CMS-5 C		Matrix: HNO ₃	
CMS-5-125ML		Volume: 125 mL	
CMS-5-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Li	10
Al	10	Mg	10
Ca	10	Mn	10
Co	10	Na	10
Cr	10	Ni	10
Cs	10	Rb	10
Cu	10	Sr	10
Fe	10	Zn	10
K	10		

For direct use in ICP analysis or as stock concentrate.

C Common Multi-Element Standard

MULTI-ELEMENT STANDARDS

Multi-Element Standards

Identical or near identical formulations

These elements are grouped for ease of use. Intended for ICP-MS and ICP-OES, they can be used individually or in any combination upon dilution into 1% HNO₃. Custom ICP-MS/OES calibration standards are available upon request.

69-Element Group

Rare Earth ICP-MS Standard			
CCS-1 C		Matrix: HNO ₃	
CCS-1-125ML		Volume: 125 mL	
CCS-1-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ce	100	Pr	100
Dy	100	Sc	100
Er	100	Sm	100
Eu	100	Tb	100
Gd	100	Th	100
Ho	100	Tm	100
La	100	U	100
Lu	100	Y	100
Nd	100	Yb	100

For ICP analysis of all rare earth elements plus U and Th. Uranium is isotopically depleted. Can be diluted with CCS-4 and CCS-6 to working concentrations.

Precious Metals ICP-MS Standard			
CCS-2 C		Matrix: HCl	
CCS-2-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	100	Pt	100
Ir	100	Rh	100
Pd	100	Ru	100

For simultaneous ICP analysis of precious metals. Can be diluted with CCS-1 or CCS-5 to working concentrations. For dilution with CCS-6 please see Silver Chemical Stability article for more information about Ag stability in HCl.

Alkali, Alkaline, Non-Transition ICP-MS Standard			
CCS-4 C		Matrix: HNO ₃	
CCS-4-125ML		Volume: 125 mL	
CCS-4-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	100	In	100
As	100	K	100
Ba	100	Li	100
Be	100	Mg	100
Bi	100	Na	100
Ca	100	Rb	100
Cs	100	Se	100
Ga	100	Sr	100

For use as stock concentrate for ICP analysis. Can be diluted with CCS-1 and CCS-6 to working concentrations.

Fluoride Soluble ICP-MS Standard			
CCS-5 C		Matrix: HNO ₃ /HF	
CCS-5-125ML		Volume: 125 mL	
CCS-5-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
B	100	Sb	100
Ge	100	Si	100
Hf	100	Sn	100
Mo	100	Ta	100
Nb	100	Ti	100
P	100	W	100
Re	100	Zr	100
S	100		

For ICP analysis of elements that tolerate or require HF for stability. Can be diluted with CCS-2 and CCS-6 to working concentrations. Can be diluted with CCS-4 to lower working concentrations (<10 ppm recommended).

C Common Multi-Element Standard

These elements are grouped for ease of use. Intended for ICP-MS and ICP-OES, they can be used individually or in any combination upon dilution into 1% HNO₃. Custom ICP-MS/OES calibration standards are available upon request.

69-Element Group

Transition ICP-MS Standard			
CCS-6 C		Matrix: HNO ₃	
CCS-6-125ML CCS-6-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	Mn	100
Cd	100	Ni	100
Co	100	Pb	100
Cr	100	Tl	100
Cu	100	V	100
Fe	100	Zn	100
Hg	100		

For use as stock concentrate for ICP analysis. Can be diluted with CCS-1 and CCS-4 to working concentrations. Contains mercury (Hg); please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

ICP-MS Refractory Elements Standard			
IV-ICPMS-71B C		Matrix: HNO ₃ / HF	
IV-ICPMS-71B-125ML IV-ICPMS-71B-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ge	10	Sn	10
Hf	10	Ta	10
Mo	10	Te	10
Nb	10	Ti	10
Sb	10	W	10
Si	10	Zr	10

Can be diluted to working concentrations without additional HF for stability.

71-Element Group

ICP-MS Complete Standard			
IV-ICPMS-71A C		Matrix: HNO ₃	
IV-ICPMS-71A-125ML IV-ICPMS-71A-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Lu	10
Al	10	Mg	10
As	10	Mn	10
B	10	Na	10
Ba	10	Nd	10
Be	10	Ni	10
Ca	10	P	10
Cd	10	Pb	10
Ce	10	Pr	10
Co	10	Rb	10
Cr	10	S	10
Cs	10	Se	10
Cu	10	Sm	10
Dy	10	Sr	10
Er	10	Th	10
Eu	10	Tl	10
Fe	10	Tm	10
Ga	10	U	10
Gd	10	V	10
Ho	10	Yb	10
K	10	Zn	10
La	10		

Uranium is isotopically depleted. Can be diluted with other multi-element standards to working concentrations.

C Common Multi-Element Standard

MULTI-ELEMENT STANDARDS

Multi-Element Standards

Identical or near identical formulations |

These elements are grouped for ease of use. Intended for ICP-MS and ICP-OES, they can be used individually or in any combination upon dilution into 1% HNO₃. Custom ICP-MS/OES calibration standards are available upon request.

71-Element Group

ICP-MS Precious Metals Standard			
IV-ICPMS-71C C		Matrix: HCl	
IV-ICPMS-71C-125ML		Volume: 125 mL	
IV-ICPMS-71C-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	10	Pt	10
Ir	10	Re	10
Os	10	Rh	10
Pd	10	Ru	10

Contains osmium (Os); avoid dilutions with oxidizing acids such as concentrated HNO₃. For dilutions including Ag please see Silver Chemical Stability article for more information about Ag stability in HCl.

ICP-MS Internal Standard			
IV-ICPMS-71D C		Matrix: HNO ₃	
IV-ICPMS-71D-125ML		Volume: 125 mL	
IV-ICPMS-71D-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	10	Sc	10
In	10	Tb	10
⁶ Li	10	Y	10

Covers mass range from 6-Li to 209-Bi. Certified reference material that may also be used for calibration. Can be diluted with other multi-element standards to working concentrations.

Lithium ICP-MS Standard	
MSLI-10PPM C	Matrix: HNO ₃
MSLI-10PPM-125ML	Volume: 125 mL
Analyte	µg/mL
Li	10

Mercury ICP-MS Standard	
MSHG-10PPM C	Matrix: HCl
MSHG-10PPM-125ML	Volume: 125 mL
MSHG-10PPM-500ML	Volume: 500 mL
Analyte	µg/mL
Hg	10

Tellurium ICP-MS Standard	
MSTEN-100PPM C	Matrix: HNO ₃
MSTEN-100PPM-125ML	Volume: 125 mL
Analyte	µg/mL
Te	100

C Common Multi-Element Standard

AGI Tuning Solution			
AGI-TS-1 		Matrix: HNO ₃	
AGI-TS-1-125ML		Volume: 125 mL	
AGI-TS-1-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ce	10	Tl	10
Co	10	Y	10
Li	10		

ICP-MS stock tuning solution designed for dilution to working concentrations. Covers mass range from Li to Tl. Certified reference material that may also be used for calibration. Agilent P/N 5188-6564.

CIROS Tuning Solution			
CIROS-OES-TS 		Matrix: HCl / HNO ₃	
CIROS-OES-TS-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Fe	10	P	10
K	10	S	50
La	10	Sc	10
Mg	5	Ti	10
Mn	5		

For reprofiling optics of Spectro Ciros ICP-OES.

GENESIS Calibration Standard			
GENESIS-ICAL 		Matrix: HNO ₃ / HCl / HF	
GENESIS-ICAL-125ML		Volume: 125 mL	
GENESIS-ICAL-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Be	2	Na	5
Ca	1	Ni	10
Ce	10	P	10
Cu	10	S	50
Eu	10	Sc	5
Fe	10	Si	10
In	10	Sr	2
K	10	Ti	10
Li	2	V	10
Mn	5	Y	10
Mo	5	Zr	10

For reprofiling optics of Spectro Ciros ICP-OES.

Trace Metals in Water- SRM1643			
IV-STOCK-1643 		Matrix: HNO ₃	
IV-STOCK-1643-125ML		Volume: 125 mL	
IV-STOCK-1643-500ML		Volume: 500 mL	
Analyte	µg/L*	Analyte	µg/L*
Ag	1	Mg	8,000
Al	142	Mn	39
As	60	Mo	121
B	158	Na	21,000
Ba	544	Ni	62
Be	14	Pb	20
Bi	14	Rb	14
Ca	32,000	Re	113
Cd	7	Sb	58
Co	27	Se	12
Cr	20	Sr	323
Cu	23	Te	1
Fe	98	Tl	7
K	2,000	V	38
Li	17	Zn	79

*Parts per billion

For quality control and method evaluation of fresh water trace element analyses. Ready to use without dilution.

Instrument Check Standard			
PE-CHK-1 		Matrix: HNO ₃ / HF	
PE-CHK-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Mn	10
Al	10	Ni	10
As	10	Pb	10
Ba	10	Sb	10
Be	10	Se	10
Cd	10	Tl	10
Co	10	V	10
Cr	10	Zn	10
Cu	10		

For daily instrument calibration.

 Agilent/Varian

 NIST

 Perkin Elmer

 Spectro

MULTI-ELEMENT STANDARDS

Multi-Element Standards

Identical or near identical formulations

Tuning Solution			
PE-TS-1 		Matrix: HNO ₃	
PE-TS-1-125ML		Volume: 125 mL	
PE-TS-1-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba	10	Mg	10
Be	10	Pb	10
Ce	10	Rh	10
Co	10	Tl	10
In	10	U	10
Li	10	Y	10

For instrument set-up and calibration. Covers mass range from Li to U (isotopically depleted).

ICP-MS Tuning Solution – Tune B iCAP			
THERMO-4AREV 		Matrix: HNO ₃ /HCl	
THERMO-4AREV-500ML		Volume: 500 mL	
THERMO-4AREV-1L		Volume: 1 L	
Analyte	µg/L*	Analyte	µg/L*
Ba	1	In	1
Bi	1	Li	1
Ce	1	U	1
Co	1		

*Parts per billion

Tuning solution for Thermo iCAP Q ICP-MS. Equivalent to Thermo P/N 1323770.

ICP-MS Tuning Solution – iCAP Q			
THERMO-5A 		Matrix: HNO ₃	
THERMO-5A-125ML		Volume: 125 mL	
THERMO-5A-250ML		Volume: 250 mL	
Analyte	µg/L*	Analyte	µg/L*
Ag	6	Mg	10
Al	10	Mn	6
Ba	4	Ni	15
Be	35	Rh	3
Bi	3	Sc	8
Ce	3	Sr	5
Co	8	Ta	3
Cs	3	Tb	3
Cu	15	Tl	4
Ga	10	U	3
Ho	3	Y	3
In	3	Zn	20
Li	8		

*Parts per billion

Calibration standard for Thermo iCAP Q ICP-MS. Equivalent to Thermo P/N 1323760.

-  Common Multi-Element Standard
-  Perkin Elmer
-  Thermo Scientific

Tune F-X-Series Tuning Solution			
TUNE F-X-SERIES 		Matrix: HNO ₃ /HF	
TUNE F-X-SERIES-125ML		Volume: 125 mL	
Analyte	ng/mL*	Analyte	ng/mL*
Ag	40	Na	40
Al	50	Nb	20
As	250	Nd	45
B	200	Ni	150
Ba	50	P	1000
Be	500	Pb	10
Bi	5	Pd	100
Ca	1000	Pr	10
Cd	100	Rb	30
Ce	10	Re	15
Co	35	Sb	40
Cr	40	Sc	30
Cs	15	Se	1250
Cu	150	Si	1000
Dy	25	Sm	45
Er	15	Sn	45
Eu	10	Sr	20
Fe	20	Ta	5
Ga	45	Tb	5
Gd	45	Te	500
Ge	150	Th	5
Hf	15	Ti	500
Ho	5	Tl	10
In	10	Tm	5
K	35	U	5
La	10	V	40
Li	100	W	25
Lu	5	Y	15
Mg	50	Yb	25
Mn	20	Zn	150
Mo	100	Zr	35

*Parts per billion

For detector cross-calibration on Thermo X-Series ICP-MS.

 Agilent/Varian Thermo Scientific

Calibration Standard			
VAR-CAL-1 		Matrix: HNO ₃ / HF	
VAR-CAL-1-125ML		Volume: 125 mL	
VAR-CAL-1-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Mo	100	Sn	100
Sb	100	Ti	100

General ICP-OES calibration standard. Designed to be mixed with VAR-CAL-2 at working concentrations.

Calibration Standard			
VAR-CAL-2 		Matrix: HNO ₃	
VAR-CAL-2-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	Mn	100
Al	100	Ni	100
As	100	Pb	100
Ba	100	Se	100
Be	100	Th	100
Cd	100	Tl	100
Co	100	U	100
Cr	100	V	100
Cu	100	Zn	100

General ICP-OES calibration standard. Designed to be mixed with VAR-CAL-1 at working concentrations.

Calibration Standard			
VAR-CAL-7 		Matrix: HNO ₃ /HF	
VAR-CAL-7-125ML		Volume: 125 mL	
VAR-CAL-7-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	5	Mn	5
As	5	Mo	5
Ba	5	Ni	5
Cd	5	Pb	5
Co	5	Se	5
Cr	5	Sr	5
Cu	5	Zn	5
K	50		

ICP-OES calibration standard.

MULTI-ELEMENT STANDARDS

Multi-Element Standards

Identical or near identical formulations

ICP Internal Standard			
VAR-IS-1		Matrix: HNO ₃	
VAR-IS-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	100	Sc	100
In	100	Tb	100
⁶ Li	100	Y	100

For use as ICP-MS multi-element internal standard. Covers mass range from 6-Li to 209-Bi

 Agilent/Varian

Tuning Solution			
VAR-TS-MS		Matrix: HNO ₃	
VAR-TS-MS-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba	10	Mg	10
Be	10	Pb	10
Ce	10	Th	10
Co	10	Tl	10
In	10		

For use as ICP-MS tuning solution. Covers mass range from 9-Be to 232-Th. Certified reference material that may also be used for calibration.

HIGH-PURITY IONIZATION BUFFERS

Ionization Buffers

Ionization buffers are 99.999+% pure. They are analyzed using both axial-view ICP-OES and ICP-MS for 70+ impurities. Custom ionization buffers are available upon request.

1% Cesium Ionization Buffer	
CSN-ISB	Matrix: HNO ₃
CSN-ISB-500ML	Volume: 500 mL
Analyte	µg/mL
Cs	10,000
High Purity buffer; ideal for Axial View ICP-OES	

For stabilizing the degree of ionization in flame AA and ICP-OES analysis.

2% Lithium Ionization Buffer	
LINB2	Matrix: HNO ₃
LINB2-125ML	Volume: 125 mL
Analyte	µg/mL
Li	20,000

For stabilizing the degree of ionization in flame AA and ICP-OES analysis. Not to be used as a calibration standard, for analytical reagent use only.

5% Cesium Ionization Buffer	
CSN-ISB5	Matrix: HNO ₃
CSN-ISB5-500ML	Volume: 500 mL
Analyte	µg/mL
Cs	50,000

For stabilizing the degree of ionization in flame AA and ICP-OES analysis.

Manufactured with high-purity starting materials and reagents. Products can be used for calibration of analytical instruments, validation of analytical methods, or for other applications deemed fit for purpose by the end-user. Can be diluted with other multi-element standards to working concentrations. Certificate of Analysis includes lot specific trace metal impurity analysis.

For the pharmaceutical industry, Inorganic Ventures has developed CRMs to comply with the United States Pharmacopeia (USP) general chapters on elemental impurity USP <232> limits and USP <233> procedures.

These methods are for testing inorganic impurities in pharmaceutical products by ICP. The International Conference on Harmonization (ICH) Working Group on Elemental Impurities is in the process of developing a harmonized approach for controlling these impurities as well.

USP <232> Precious Metals Elemental Impurities			
IV-STOCK-38 		Matrix: HCl	
IV-STOCK-38-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ir	100	Pt	100
Os	100	Rh	100
Pd	100	Ru	100

USP <232> Oral Elemental Impurities			
IV-STOCK-40 		Matrix: HNO ₃	
IV-STOCK-40-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	1.5	Mo	100
Cd	25	Ni	500
Cu	1000	Pb	5
Hg	15	V	100

USP <232> Parenteral Elemental Impurities			
IV-STOCK-41 		Matrix: HNO ₃	
IV-STOCK-41-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	1.5	Mo	10
Cd	2.5	Ni	50
Cu	100	Pb	5
Hg	1.5	V	10

USP <232> Drug Substance and Excipients			
IV-STOCK-60 		Matrix: HCl	
IV-STOCK-60-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	1.5	Os	10
Cd	0.5	Pb	0.5
Cr	1100	Pd	10
Cu	300	Pt	10
Hg	3	Rh	10
Ir	10	Ru	10
Mo	300	V	10
Ni	20		

USP <232> / ICH Q3D Class 1 Oral Elemental Impurities			
IV-STOCK-65  		Matrix: HNO ₃	
IV-STOCK-65-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	15	Hg	30
Cd	5	Pb	5

USP <232> / ICH Q3D Class 2A Oral Elemental Impurities			
IV-STOCK-66 		Matrix: HNO ₃	
IV-STOCK-66-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Co	50	V	100
Ni	200		

 Agilent/Varian

 USP Method <232>

MULTI-ELEMENT STANDARDS

USP Method <232> – Elemental Impurities Compliance Standards

USP <232> / ICH Q3D Class 2B Oral Elemental Impurities			
IV-STOCK-67 U		Matrix: HCl	
IV-STOCK-67-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	100	Rh	100
Ir	100	Ru	100
Os	100	Se	150
Pd	100	Tl	8
Pt	100		

USP <232> / ICH Q3D Class 2B Oral Elemental Impurities			
IV-STOCK-68 U		Matrix: HNO ₃	
IV-STOCK-68-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag*	150		

* Silver has been separated from the other Class 2B elements due to long-term stability concerns. However, IV-STOCK-68 can be combined with IV-STOCK-67 at working levels. Contact Technical Support or visit our Technical Forum for more information regarding Ag in HCl matrices.

USP <232> / ICH Q3D Class 3 Oral Elemental Impurities			
IV-STOCK-69 U		Matrix: HNO ₃ /tr HF	
IV-STOCK-69-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba	140	Mo	300
Cr	1100	Sb	120
Cu	300	Sn	600
Li	55		

USP <232> / ICH Q3D Oral Elemental Impurities			
IV-STOCK-70 U		Matrix: HCl	
IV-STOCK-70-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	15	Ni	20
As	1.5	Os	10
Au	10	Pb	0.5
Ba	140	Pd	10
Cd	0.5	Pt	10
Co	5	Rh	10
Cr	1100	Ru	10
Cu	300	Sb	120
Hg	3	Se	15
Ir	10	Sn	600
Li	55	Tl	0.8
Mo	300	V	10

USP <232> / ICH Q3D Oral Elemental Impurities			
IV-STOCK-78 U		Matrix: HCl	
IV-STOCK-78-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	15	Ni	20
As	1.5	Os	10
Au	30	Pb	0.5
Ba	140	Pd	10
Cd	0.5	Pt	10
Co	5	Rh	10
Cr	1100	Ru	10
Cu	300	Sb	120
Hg	3	Se	15
Ir	10	Sn	600
Li	55	Tl	0.8
Mo	300	V	10

USP <232> / ICH Q3D Class 2B Oral Elemental Impurities			
IV-STOCK-79 U		Matrix: HCl	
IV-STOCK-79-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Au	300	Rh	100
Ir	100	Ru	100
Os	100	Se	150
Pd	100	Tl	8
Pt	100		

U USP Method <232>

Bringing Confidence to the Cannabis Industry

In an industry where regulations and testing requirements vary by state, put your trust in Inorganic Ventures. We aim to squash inter-laboratory variations in the Cannabis industry and make it easier for your lab to get accurate results time and time again.

We offer single and multi-element standards to ensure your final product is free of dangerous heavy metals that may arise from soil contamination in agricultural production or manufacturing processes.

- Extensive documentation: Certificate of Analysis (CoA) and Safety Data Sheet (SDS)
- TCT packaging prevents transpiration, guarantees a 5-year shelf life, and allows for storage outside of normal lab conditions.
- 100% satisfaction guarantee



10 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Arsenic, As	HNO ₃	125 mL	MSAS-10PPM-125ML
Cadmium, Cd	HNO ₃	125 mL	MSCD-10PPM-125ML
Lead, Pb	HNO ₃	125 mL	MSPB-10PPM-125ML
Mercury, Hg	HCl	125 mL 500 mL	MSHG-10PPM-125ML MSHG-10PPM-500ML
Mercury, Hg	HNO ₃	125 mL 500 mL	MSHGN-10PPM-125ML MSHGN-10PPM-500ML

100 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Arsenic, As	HNO ₃	125 mL	MSAS-100PPM-125ML
Cadmium, Cd	HNO ₃	125 mL	MSCD-100PPM-125ML
Lead, Pb	HNO ₃	125 mL 500 mL	MSPB-100PPM-125ML MSPB-100PPM-500ML
Mercury, Hg	HCl	125 mL	MSHG-100PPM-125ML
Mercury, Hg	HNO ₃	125 mL	MSHGN-100PPM-125ML

MULTI-ELEMENT STANDARDS

Heavy Metals – Single Element Standards

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Arsenic, As	HNO ₃	30 mL	CGAS1-30ML
		125 mL	CGAS1-125ML
		500 mL	CGAS1-500ML
Arsenic ⁺³ , As ⁺³	HCl / NaOH / NaHCO ₃	30 mL	CGAS(3)1-30ML
		125 mL	CGAS(3)1-125ML
		500 mL	CGAS(3)1-500ML
Arsenic ⁺⁵ , As ⁺⁵	H ₂ O	30 mL	CGAS(5)1-30ML
		125 mL	CGAS(5)1-125ML
		500 mL	CGAS(5)1-500ML
Cadmium, Cd	HNO ₃	30 mL	CGCD1-30ML
		125 mL	CGCD1-125ML
		500 mL	CGCD1-500ML
Lead, Pb	HNO ₃	30 mL	CGPB1-30ML
		125 mL	CGPB1-125ML
		500 mL	CGPB1-500ML
Mercury, Hg	HNO ₃	30 mL	CGHG1-30ML
		125 mL	CGHG1-125ML
		500 mL	CGHG1-500ML

10,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Arsenic, As	HNO ₃	30 mL	CGAS10-30ML
		125 mL	CGAS10-125ML
		500 mL	CGAS10-500ML
Cadmium, Cd	HNO ₃	125 mL	CGCD10-125ML
		500 mL	CGCD10-500ML
Lead, Pb	HNO ₃	30 mL	CGPB10-30ML
		125 mL	CGPB10-125ML
		500 mL	CGPB10-500ML
Mercury, Hg	HNO ₃	125 mL	CGHG10-125ML
		500 mL	CGHG10-500ML



Does your state require testing for elements in addition to the Big 4? Request a custom quote! As the leading manufacturer of custom inorganic standards, we've produced tens of thousands of unique custom blends for laboratories worldwide.

While our USP <232>/ICH Q3D stock products were formulated for the pharmaceutical industry, these products can be used to test heavy metals in cannabis. The Big 4 (Arsenic, Mercury, Cadmium, and Lead) are common analytes for cannabis testing.

USP <232> / ICH Q3D Class 1 Oral Elemental Impurities			
IV-STOCK-65-125ML		Volume: 125 mL	Matrix: HNO ₃
Analyte	µg/mL	Analyte	µg/mL
As	15	Hg	30
Cd	5	Pb	5


The following are custom products. They are available to order, but not are not stock items.

Custom Heavy Metal Standard			
IV-6239		Matrix: HNO ₃	
Analyte	µg/mL	Analyte	µg/mL
As	15	Hg	3
Cd	2	Pb	5

Custom Heavy Metal Standard			
IV-48592		Matrix: HNO ₃	
Analyte	µg/mL	Analyte	µg/mL
As	2	Hg	1
Cd	2	Pb	5

Don't see what you need?

Contact us with the solution part number and instrument manufacturer you're seeking, and we'll check our extensive library of solutions.



Over the years, we've developed a unique line of EPA standards. If you do not see what you are looking for, please contact us with an EPA custom request and we will get you competitive pricing guaranteed.



Contents

ILM03.0.....	59
ILM04.0.....	61
ILM05.2 & ILM05.3.....	63
Method 200.7.....	66
Method 200.8.....	76
Method 1311.....	79
Method 6020.....	80
Need a Custom CRM?.....	16

- ✓ **Five-year shelf life**
- ✓ **Traceable to NIST SRMs**
- ✓ **Produced under ISO 9001**
- ✓ **Produced under ISO 17025**
- ✓ **Produced under ISO 17034**
- ✓ **Assayed by validated Wet Chemical procedures**
- ✓ **Assayed by validated instrument procedures**

Standards for ILM03.0 are designed for use with ICP-OES. Custom EPA standards are available upon request.

Calibration Standard			
CLPP-CAL-1		Matrix: HNO ₃ Dilution 1:100	
CLPP-CAL-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	250	Fe	1,000
Al	2,000	K	5,000
Ba	2,000	Mg	5,000
Be	50	Mn	500
Ca	5,000	Na	5,000
Co	500	Ni	500
Cr	200	V	500
Cu	250	Zn	500

For use as ICP calibration standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or CLPP-SPK-2 for suitable Sb reference standard.

Calibration Standard			
CLPP-CAL-3		Matrix: HNO ₃ Dilution 1:100	
CLPP-CAL-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	1,000	Se	1,000
Cd	500	Tl	1,000
Pb	1,000		

For use as ICP calibration standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or CLPP-SPK-2 for suitable Sb reference standard.

Calibration Standard	
CGSB1	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100
CGSB1-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	1,000

CICV Standards – Continuing and Initial Calibration Verification

CICV Standard†			
QCP-CICV-1		Matrix: HNO ₃ Dilution 1:100 or 1:500	
QCP-CICV-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	125	Fe	500
Al	1,000	K	2,500
Ba	1,000	Mg	2,500
Be	25	Mn	250
Ca	2,500	Na	2,500
Co	250	Ni	250
Cr	100	V	250
Cu	125	Zn	250

For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or QCP-CICV-2 for suitable Sb reference standard.

CICV Standard†	
QCP-CICV-2	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100 or 1:500
QCP-CICV-2-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	500

For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. For analyses requiring antimony (Sb). Designed to be diluted to working concentrations with QCP-CICV-1 and/or QCP-CICV-3.

CICV Standard†			
QCP-CICV-3		Matrix: HNO ₃ Dilution 1:100 or 1:500	
QCP-CICV-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	500	Se	500
Cd	250	Tl	500
Pb	500		


For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or QCP-CICV-2 for suitable Sb reference standard.

†Manufactured from in-house Second Source concentrates, whenever possible.

CRDL Standards – Contract Required Detection Limit

We can create any CRDL standard to best fit your needs.

Custom solutions are our specialty.



Soil & Water Spike Standards

Spike Standard*			
CLPP-SPK-1		Matrix: HNO ₃ Dilution 1:1,000	
CLPP-SPK-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	50	Cu	250
Al	2,000	Fe	1,000
Ba	2,000	Mn	500
Be	50	Ni	500
Co	500	V	500
Cr	200	Zn	500

For use as ICP-OES soil or water spike standard in EPA Contract Laboratory Program (CLP) methods.

Spike Standard*	
CLPP-SPK-2	Matrix: HNO ₃ /Tartaric Acid Dilution 1:1,000
CLPP-SPK-2-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	500

For use as Sb spike standard in EPA Contract Laboratory Program (CLP) methods.

*Instructions included.

Interference Check Standards

For use as ICP-OES soil or water spike standard in EPA Contract Laboratory Program (CLP) methods.

Interference Check Standard			
CLPP-ICS-A  		Matrix: HNO ₃ Dilution 1:10	
CLPP-ICS-A-125ML		Volume: 125 mL	
CLPP-ICS-A-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	5,000	Fe	2,000
Ca	5,000	Mg	5,000

Interference Check Standard			
CLPP-ICS-B		Matrix: HNO ₃ Dilution 1:100	
CLPP-ICS-B-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	Cu	50
Ba	50	Mn	50
Be	50	Ni	100
Cd	100	Pb	100
Co	50	V	50
Cr	50	Zn	100

 Agilent/Varian

 Perkin Elmer

Standards for ILM04.0 are designed for use with ICP-OES. Custom EPA standards are available upon request.

Calibration Standards

Calibration Standard			
CLPP-CAL-1		Matrix: HNO ₃ Dilution 1:100	
CLPP-CAL-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	250	Fe	1,000
Al	2,000	K	5,000
Ba	2,000	Mg	5,000
Be	50	Mn	500
Ca	5,000	Na	5,000
Co	500	Ni	500
Cr	200	V	500
Cu	250	Zn	500

For use as ICP calibration standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or CLPP-SPK-2 for suitable Sb reference standard.

Calibration Standard			
CLPP-CAL-3		Matrix: HNO ₃ Dilution 1:100	
CLPP-CAL-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	1,000	Se	1,000
Cd	500	Tl	1,000
Pb	1,000		

For use as ICP calibration standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or CLPP-SPK-2 for suitable Sb reference standard.

Calibration Standard	
CGSB1	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100
CGSB1-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	1,000

†Manufactured from in-house Second Source concentrates, whenever possible.

CICV Standards – Continuing and Initial Calibration Verification

CICV Standard†			
QCP-CICV-1		Matrix: HNO ₃ Dilution 1:100 or 1:500	
QCP-CICV-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	125	Fe	500
Al	1,000	K	2,500
Ba	1,000	Mg	2,500
Be	25	Mn	250
Ca	2,500	Na	2,500
Co	250	Ni	250
Cr	100	V	250
Cu	125	Zn	250

For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or QCP-CICV-2 for suitable Sb reference standard.

CICV Standard†	
QCP-CICV-2	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100 or 1:500
QCP-CICV-2-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	500

For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. For analyses requiring antimony (Sb). Designed to be diluted to working concentrations with QCP-CICV-1 and/or QCP-CICV-3.

CICV Standard†			
QCP-CICV-3		Matrix: HNO ₃ Dilution 1:100 or 1:500	
QCP-CICV-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	500	Se	500
Cd	250	Tl	500
Pb	500		


For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or QCP-CICV-2 for suitable Sb reference standard.

ILMO4.0

CRDL Standards – Contract Required Detection Limit

We can create any CRDL standard to best fit your needs.

Custom solutions are our specialty.



Soil & Water Spike Standards

Spike Standard*			
CLPP-SPK-1		Matrix: HNO ₃ Dilution 1:1,000	
CLPP-SPK-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	50	Cu	250
Al	2,000	Fe	1,000
Ba	2,000	Mn	500
Be	50	Ni	500
Co	500	V	500
Cr	200	Zn	500

*Instructions included.

For use as ICP-OES soil or water spike standard in EPA Contract Laboratory Program (CLP) methods.



Don't see what you need?

Contact us with the solution part number and instrument manufacturer you're seeking, and we'll check our extensive library of solutions.

Interference Check Standards

Interference Check Standard A			
CLPP-ICS-A		Matrix: HNO ₃ Dilution 1:10	
CLPP-ICS-A-125ML		Volume: 125 mL	
CLPP-ICS-A-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	5,000	Fe	2,000
Ca	5,000	Mg	5,000

For use as ICP-OES interference check standard in EPA Contract Laboratory Program (CLP) methods.

Interference Check Standard B4			
CLPP-ICS-B4		Matrix: HNO ₃ Dilution 1:100	
CLPP-ICS-B4-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	20	Mn	50
As	10	Ni	100
Ba	50	Pb	5
Be	50	Sb	60
Cd	100	Se	5
Co	50	Tl	10
Cr	50	V	50
Cu	50	Zn	100

For use as ICP-OES interference check standard in EPA Contract Laboratory Program (CLP) methods.

Agilent/Varian

Perkin Elmer

See individual products for recommended instrumentation and revision. Custom EPA standards are available upon request.

Calibration Standards

Calibration Standard	
CGSB1	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100
CGSB1-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	1,000

Calibration Standard			
CLPP-CAL-1		Matrix: HNO ₃ Dilution 1:100	
CLPP-CAL-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	250	Fe	1,000
Al	2,000	K	5,000
Ba	2,000	Mg	5,000
Be	50	Mn	500
Ca	5,000	Na	5,000
Co	500	Ni	500
Cr	200	V	500
Cu	250	Zn	500

For use as ICP calibration standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or CLPP-SPK-2 for suitable Sb reference standard.

Calibration Standard			
CLPP-CAL-3		Matrix: HNO ₃ Dilution 1:100	
CLPP-CAL-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	1,000	Se	1,000
Cd	500	Tl	1,000
Pb	1,000		

For use as ICP calibration standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or CLPP-SPK-2 for suitable Sb reference standard.

CICV Standards – Continuing and Initial Calibration Verification

CICV Standard†			
QCP-CICV-1		Matrix: HNO ₃ Dilution 1:100 or 1:500	
QCP-CICV-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	125	Fe	500
Al	1,000	K	2,500
Ba	1,000	Mg	2,500
Be	25	Mn	250
Ca	2,500	Na	2,500
Co	250	Ni	250
Cr	100	V	250
Cu	125	Zn	250

For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or QCP-CICV-2 for suitable Sb reference standard.

CICV Standard†	
QCP-CICV-2	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100 or 1:500
QCP-CICV-2-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	500

For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. For analyses requiring antimony (Sb). Designed to be diluted to working concentrations with QCP-CICV-1 and/or QCP-CICV-3.

CICV Standard†			
QCP-CICV-3		Matrix: HNO ₃ Dilution 1:100 or 1:500	
QCP-CICV-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	500	Se	500
Cd	250	Tl	500
Pb	500		

For use as initial/continuing calibration verification standard in EPA Contract Laboratory Program (CLP) methods. Does not contain antimony (Sb). Please see CGSB1 or QCP-CICV-2 for suitable Sb reference standard.

†Manufactured from in-house Second Source concentrates, whenever possible.

Interference Check Standards

Interference Check Standard A			
CLPP-ICS-A 		Matrix: HNO ₃ Dilution 1:10	
CLPP-ICS-A-125ML		Volume: 125 mL	
CLPP-ICS-A-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	5,000	Fe	2,000
Ca	5,000	Mg	5,000

For use as ICP-OES interference check standard in EPA Contract Laboratory Program (CLP) methods.

Interference Check Standard B4			
CLPP-ICS-B4 		Matrix: HNO ₃ Dilution 1:100	
CLPP-ICS-B4-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	20	Mn	50
As	10	Ni	100
Ba	50	Pb	5
Be	50	Sb	60
Cd	100	Se	5
Co	50	Tl	10
Cr	50	V	50
Cu	50	Zn	100

For use as ICP-OES interference check standard in EPA Contract Laboratory Program (CLP) methods.



Now, YOU control the expiration date.

With TCT, concerns about shipping and storage conditions are eliminated, as transpiration is no longer an issue. This saves you money and simplifies research by removing the need to constantly inventory and restock CRMs.

Find out more at inorganicventures.com/tct.

 Agilent/Varian

 Perkin Elmer

Soil & Water Spike Standards

Spike Standard			
CLP-MS-SPK		Matrix: HNO ₃ Dilution 1:100	
CLP-MS-SPK-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	5	Mn	50
Al	200	Ni	50
As	4	Pb	2
Ba	200	Sb	10
Be	5	Se	1
Cd	5	Tl	5
Co	50	V	50
Cr	20	Zn	50
Cu	25		

For use as ICP soil or water spike standard in EPA Contract Laboratory Program (CLP) methods

Spike Standard			
CLPP-SPK-1		Matrix: HNO ₃ Dilution 1:1,000	
CLPP-SPK-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	50	Cu	250
Al	2,000	Fe	1,000
Ba	2,000	Mn	500
Be	50	Ni	500
Co	500	V	500
Cr	200	Zn	500

For use as ICP-OES soil or water spike standard in EPA Contract Laboratory Program (CLP) methods

Internal Standards & Tuning Solutions

Internal Standard			
6020ISS		Matrix: HNO ₃ Dilution 1:100	
6020ISS-125ML		Volume: 125 mL	
6020ISS-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	10	Rh	10
Ho	10	Sc	10
In	10	Tb	10
⁶ Li	10	Y	10

Internal standard for ICP-MS in all versions of EPA Method 6020.

Tuning Solution			
6020TS		Matrix: HNO ₃ Dilution 1:100	
6020TS-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Co	10	Li	10
In	10	Tl	10

For use as general tuning solution suitable for numerous ICP-MS designs and models. Covers mass range from Li to Tl. Certified reference material that may also be used for calibration.

Tuning Solution			
2008TS		Matrix: HNO ₃ Dilution 1:100 to 1:1,000	
2008TS-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Be	10	Mg	10
Co	10	Pb	10
In	10		

For use as ICP-MS tuning solution in EPA Method 200.8.

 Agilent/Varian

Blank & Rinse Solutions

Blank & Rinse solutions are prepared using double-distilled reagents and 18 megohm (MΩ) deionized water. They come packaged in ultra-clean LDPE bottles and are ready to use. Custom solutions are available upon request.

2% (v/v) Nitric Acid Rinse	
CLP-MS-RINSE	Matrix: HNO ₃
CLP-MS-RINSE-125ML	Volume: 125 mL
CLP-MS-RINSE-500ML	Volume: 500 mL

See pg. 108 for more Blank & Rinse Solution options.

For use as ultra-pure nitric acid ICP rinse or blank solution. Suitable for EPA Contract Laboratory Program (CLP) methods. Manufactured using ultra-high purity HNO₃, >18 MΩ deionized water, and packaged in specially cleaned LDPE bottles. Certificate includes trace metal impurity values representative of typical analyses.

200.7 Calibration

Standards for Method 200.7 are designed for use with ICP-OES. Custom EPA standards are available upon request. Standards are designed for Method 200.7, Method 3120, Method 6010A Rev. 1 and Method 200.7 CLP-M.

Calibration Standard		
CLPP-SPK-2	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100	
CLPP-SPK-2-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
Sb	500	206.833

For use as Sb spike standard in EPA Contract Laboratory Program (CLP) methods.

Calibration Standard		
WW-CAL-2	Matrix: HNO ₃ / HF Dilution 1:100	
WW-CAL-2-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
K	2,000	766.491
Li	500	670.784
Mo	1,000	203.844
Na	1,000	588.995
Ti	1,000	334.941

For use as ICP-OES calibration standard II in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Calibration Standard		
WW-CAL-1A	Matrix: HNO ₃ Dilution 1:100	
WW-CAL-1A-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
Ag	50	328.068
As	1,000	193.759
B	100	249.678
Ba	100	493.409
Ca	1,000	315.887
Cd	200	226.502
Cu	200	324.754
Mn	200	257.610
Se	500	196.090
Sr*	100	421.552

For use as ICP-OES calibration standard I in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions

*NOTE: Sr does not exhibit spectral interference problems with any of the EPA Method 200.7 analytes.

Calibration Standard		
WW-CAL-3	Matrix: HNO ₃ Dilution 1:100	
WW-CAL-3-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
Ce	200	413.765
Co	200	228.616
P	1,000	214.914
V	200	292.402

For use as ICP-OES calibration standard III in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

200.7 Calibration

Calibration Standard		
WW-CAL-4A		Matrix: HNO ₃ Dilution 1:100
WW-CAL-4A-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Al	1,000	308.215
Cr	500	205.552
Hg	200	194.227
Zn	500	213.856

For use as ICP-OES calibration standard IV (Part A) in EPA Method 200.7. Designed to be mixed with WW-CAL-4B at working concentrations. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Calibration Standard		
WW-CAL-4B		Matrix: HNO ₃ / HF Dilution 1:100
WW-CAL-4B-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
SiO₂	1,000	251.611
Sn	400	189.980

For use as ICP-OES calibration standard IV (Part B) in EPA Method 200.7. Designed to be mixed with WW-CAL-4A at working concentrations. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Calibration Standard		
WW-CAL-5		Matrix: HNO ₃ Dilution 1:100
WW-CAL-5-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Be	100	313.042
Fe	1,000	259.940
Mg	1,000	279.079
Ni	200	231.604
Pb	1,000	220.353
Tl	500	190.864

For use as ICP-OES calibration standard V in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

 Perkin Elmer

200.7 Interference Checks

Interference Check Standard	
CGSB1	Matrix: HNO ₃ /Tartaric Acid Dilution 1:100
CGSB1-125ML	Volume: 125 mL
Analyte	µg/mL
Sb	1,000

Interference Check Standard			
2007ICS-1		Matrix: HNO ₃ / HF Dilution 1:100	
2007ICS-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
B	500	Si	230
Mo	300	Ti	1,000

For use as ICP-OES interference check standard in EPA Method 200.7.

Interference Check Standard			
2007ICS-3		Matrix: HNO ₃ Dilution 1:100	
2007ICS-3-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	300	K	20,000
As	1,000	Mn	200
Ba	300	Ni	300
Be	100	Pb	1,000
Cd	300	Se	500
Co	300	Tl	1,000
Cr	300	V	300
Cu	300	Zn	300

For use as ICP-OES interference check standard in EPA Method 200.7.

Interference Check Standard			
2007ICS-4		Matrix: HNO ₃ Dilution 1:50	
2007ICS-4-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	3,000	Mg	7,500
Ca	15,000	Na	2,500
Fe	12,500		

For use as ICP-OES interference check standard in EPA Method 200.7.

200.7 Quality Controls

Quality Control Standard†		
QCP-QCS-1		Matrix: HNO ₃ Dilution 1:100
QCP-QCS-1-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Ag	25	328.068
Al	100	308.215
As	200	193.759
B	100	249.678
Ba	100	493.409
Be	100	313.042
Ca	100	315.887
Cd	100	226.502
Ce	100	413.765
Co	100	228.616
Cr	100	205.552
Cu	100	324.754
Fe	100	259.940
Hg	200	194.227
K	500	766.491
Li	100	670.784
Mg	100	279.079
Mn	100	257.610
Na	100	588.995
Ni	100	231.604
P	500	214.914
Pb	200	220.353
Se	100	196.090
Sr	100	421.552
Tl	500	190.864
V	100	292.402
Zn	100	213.856

For use as ICP-OES QC standard in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Quality Control Standard†		
QCP-QCS-2		Matrix: HNO ₃ / HF Dilution 1:100
QCP-QCS-2-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Mo	100	203.844
Sb	200	206.833
SiO ₂	500	251.611
Sn	500	189.980
Ti	100	334.941

For use as ICP-OES QC standard in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Quality Control Standard†			
IV-7		Matrix: HNO ₃ / HF Dilution 1:100	
IV-7-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	K	1,000
Al	100	Na	100
B	100	Si	50
Ba	100		

For use as a certified reference standard in ICP applications. Can be diluted with other standards to working concentrations.

JY HORIBA Jobin Yvon



Don't see what you need?

Contact us with the solution part number and instrument manufacturer you're seeking, and we can check our extensive library of solutions.

†Manufactured from in-house Second Source concentrates, whenever possible.

200.7 Quality Controls

For use as a certified reference standard in ICP applications. Can be diluted with other standards to working concentrations.

Quality Control Standard†			
IV-19		Matrix: HNO ₃ / HF Dilution 1:100	
IV-19-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	100	Mo	100
Be	100	Ni	100
Ca	100	Pb	100
Cd	100	Sb	100
Co	100	Se	100
Cr	100	Ti	100
Cu	100	Tl	100
Fe	100	V	100
Mg	100	Zn	100
Mn	100		

Quality Control Standard†			
IV-21		Matrix: HNO ₃ / HF Dilution 1:100	
IV-21-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	100	Mo	100
Be	100	Ni	100
Ca	100	Pb	100
Cd	100	Sb	100
Co	100	Se	100
Cr	100	Sr	100
Cu	100	Ti	100
Fe	100	Tl	100
Li	100	V	100
Mg	100	Zn	100
Mn	100		

JY HORIBA Jobin Yvon

M Merck/MilliporeSigma

PE Perkin Elmer

†Manufactured from in-house Second Source concentrates, whenever possible.

Quality Control Standard†			
IV-26		Matrix: HNO ₃ / HF Dilution 1:100	
IV-26-125ML IV-26-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	Mg	100
Al	100	Mn	100
As	100	Mo	100
B	100	Na	100
Ba	100	Ni	100
Be	100	Pb	100
Ca	100	Sb	100
Cd	100	Se	100
Co	100	Si	50
Cr	100	Ti	100
Cu	100	Tl	100
Fe	100	V	100
K	1,000	Zn	100

Quality Control Standard†			
IV-28		Matrix: HNO ₃ / HF Dilution 1:100	
IV-28-125ML IV-28-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	Mg	100
Al	100	Mn	100
As	100	Mo	100
B	100	Na	100
Ba	100	Ni	100
Be	100	Pb	100
Ca	100	Sb	100
Cd	100	Se	100
Co	100	Si	50
Cr	100	Sr	100
Cu	100	Ti	100
Fe	100	Tl	100
K	1,000	V	100
Li	100	Zn	100

Rev. 3.3 & 4.4 Calibrations – Standards may be used for either revision.

Calibration Standard		
CLPP-SPK-2		Matrix: HNO ₃ /Tartaric Acid Dilution 1:100
CLPP-SPK-2-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Sb	500	206.833

For use as Sb spike standard in EPA Contract Laboratory Program (CLP) methods.

Calibration Standard		
WW-CAL-1A		Matrix: HNO ₃ Dilution 1:100
WW-CAL-1A-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Ag	50	328.068
As	1,000	193.759
B	100	249.678
Ba	100	493.409
Ca	1,000	315.887
Cd	200	226.502
Cu	200	324.754
Mn	200	257.610
Se	500	196.090
Sr	100	421.552

NOTE: Sr does not exhibit spectral interference problems with any of the EPA Method 200.7 analytes.

For use as ICP-OES calibration standard I in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Calibration Standard		
WW-CAL-2		Matrix: HNO ₃ / HF Dilution 1:100
WW-CAL-2-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
K	2,000	766.491
Li	500	670.784
Mo	1,000	203.844
Na	1,000	588.995
Ti	1,000	334.941

For use as ICP-OES calibration standard II in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Calibration Standard		
WW-CAL-3		Matrix: HNO ₃ Dilution 1:100
WW-CAL-3-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Ce	200	413.765
Co	200	228.616
P	1,000	214.914
V	200	292.402

For use as ICP-OES calibration standard III in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Calibration Standard		
WW-CAL-4A		Matrix: HNO ₃ Dilution 1:100
WW-CAL-4A-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Al	1,000	308.215
Cr	500	205.552
Hg	200	194.227
Zn	500	213.856

For use as ICP-OES calibration standard IV (Part A) in EPA Method 200.7. Designed to be mixed with WW-CAL-4B at working concentrations. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Calibration Standard		
WW-CAL-4B		Matrix: HNO ₃ / HF Dilution 1:100
WW-CAL-4B-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
SiO₂	1,000	251.611
Sn	400	189.980

For use as ICP-OES calibration standard IV (Part B) in EPA Method 200.7. Designed to be mixed with WW-CAL-4A at working concentrations. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Rev. 3.3 & 4.4 Calibrations – Standards may be used for either revision.

Calibration Standard		
WW-CAL-5		Matrix: HNO ₃ Dilution 1:100
WW-CAL-5-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Be	100	313.042
Fe	1,000	259.940
Mg	1,000	279.079
Ni	200	231.604
Pb	1,000	220.353
Tl	500	190.864

For use as ICP-OES calibration standard V in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Rev. 3.3 & 4.4 Instrument Performance Checks – Standards may be used for either revision.

Instrument Performance Check		
WW-IPC-1		Matrix: HNO ₃ Dilution 1:100
WW-IPC-1-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Ag	25	328.068
Al	200	308.215
As	200	193.759
B	200	249.678
Ba	200	493.409
Be	200	313.042
Ca	200	315.887
Cd	200	226.502
Ce	200	413.765
Co	200	228.616
Cr	200	205.552
Cu	200	324.754
Fe	200	259.940
Hg	200	194.227
K	1,000	766.491
Li	200	670.784
Mg	200	279.079
Mn	200	257.610
Na	200	588.995
Ni	200	231.604
P	1,000	214.914
Pb	200	220.353
Se	200	196.090
Sr	200	421.552
Tl	200	190.864
V	200	292.402
Zn	200	213.856

Performance Check solution for EPA Method 200.7. Designed to be mixed with WW-IPC-2 at working concentrations. Based upon Revision 4.4 and suitable for all 200.7 versions.

Method 200.7

Instrument Performance Check		
WW-IPC-2		Matrix: HNO ₃ / HF Dilution 1:100
WW-IPC-2-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Mo	200	203.844
Sb	200	206.833
SiO₂	1,000	251.611
Sn	200	189.980
Ti	200	334.941

Performance Check solution for EPA Method 200.7. Designed to be mixed with WW-IPC-1 or WW-IPC-3 at working concentrations. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions

Instrument Performance Check		
WW-IPC-3		Matrix: HNO ₃ Dilution 1:100
WW-IPC-3-125ML		Volume: 125 mL
Analyte	µg/mL	λ(nm)
Ag	25	328.068
Al	200	308.215
As	200	193.759
B	200	249.678
Ba	200	493.409
Be	200	313.042
Ca	200	315.887
Cd	200	226.502
Co	200	228.616
Cr	200	205.552
Cu	200	324.754
Fe	200	259.940
K	1,000	766.491
Li	200	670.784
Mg	200	279.079
Mn	200	257.610
Na	200	588.995
Ni	200	231.604
P	1,000	214.914
Pb	200	220.353
Se	200	196.090
Sr	200	421.552
Tl	200	190.864
V	200	292.402
Zn	200	213.856

Performance Check solution for EPA Method 200.7. Designed to be mixed with WW-IPC-2 at working concentrations. Does not contain mercury (Hg). Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.



You can't put time in a bottle.
But we can save it in a bag.

Transpiration Control Technology (TCT)

Based on years of study and data evaluation, we improved the way we deliver our quality products. inorganicventures.com/tct.

Rev. 3.3 & 4.4 Laboratory Fortified Stocks – Standards may be used for either revision.

Laboratory Fortified Stock Solution		
WW-LFS-1	Matrix: HNO ₃ Dilution 1:100	
WW-LFS-1-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
Ag	7.5	328.068
Al	200	308.215
As	80	193.759
B	30	249.678
Ba	20	493.409
Be	20	313.042
Ca	100	315.887
Cd	20	226.502
Ce	200	413.765
Co	20	228.616
Cr	40	205.552
Cu	30	324.754
Fe	300	259.940
Hg	70	194.227
K	1,000	766.491
Li	20	670.784
Mg	200	279.079
Mn	20	257.610
Na	300	588.995
Ni	50	231.604
P	600	214.914
Pb	100	220.353
Se	200	196.090
Sr	20	421.552
Tl	200	190.864
V	30	292.402
Zn	20	213.856

Laboratory Fortified Stock Solution for EPA Method 200.7. Designed to be mixed with WW-LFS-2 at working concentrations. Suitable for use with all 200.7 versions.

Laboratory Fortified Stock Solution		
WW-LFS-2	Matrix: HNO ₃ / HF Dilution 1:100	
WW-LFS-2-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
Mo	40	203.844
Sb	80	206.833
SiO ₂	200	251.611
Sn	70	189.980
Ti	20	334.941

Laboratory Fortified Stock Solution for EPA Method 200.7. Designed to be mixed with WW-LFS-1 at working concentrations. Suitable for use with all 200.7 versions.

Rev. 3.3 & 4.4 Quality Controls – Standards may be used for either revision.

Quality Control Standard†		
QCP-QCS-1	Matrix: HNO ₃ Dilution 1:100	
QCP-QCS-1-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
Ag	25	328.068
Al	100	308.215
As	200	193.759
B	100	249.678
Ba	100	493.409
Be	100	313.042
Ca	100	315.887
Cd	100	226.502
Ce	100	413.765
Co	100	228.616
Cr	100	205.552
Cu	100	324.754
Fe	100	259.940
Hg	200	194.227
K	500	766.491
Li	100	670.784
Mg	100	279.079
Mn	100	257.610
Na	100	588.995
Ni	100	231.604
P	500	214.914
Pb	200	220.353
Se	100	196.090
Sr	100	421.552
Tl	500	190.864
V	100	292.402
Zn	100	213.856

For use as ICP-OES QC standard in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Quality Control Standard†		
QCP-QCS-2	Matrix: HNO ₃ / HF Dilution 1:100	
QCP-QCS-2-125ML	Volume: 125 mL	
Analyte	µg/mL	λ(nm)
Mo	100	203.844
Sb	200	206.833
SiO ₂	500	251.611
Sn	500	189.980
Ti	100	334.941

For use as ICP-OES QC standard in EPA Method 200.7. Based upon Revisions 3.3 and 4.4 and suitable for all 200.7 versions.

Quality Control Standard†			
IV-7		Matrix: HNO ₃ / HF Dilution 1:100	
IV-7-125ML	Volume: 125 mL		
Analyte	µg/mL	Analyte	µg/mL
Ag	100	K	1,000
Al	100	Na	100
B	100	Si	50
Ba	100		

For use as a certified reference standard in ICP applications. Can be diluted with other standards to working concentrations.

†Manufactured from in-house Second Source concentrates, whenever possible.

Rev. 3.3 & 4.4 Quality Controls – Standards may be used for either revision.

For use as a certified reference standard in ICP applications. Can be diluted with other standards to working concentrations.

Quality Control Standard [†]			
IV-19		Matrix: HNO ₃ / HF Dilution 1:100	
IV-19-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	100	Mo	100
Be	100	Ni	100
Ca	100	Pb	100
Cd	100	Sb	100
Co	100	Se	100
Cr	100	Ti	100
Cu	100	Tl	100
Fe	100	V	100
Mg	100	Zn	100
Mn	100		

Quality Control Standard [†]			
IV-21		Matrix: HNO ₃ / HF Dilution 1:100	
IV-21-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
As	100	Mo	100
Be	100	Ni	100
Ca	100	Pb	100
Cd	100	Sb	100
Co	100	Se	100
Cr	100	Sr	100
Cu	100	Ti	100
Fe	100	Tl	100
Li	100	V	100
Mg	100	Zn	100
Mn	100		

JY HORIBA Jobin Yvon**M** Merck/MilliporeSigma**PE** Perkin Elmer[†]Manufactured from in-house Second Source concentrates, whenever possible.

Quality Control Standard [†]			
IV-26		Matrix: HNO ₃ / HF Dilution 1:100	
IV-26-125ML IV-26-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	Mg	100
Al	100	Mn	100
As	100	Mo	100
B	100	Na	100
Ba	100	Ni	100
Be	100	Pb	100
Ca	100	Sb	100
Cd	100	Se	100
Co	100	Si	50
Cr	100	Ti	100
Cu	100	Tl	100
Fe	100	V	100
K	1,000	Zn	100

Quality Control Standard [†]			
IV-28		Matrix: HNO ₃ / HF Dilution 1:100	
IV-28-125ML IV-28-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	100	Mg	100
Al	100	Mn	100
As	100	Mo	100
B	100	Na	100
Ba	100	Ni	100
Be	100	Pb	100
Ca	100	Sb	100
Cd	100	Se	100
Co	100	Si	50
Cr	100	Sr	100
Cu	100	Ti	100
Fe	100	Tl	100
K	1,000	V	100
Li	100	Zn	100

Standards for Method 200.8 are designed for use with ICP-MS. Custom EPA standards are available upon request.

Rev. 4.4 & 5.4 Calibration – See individual products for recommended revisions.

Calibration Standard			
2008CAL-1		Matrix: HNO ₃ / HF Dilution 1:100	
2008CAL-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Mo	20	Sb	20

Designed for Rev. 4.4 and 5.4.

For use as ICP-MS calibration standard in EPA Method 200.8.

Calibration Standard			
2008CAL-2		Matrix: HNO ₃ Dilution 1:100	
2008CAL-2-125ML		Volume: 125 mL	
2008CAL-2-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	20	Mn	20
Al	20	Ni	20
As	20	Pb	20
Ba	20	Se	20
Be	20	Th	20
Cd	20	Tl	20
Co	20	U	20
Cr	20	V	20
Cu	20	Zn	20

Designed for Rev. 4.4.

For use as ICP-MS calibration standard in EPA Method 200.8.

Calibration Standard	
WW-MSCAL-1 	Matrix: HNO ₃ Dilution 1:1,000
WW-MSCAL-1-125ML	Volume: 125 mL
Analyte	µg/mL
Hg	5

Designed for Rev. 5.4.

For use as ICP-MS calibration standard in EPA Method 200.8. Can be combined with WW-MSCAL-2 at working concentrations if Au is used to stabilize Hg. Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

Calibration Standard			
WW-MSCAL-2		Matrix: HNO ₃ Dilution 1:100	
WW-MSCAL-2-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	20	Mn	20
Al	20	Ni	20
As	20	Pb	20
Ba	20	Se	100
Be	20	Th	20
Cd	20	Tl	20
Co	20	U	20
Cr	20	V	20
Cu	20	Zn	20

Designed for Rev. 5.4.

For use as ICP-MS calibration standard in EPA Method 200.8. Uranium is isotopically depleted. Does not contain mercury (Hg). Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

Standards for Method 200.8 are designed for use with ICP-MS. Custom EPA standards are available upon request.

Rev. 4.4 & 5.4 Calibration – See individual products for recommended revisions.

Mercury Standard	
MSHG-1PPM	Matrix: HCl
MSHG-1PPM-125ML MSHG-1PPM-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Hg	1

Mercury Standard	
IV-STOCK-73	Matrix: 10% v/v HCl
IV-STOCK-73-125ML	Volume: 125 mL
Analyte	µg/L*
Hg	100

Mercury Standard	
IV-STOCK-72	Matrix: 7% v/v HNO ₃
IV-STOCK-72-125ML	Volume: 125 mL
Analyte	µg/L*
Hg	100

*Parts per billion

Rev. 4.4 & 5.4 Internal Standards

Internal Standard			
2008ISS		Matrix: HNO ₃ Dilution 1:100 to 1:1,000	
2008ISS-125ML 2008ISS-500ML		Volume: 125 mL Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	20	Tb	20
In	20	Y	20
Sc	20		

Designed for Rev. 4.4 and 5.4. Recommended working level is 200 µg/L for Rev. 4.4; 20-200 µg/L for Rev. 5.4. Use this solution with CGAUN1 for Rev. 5.4 if Hg is to be determined by direct analysis.

Mercury Preservation Solution	
CGAUN1	Matrix: HNO ₃ Dilution 1:100
CGAUN1-30ML CGAUN1-125ML CGAUN1-500ML	Volume: 30 mL Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Au	1,000

Designed for Rev. 5.4. Add an aliquot of this solution to 2008ISS, sufficient to provide a concentration of 100 µg/L in the final dilution of all blanks, calibration standards, and samples.

Rev. 4.4 & 5.4 Quality Controls

Quality Control Standard†			
QCP-QCS-3		Matrix: HNO ₃ Dilution 1:100	
QCP-QCS-3-125ML		Volume: 125 mL	
QCP-QCS-3-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Mn	10
Al	10	Mo	10
As	10	Na	10
Ba	10	Ni	10
Be	10	Pb	10
Ca	10	Sb	10
Cd	10	Se	50
Co	10	Th	10
Cr	10	Tl	10
Cu	10	U	10
Fe	10	V	10
K	10	Zn	10
Mg	10		

Designed for Rev. 4.4 and 5.4.

For use as ICP-MS quality control standard in EPA Method 200.8. Uranium is isotopically depleted. Does not contain mercury (Hg). Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

Quality Control Standard†	
QCP-QCS-4	Matrix: HNO ₃ Dilution 1:100
QCP-QCS-4-125ML	Volume: 125 mL
Analyte	µg/mL
Hg	5

Designed for Rev. 4.4 and 5.4.

For use as ICP-MS mercury (Hg) quality control standard in EPA Method 200.8. Can be combined with QCP-QCS-3 at working concentrations if Au is used to stabilize Hg, though trace chloride from Au may cause Ag stability problems. Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions, and the Silver Chemical Stability article for more information about Ag stability in the presence of chloride.

†Manufactured from in-house Second Source concentrates, whenever possible.

Rev. 4.4 & 5.4 Tuning

Tuning Solution			
2008TS 		Matrix: HNO ₃ Dilution 1:100 to 1:1,000	
2008TS-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Be	10	Mg	10
Co	10	Pb	10
In	10		

Designed for Rev. 4.4 and 5.4.

For use as ICP-MS tuning solution in EPA Method 200.8.



For use in EPA Toxicity Characteristic Leachate Procedure (TCLP). Custom EPA standards are available upon request.

TCLP Hg Standard	
TCLP-AA-HG	Matrix: HNO ₃ Dilution: As required
TCLP-AA-HG-125ML	Volume: 125 mL
Analyte	µg/mL
Hg	20


For use in EPA Toxicity Characteristic Leachate Procedure (TCLP).

TCLP Standard			
TCLP-1REV		Matrix: HNO ₃ Dilution: As required	
TCLP-1REV-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	25	Cr	25
As	25	Pb	25
Ba	500	Se	5
Cd	5		

For use in EPA Toxicity Characteristic Leachate Procedure (TCLP).

Don't see what you need?

Contact us with the solution part number and instrument manufacturer you're seeking, and we'll check our extensive library of solutions.



Standards for Method 6020 are designed for use with ICP-MS. Custom EPA standards are available upon request.

CLP-M Version 8

Calibration Standard			
6020CAL-1		Matrix: HNO ₃ / HF Dilution 1:100	
6020CAL-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	20	K	20
Al	20	Mg	20
As	20	Mn	20
Ba	20	Na	20
Be	20	Ni	20
Ca	20	Pb	20
Cd	20	Sb	20
Co	20	Se	20
Cr	20	Tl	20
Cu	20	V	20
Fe	20	Zn	20

For use as ICP-MS calibration standard in EPA Method 6020. Does not contain mercury (Hg). Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

Interference Check Standard			
6020ICS-8A		Matrix: HNO ₃ Dilution 1:10	
6020ICS-8A-125ML		Volume: 125 mL	
6020ICS-8A-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	1,000	Mg	1,000
C	2,000	Mo	20
Ca	3,000	Na	2,500
Cl-	18,000	P	1,000
Fe	2,500	S	1,000
K	1,000	Ti	20

For evaluating ICP-MS interferences and corrections in EPA Method 6020A or 6020B. Based upon 6020A Revision 1 Solution A.

Internal Standard			
6020ISS		Matrix: HNO ₃ Dilution 1:100	
6020ISS-125ML		Volume: 125 mL	
6020ISS-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	10	Rh	10
Ho	10	Sc	10
In	10	Tb	10
⁶ Li	10	Y	10

Internal standard for ICP-MS in all versions of EPA Method 6020.

Spike Standard – Soil			
6020SPK-S		Matrix: HNO ₃ Dilution 1:100	
6020SPK-S-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Ni	25
As	10	Pb	20
Ba	50	Sb	20
Be	5	Se	5
Cd	10	Tl	5
Co	20	V	30
Cr	50	Zn	50
Cu	50		

Matrix spike for solid samples.

Standards for Method 6020 are designed for use with ICP-MS. Custom EPA standards are available upon request.

CLP-M Version 8

Spike Standard – Water			
6020SPK-W		Matrix: HNO ₃ Dilution 1:100	
6020SPK-W-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	5	Mn	20
As	10	Ni	20
Ba	50	Pb	10
Be	5	Sb	20
Cd	5	Se	5
Co	20	Tl	5
Cr	20	V	20
Cu	20	Zn	50
Fe	100		

Matrix spike for aqueous samples.

CLP-M Version 9

Calibration Standard			
6020CAL-1		Matrix: HNO ₃ / HF Dilution 1:100	
6020CAL-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	20	K	20
Al	20	Mg	20
As	20	Mn	20
Ba	20	Na	20
Be	20	Ni	20
Ca	20	Pb	20
Cd	20	Sb	20
Co	20	Se	20
Cr	20	Tl	20
Cu	20	V	20
Fe	20	Zn	20

For use as ICP-MS calibration standard in EPA Method 6020. Does not contain mercury (Hg). Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

Tuning Solution			
6020TS		Matrix: HNO ₃ Dilution 1:100	
6020TS-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Co	10	Li	10
In	10	Tl	10

For use as general tuning solution suitable for numerous ICP-MS designs and models. Covers mass range from Li to Tl. Certified reference material that may also be used for calibration.

Interference Check Standard			
6020ICS-9A		Matrix: HNO ₃ Dilution 1:10	
6020ICS-9A-125ML		Volume: 125 mL	
6020ICS-9A-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	1,000	Mg	1,000
C	2,000	Mo	20
Ca	3,000	Na	2,500
Cl-	21,215	P	1,000
Fe	2,500	S	1,000
K	1,000	Ti	20

For evaluating ICP-MS interferences and corrections in EPA Method 6020A or 6020B. Based upon 6020A Revision 1 Solution A. Suitable for analyses with higher chloride contents.

Interference Check Standard			
6020ICS-9B		Matrix: HNO ₃ Dilution 1:100	
6020ICS-9B-125ML			
Volume: 125 mL			
Analyte	µg/mL	Analyte	µg/mL
Ag	5	Mn	20
As	10	Ni	20
Cd	10	Se	10
Co	20	V	20
Cr	20	Zn	10
Cu	20		

For evaluating ICP-MS interferences and corrections in EPA Method 6020A or 6020B. Based upon 6020A Revision 1 Solution B. Does not contain mercury (Hg). Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

Internal Standard			
6020ISS		Matrix: HNO ₃ Dilution 1:100	
6020ISS-125ML 6020ISS-500ML			
Volume: 125 mL		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	10	Rh	10
Ho	10	Sc	10
In	10	Tb	10
⁶ Li	10	Y	10

Internal standard for ICP-MS in all versions of EPA Method 6020.

Spike Standard – Soil			
6020SPK-S		Matrix: HNO ₃ Dilution 1:100	
6020SPK-S-125ML			
Volume: 125 mL			
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Ni	25
As	10	Pb	20
Ba	50	Sb	20
Be	5	Se	5
Cd	10	Tl	5
Co	20	V	30
Cr	50	Zn	50
Cu	50		

Matrix spike for solid samples.

Spike Standard – Water			
6020SPK-W		Matrix: HNO ₃ Dilution 1:100	
6020SPK-W-125ML			
Volume: 125 mL			
Analyte	µg/mL	Analyte	µg/mL
Ag	5	Mn	20
As	10	Ni	20
Ba	50	Pb	10
Be	5	Sb	20
Cd	5	Se	5
Co	20	Tl	5
Cr	20	V	20
Cu	20	Zn	50
Fe	100		

Matrix spike for aqueous samples.

Tuning Solution			
6020TS		Matrix: HNO ₃ Dilution 1:100	
6020TS-125ML			
Volume: 125 mL			
Analyte	µg/mL	Analyte	µg/mL
Co	10	Li	10
In	10	Tl	10

For use as general tuning solution suitable for numerous ICP-MS designs and models. Covers mass range from Li to Tl. Certified reference material that may also be used for calibration.



REV. 0

Calibration Standard			
6020CAL-1		Matrix: HNO ₃ / HF Dilution 1:100	
6020CAL-1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	20	K	20
Al	20	Mg	20
As	20	Mn	20
Ba	20	Na	20
Be	20	Ni	20
Ca	20	Pb	20
Cd	20	Sb	20
Co	20	Se	20
Cr	20	Tl	20
Cu	20	V	20
Fe	20	Zn	20

For use as ICP-MS calibration standard in EPA Method 6020. Does not contain mercury (Hg). Please see the Mercury Chemical Stability article for more information regarding accurate Hg analyses in multi-element solutions.

Internal Standard			
6020ISS		Matrix: HNO ₃ Dilution 1:100	
6020ISS-125ML		Volume: 125 mL	
6020ISS-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Bi	10	Rh	10
Ho	10	Sc	10
In	10	Tb	10
⁶ Li	10	Y	10

Internal standard for ICP-MS in all versions of EPA Method 6020.

Interference Check Standard			
6020ICS-0A		Matrix: HNO ₃ Dilution 1:10	
6020ICS-0A-125ML		Volume: 125 mL	
6020ICS-0A-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Al	1,000	Mg	1,000
C	2,000	Mo	20
Ca	1,000	Na	1,000
Cl-	10,000	P	1,000
Fe	1,000	S	1,000
K	1,000	Ti	20

For evaluating ICP-MS interferences and corrections in EPA Method 6020. Based upon Revision 0 Solution A.

Interference Check Standard			
6020ICS-0B		Matrix: HNO ₃ Dilution 1:100	
6020ICS-0B-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	2	Cu	2
As	2	Mn	2
Cd	2	Ni	2
Co	2	Zn	2
Cr	2		

For evaluating ICP-MS interferences and corrections in EPA Method 6020. Based upon Revision 0 Solution B.

Spike Standard – Soil			
6020SPK-S		Matrix: HNO ₃ Dilution 1:100	
6020SPK-S-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	10	Ni	25
As	10	Pb	20
Ba	50	Sb	20
Be	5	Se	5
Cd	10	Tl	5
Co	20	V	30
Cr	50	Zn	50
Cu	50		

Matrix spike for solid samples.

REV. 0

Spike Standard – Water			
6020SPK-W		Matrix: HNO ₃ Dilution 1:100	
6020SPK-W-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ag	5	Mn	20
As	10	Ni	20
Ba	50	Pb	10
Be	5	Sb	20
Cd	5	Se	5
Co	20	Tl	5
Cr	20	V	20
Cu	20	Zn	50
Fe	100		

Matrix spike for aqueous samples.

Tuning Solution			
6020TS		Matrix: HNO ₃ Dilution 1:100	
6020TS-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Co	10	Li	10
In	10	Tl	10

For use as general tuning solution suitable for numerous ICP-MS designs and models. Covers mass range from Li to Tl. Certified reference material that may also be used for calibration.



Now, YOU control the expiration date.

With TCT, concerns about shipping and storage conditions are eliminated, as transpiration is no longer an issue. This saves you money and simplifies research by removing the need to constantly inventory and restock CRMs.

Find out more at inorganicventures.com/tct.

ION CHROMATOGRAPHY

If you've been searching for an atypical Ion Chromatography standard, look no further.

Over the years, we've developed the most complete line of IC standards on the market. Our technicians have stabilized more than a dozen rare anion and cation standards that you won't find anywhere else.



Contents

Anion Standards	86
Cation Standards	88
Multi-Ion Standards	89
Eluent Concentrates.....	90
EPA Standards	91
Need a Custom CRM?	16

- ✓ **Five-year shelf life**
- ✓ **Traceable to NIST SRMs**
- ✓ **Produced under ISO 9001**
- ✓ **Produced under ISO 17025**
- ✓ **Produced under ISO 17034**
- ✓ **Assayed by validated Wet Chemical procedures**
- ✓ **Assayed by validated IC procedures**

1,000 µg/mL Anions

Custom anion standards are available upon request.

1,000 µg/mL

ANALYTE	MATRIX	STARTING MATERIAL	VOLUME	CATALOG #
Acetate, C ₂ H ₃ O ₂ ⁻	H ₂ O	Sodium acetate	125 mL 500 mL	ICOAC1-125ML ICOAC1-500ML
Adipate, C ₆ H ₈ O ₄ ⁻²	H ₂ O	Adipic acid	125 mL	ICADP1-125ML
Benzoate, C ₆ H ₅ CO ₂ ⁻	H ₂ O	Benzoic acid	125 mL	ICBEN1-125ML
Bromate, BrO ₃ ⁻	H ₂ O	KBrO ₃	125 mL 500 mL	ICBRO31-125ML ICBRO31-500ML
Bromide, Br ⁻	H ₂ O	KBr	125 mL 500 mL	ICBR1-125ML ICBR1-500ML
Butyrate, C ₄ H ₇ O ₂ ⁻	H ₂ O	Butyric acid	125 mL	ICBTR1-125ML
Carbonate, CO ₃ ⁻²	H ₂ O	Na ₂ CO ₃	125 mL 500 mL	ICCO31-125ML ICCO31-500ML
Chlorate, ClO ₃ ⁻	H ₂ O	KClO ₃	125 mL 500 mL	ICCLO31-125ML ICCLO31-500ML
Chloride, Cl ⁻	H ₂ O	KCl	125 mL 500 mL	ICCL1-125ML ICCL1-500ML
Chlorite, ClO ₂ ⁻	H ₂ O	NaClO ₂	125 mL 500 mL	ICCLO21-125ML ICCLO21-500ML
Chromate, CrO ₄ ⁻²	H ₂ O	(NH ₄) ₂ Cr ₂ O ₇	125 mL	ICCRO41-125ML
Citrate, C ₆ H ₅ O ₇ ⁻³	H ₂ O	Citric acid	125 mL 500 mL	ICCIT1-125ML ICCIT1-500ML
Cyanide, NaCN	H ₂ O	Sodium cyanide	20 mL	CN-1000-25-20ML
Fluoride, F ⁻	H ₂ O	NaF	125 mL 500 mL	ICF1-125ML ICF1-500ML
Formate, HCO ₂ ⁻	H ₂ O	Sodium formate	125 mL 500 mL	ICHCO1-125ML ICHCO1-500ML
Glutarate, C ₅ H ₆ O ₄ ⁻²	H ₂ O	Glutaric acid	125 mL	ICGTR1-125ML
Glycolate, C ₂ H ₃ O ₃ ⁻	H ₂ O	Glycolic acid	125 mL	ICGLY1-125ML
Iodide, I ⁻	H ₂ O / stabilizer	NH ₄ I	125 mL 500 mL	ICI1-125ML ICI1-500ML
Lactate, C ₃ H ₅ O ₃ ⁻	H ₂ O	Lactic acid	125 mL	ICLCT1-125ML
Malate, C ₄ H ₄ O ₅ ⁻²	H ₂ O	Malic acid	125 mL	ICMLA1-125ML
Maleate, C ₄ H ₂ O ₄ ⁻²	H ₂ O	Maleic acid	125 mL	ICMLE1-125ML
Malonate, C ₃ H ₂ O ₄ ⁻²	H ₂ O	Malonic acid	125 mL	ICMLO1-125ML
Methanesulfonate, CH ₃ SO ₃ ⁻	H ₂ O	Methanesulfonic acid	125 mL	ICMSA1-125ML
Nitrate, NO ₃ ⁻	H ₂ O	NaNO ₃	125 mL 500 mL	ICNO31-125ML ICNO31-500ML
Nitrate as Nitrogen	H ₂ O	NaNO ₃	125 mL 500 mL	ICNNO31-125ML ICNNO31-500ML
Nitrilotriacetate, NC ₆ H ₆ O ₆ ⁻³	H ₂ O	Nitrilotriacetic acid	125 mL	ICNTA1-125ML
Nitrite, NO ₂ ⁻	H ₂ O	NaNO ₂	125 mL 500 mL	ICNO21-125ML ICNO21-500ML
Nitrite as Nitrogen	H ₂ O	NaNO ₂	125 mL 500 mL	ICNNO21-125ML ICNNO21-500ML

ION CHROMATOGRAPHY ANION STANDARDS

1,000 µg/mL Anions

Custom anion standards are available upon request.

1,000 µg/mL

ANALYTE	MATRIX	STARTING MATERIAL	VOLUME	CATALOG #
Oxalate, $C_2O_4^{-2}$	H ₂ O	Sodium oxalate	125 mL	ICXA1-125ML
			500 mL	ICXA1-500ML
Perchlorate, ClO_4^-	H ₂ O	KClO ₄	125 mL	ICCL041-125ML
			500 mL	ICCL041-500ML
Phosphate, PO_4^{-3}	H ₂ O	NH ₄ H ₂ PO ₄	125 mL	ICP041-125ML
			500 mL	ICP041-500ML
Phosphate as Phosphorus	H ₂ O	NH ₄ H ₂ PO ₄	125 mL	ICPPO41-125ML
			500 mL	ICPPO41-500ML
Phthalate, $C_6H_4(CO_2)_2^{-2}$	H ₂ O	Potassium hydrogen phthalate	125 mL	ICKHP1-125ML
Propionate, $C_2H_5CO_2^-$	H ₂ O	Sodium propionate	125 mL	ICOPR1-125ML
Succinate, $C_4H_4O_4^{-2}$	H ₂ O	Succinic acid	125 mL	ICSCC1-125ML
Sulfate, SO_4^{-2}	H ₂ O	K ₂ SO ₄	125 mL	ICSO41-125ML
			500 mL	ICSO41-500ML
Tartrate, $C_4H_4O_6^{-2}$	H ₂ O	Tartaric acid	125 mL	ICTRTR1-125ML
Thiocyanate, SCN ⁻	H ₂ O	KSCN	125 mL	ICSCN1-125ML
Thiosulfate, $S_2O_3^{-2}$	H ₂ O	Sodium thiosulfate	125 mL	ICS2031-125ML
			500 mL	ICS2031-500ML

10,000 µg/mL

Custom anion standards are available upon request.

10,000 µg/mL

ANALYTE	MATRIX	STARTING MATERIAL	VOLUME	CATALOG #
Chloride, Cl ⁻	H ₂ O	KCl	125 mL	ICCL10-125ML
			500 mL	ICCL10-500ML
Sulfate, SO_4^{-2}	H ₂ O	K ₂ SO ₄	125 mL	ICSO410-125ML
			500 mL	ICSO410-500ML

100 ppm

Custom anion standards are available upon request.

100 ppm

ANALYTE	MATRIX	STARTING MATERIAL	VOLUME	CATALOG #
Nitrite, NO ₂ ⁻	H ₂ O	100	125 mL	ICNO2-100PPM-125ML

Custom cation standards are available upon request.

1,000 µg/mL

ANALYTE	MATRIX	STARTING MATERIAL	VOLUME	CATALOG #
3-Methoxypropylamine $\text{CH}_3\text{O}(\text{CH}_2)_3\text{NH}_2$	HCl	3-Methoxypropylamine	125 mL	ICMPA1-125ML
Ammonium, NH_4^+	H_2O	NH_4Cl	125 mL 500 mL	ICNH41-125ML ICNH41-500ML
Ammonium as Nitrogen	H_2O	NH_4Cl	125 mL 500 mL	ICNNH41-125ML ICNNH41-500ML
Barium, Ba^{+2}	HNO_3	$\text{Ba}(\text{NO}_3)_2$	125 mL	ICBA1-125ML
Calcium, Ca^{+2}	HNO_3	CaO	125 mL 500 mL	ICCA1-125ML ICCA1-500ML
Cesium, Cs^+	HNO_3	CsNO_3	125 mL	ICCS1-125ML
Diethanolamine, $(\text{HOCH}_2\text{CH}_2)_2\text{NH}$	H_2O	Diethanolamine	125 mL	ICDEA1-125ML
Dimethylamine, $\text{NH}(\text{CH}_3)_2$	HCl	Dimethylamine	125 mL	ICDMA1-125ML
Lithium, Li^+	HNO_3	Li_2CO_3	125 mL	ICLI1-125ML
Magnesium, Mg^{+2}	HNO_3	Mg metal	125 mL 500 mL	ICMG1-125ML ICMG1-500ML
Monoethanolamine, $\text{HOCH}_2\text{CH}_2\text{NH}_2$	H_2O	Monoethanolamine	125 mL 500 mL	ICMEA1-125ML ICMEA1-500ML
Monomethylamine, NH_2CH_3	HCl	Monomethylamine	125 mL	ICMMA1-125ML
Potassium, K^+	HNO_3	KNO_3	125 mL 500 mL	ICK1-125ML ICK1-500ML
Rubidium, Rb^+	HNO_3	RbNO_3	125 mL	ICRB1-125ML
Sodium, Na^+	HNO_3	Na_2CO_3	125 mL 500 mL	ICNA1-125ML ICNA1-500ML
Strontium, Sr^{+2}	HNO_3	SrCO_3	125 mL	ICSR1-125ML
Tetramethylammonium, $\text{N}^+(\text{CH}_3)_4$	H_2O	Tetramethylammonium hydroxide	125 mL	ICTMAH1-125ML
Triethanolamine, $(\text{HOCH}_2\text{CH}_2)_3\text{N}$	H_2O	Triethanolamine	125 mL	ICTEA1-125ML
Trimethylamine, $(\text{CH}_3)_3\text{N}$	HCl	Trimethylamine	125 mL	ICTMA1-125ML

Anion Calibration Standard			
IC-FAS-1A D I		Matrix: H ₂ O	
IC-FAS-1A-125ML		Volume: 125 mL	
IC-FAS-1A-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Br ⁻	100	NO ₂ ⁻	100
Cl ⁻	30	PO ₄ ⁻³	150
F ⁻	20	SO ₄ ⁻²	150
NO ₃ ⁻	100		

For anion analysis of water samples by ion chromatography (IC). Contains 7 "common anions" as defined by EPA and Standard Methods.

Cation Calibration Standard			
IC-SCS1 D I		Matrix: HNO ₃	
IC-SCS1-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Ca ⁺²	1,000	Mg ⁺²	200
K ⁺	200	Na ⁺	200
Li ⁺	50	NH ₄ ⁺	400

Used for daily calibration.

For cation analysis of water samples by ion chromatography (IC).

Cation Calibration Standard			
IV-STOCK-7 M I		Matrix: HNO ₃	
IV-STOCK-7-125ML		Volume: 125 mL	
IV-STOCK-7-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Ba ⁺²	100	Mn ⁺²	100
Ca ⁺²	100	Na ⁺	100
K ⁺	100	NH ₄ ⁺	100
Li ⁺	100	Sr ⁺²	100
Mg ⁺²	100		

Used for daily calibration.

For use as a certified reference standard in ion chromatography (IC) applications.

- D** DIONEX
- M** Merck/MilliporeSigma
- I** Common Multi-Ion Standards

Anion Calibration Standard			
IV-STOCK-59 I		Matrix: H ₂ O	
IV-STOCK-59-125ML		Volume: 125 mL	
IV-STOCK-59-500ML		Volume: 500 mL	
Analyte	µg/mL	Analyte	µg/mL
Br ⁻	1,000	NO ₂ ⁻	1,000
Cl ⁻	1,000	PO ₄ ⁻³	1,000
F ⁻	1,000	SO ₄ ⁻²	1,000
NO ₃ ⁻	1,000		

For use as a certified reference standard in ion chromatography (IC) applications.

Anion Mix A	
IV-STOCK-61	Matrix: H ₂ O
IV-STOCK-61-125ML	Volume: 125 mL
IV-STOCK-61-500ML	Volume: 500 mL
Analyte	Range
Br ⁻	20
F ⁻	10
NO ₂ ⁻	20
SO ₄ ⁻²	30
Cl ⁻	20
NO ₃ ⁻	20
PO ₄ ⁻³	30

For use as a certified reference standard in ion chromatography (IC) applications.

Cation Mix B	
IV-STOCK-62	Matrix: H ₂ O
IV-STOCK-62-125ML	Volume: 125 mL
Analyte	Range
Ca	2.0
K	2.5
Li	0.2
Mg	2.0
Na	1.5
NH ₄ ⁺	1.5

For use as a certified reference standard in ion chromatography (IC) applications.

Anion Mix 4	
IV-STOCK-63	Matrix: H ₂ O
IV-STOCK-63-125ML	Volume: 125 mL
Analyte	Range
Br⁻	40
F⁻	20
NO₂⁻	40
Cl⁻	40
NO₃⁻	40
SO₄⁻²	40

For use as a certified reference standard in ion chromatography (IC) applications.

Anion Mix 5	
IV-STOCK-64	Matrix: H ₂ O
IV-STOCK-64-125ML IV-STOCK-64-500ML	Volume: 125 mL Volume: 500 mL
Analyte	Range
Br⁻	50
Cl⁻	50
F	25
NO₃⁻	50
NO₂⁻	50
PO₄⁻³	50
SO₄⁻²	50

For use as a certified reference standard in ion chromatography (IC) applications.

Custom eluent concentrates are available upon request. All Eluents supplied with Product Information Sheet.

0.18 M Sodium Carbonate/0.17 M Sodium Bicarbonate D		
ELUENT1817-100ML	Volume: 100 mL	Matrix: H ₂ O
ELUENT1817-500ML	Volume: 500 mL	Dilution: 1:100
For preparation of 1.8 mM CO ₃ ⁻² / 1.7 mM HCO ₃ ⁻ eluent.		

This solution is a reagent and is not intended to be used as a certified reference material. Concentrate (100x); To prepare 1.8 mM carbonate/1.7 mM bicarbonate eluent for IC applications; ISO 17034 Reference Material.

0.35 M Sodium Carbonate/0.10 M Sodium Bicarbonate D		
ELUENT3510-100ML	Volume: 100 mL	Matrix: H ₂ O
ELUENT3510-500ML	Volume: 500 mL	Dilution: 1:100
For preparation of 3.5 mM CO ₃ ⁻² / 1.0 mM HCO ₃ ⁻ eluent.		

This solution is a reagent and is not intended to be used as a certified reference material. Concentrate (100x); To prepare 3.5 mM carbonate/1.0 mM bicarbonate eluent for IC applications; ISO 17034 Reference Material.

0.45 M Sodium Carbonate/0.14 M Sodium Bicarbonate D		
ELUENT4514-500ML	Volume: 500 mL	Matrix: H ₂ O Dilution: 1:100
For preparation of 3.5 mM CO ₃ ⁻² / 1.0 mM HCO ₃ ⁻ eluent.		

This solution is a reagent and is not intended to be used as a certified reference material. Concentrate (100x); To prepare 3.5 mM carbonate/1.0 mM bicarbonate eluent for IC applications; ISO 17034 Reference Material.

0.5 M Sodium Bicarbonate D		
BICARB-100ML	Volume: 100 mL	Matrix: H ₂ O
BICARB-500ML	Volume: 500 mL	Dilution: 1:100
For preparation of various CO ₃ ⁻² / HCO ₃ ⁻ eluents.		

This solution is a reagent and is not intended to be used as a certified reference material. Concentrate; To prepare carbonate/bicarbonate eluent for IC applications; ISO 17034 Reference Material.

D DIONEX

0.5 M Sodium Carbonate

D

CARB-100ML CARB-500ML	Volume: 100 mL Volume: 500 mL	Matrix: H ₂ O Dilution: 1:100
For preparation of various CO ₃ ²⁻ / HCO ₃ ⁻ eluents.		

This solution is a reagent and is not intended to be used as a certified reference material. Concentrate; To prepare carbonate/bicarbonate eluent for IC applications; ISO 17034 Reference Material; Supplied with Product Information Sheet.

1.8 M Methanesulfonic Acid

MSAELUENT-100ML MSAELUENT-500ML	Volume: 100 mL Volume: 500 mL	Matrix: H ₂ O Dilution: 1:100
For preparation of 18 mM CH ₃ SO ₃ H eluent for analyzing cations.		

This solution is a reagent and is not intended to be used as a certified reference material. Concentrate (100x); To prepare 18 mM methanesulfonic acid eluent for IC applications; ISO 17034 Reference Material; Supplied with Product Information Sheet.

300.0 Rev. 2.1 Part A / 300.1 Part A Custom EPA standards are available upon request.

0.18 M Sodium Carbonate/0.17 M Sodium Bicarbonate

D

ELUENT1817-100ML ELUENT1817-500ML	Volume: 100 mL Volume: 500 mL	Matrix: H ₂ O Dilution 1:100
For preparation of 1.8 mM CO ₃ ²⁻ / 1.7 mM HCO ₃ ⁻ eluent.		

This solution is a reagent and is not intended to be used as a certified reference material. Concentrate (100x); To prepare 1.8 mM carbonate/1.7 mM bicarbonate eluent for IC applications; ISO 17034 Reference Material; Supplied with Product Information Sheet.

Calibration Standard

300-CAL-A-125ML 300-CAL-A-500ML		Volume: 125 mL Volume: 500 mL	Matrix: H ₂ O Dilution 1:10 to 1:100
Analyte	µg/mL	Analyte	µg/mL
Br-	100	Nitrite as Nitrogen	30
Cl-	30	Nitrate as Nitrogen	25
F-	20	Phosphate as Phosphorus	50
SO₄⁻²	150		

For use as ion chromatography calibration standard in EPA Method 300.0. Also suitable for use as common anions standard in EPA Method 300.1.

Dichloroacetate Standard

ICDCA-S-125ML ICDCA-S-500ML	Volume: 125 mL Volume: 500 mL	Matrix: H ₂ O
Analyte	µg/mL	
Cl₂HC₂O₂⁻	500	

For use as a surrogate analyte.

For use as a surrogate analyte in ion chromatography (IC) analysis.

D DIONEX

300.0 Rev. 2.1 Part A / 300.1 Part A Custom EPA standards are available upon request.

Laboratory Fortification Stock Standard			
300-LFS-A-125ML		Volume: 125 mL	Matrix: H ₂ O Dilution 1:100 to 1:1,000
Analyte	µg/mL	Analyte	µg/mL
Br-	1,000	Nitrite as Nitrogen	300
Cl-	300	Nitrate as Nitrogen	300
F-	200	Phosphate as Phosphorus	500
SO₄⁻²	1,500		

This standard is used to prepare the Laboratory Fortified Blank and the Laboratory Fortified Sample Matrix

For use as ion chromatography laboratory fortified blank or fortified sample matrix solution in EPA Method 300.0. Also suitable for use as common anions fortified blank or fortified sample matrix solution in EPA Method 300.1.

QC Standard/Instrument Performance Check [†]			
QCP-QCS-5-125ML		Volume: 125 mL	Matrix: H ₂ O Dilution 1:10 to 1:100
Analyte	µg/mL	Analyte	µg/mL
Br-	50	Nitrite as Nitrogen	15
Cl-	15	Nitrate as Nitrogen	10
F-	10	Phosphate as Phosphorus	25
SO₄⁻²	75	†Manufactured from in-house Second Source concentrates.	

Can be used to prepare the QC Sample or the IPC Solution.

For use as a general ion chromatography quality control standard

300.1 Part B Custom EPA standards are available upon request.

Bromate	
ICBR031	Matrix: H ₂ O
ICBR031-125ML ICBR031-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
BrO ₃ ⁻	1,000

Chlorate	
ICCL031	Matrix: H ₂ O
ICCL031-125ML ICCL031-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
ClO ₃ ⁻	1,000

Bromide	
ICBR1	Matrix: H ₂ O
ICBR1-125ML ICBR1-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Br ⁻	1,000

Dichloroacetate Standard	
ICDCA-S	Matrix: H ₂ O
ICDCA-S-125ML ICDCA-S-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Cl ₂ HC ₂ O ₂ ⁻	500

For use as a surrogate analyte.

Chlorite	
ICCL021	Matrix: H ₂ O
ICCL021-125ML ICCL021-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
ClO ₂ ⁻	1,000

NOTE: Contains less than 10ppm ClO₃⁻.

Custom EPA standards are available upon request.

1,400 µmhos/cm Conductivity at 25°C	
CON1400-25	Matrix: H ₂ O
CON1400-25-125ML CON1400-25-500ML CON1400-25-1L	Volume: 125 mL Volume: 500 mL Volume: 1 L

For the calibration of analytical instruments and validation of analytical methods as appropriate. Preserved with antimicrobial agent.

Perchlorate	
ICCL041	Matrix: H ₂ O
ICCL041-125ML ICCL041-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
ClO ₄ ⁻	1,000

If Atomic Absorption (AA) is your technique of choice, we think you'll appreciate our full line of AA standards.



Contents

Single-Element Standards	95
Modifiers, Buffers & Releasing Agents	98
Instrument Cross-Reference Table.....	34

- ✓ **Five-year shelf life**
- ✓ **Traceable to NIST SRMs**
- ✓ **Produced under ISO 9001**
- ✓ **Assayed by validated procedures**

For the calibration of analytical instruments and validation of analytical methods as appropriate.

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Aluminum, Al	HNO ₃	125 mL	AAAL1-125ML
		500 mL	AAAL1-500ML
Antimony, Sb	HNO ₃ / Tartaric Acid	125 mL	AASB1-125ML
		500 mL	AASB1-500ML
Arsenic, As	HNO ₃	125 mL	AAAS1-125ML
		500 mL	AAAS1-500ML
Barium, Ba	HNO ₃	125 mL	AABA1-125ML
		500 mL	AABA1-500ML
Beryllium, Be	HNO ₃	125 mL	AABE1-125ML
		500 mL	AABE1-500ML
Bismuth, Bi	HNO ₃	125 mL	AABI1-125ML
		500 mL	AABI1-500ML
Boron, B	NH ₄ OH	125 mL	AAB1-125ML
		500 mL	AAB1-500ML
Cadmium, Cd	HNO ₃	125 mL	AACD1-125ML
		500 mL	AACD1-500ML
Calcium, Ca	HNO ₃	125 mL	AACA1-125ML
		500 mL	AACA1-500ML
Cerium, Ce	HNO ₃	125 mL	AACE1-125ML
		500 mL	AACE1-500ML
Cesium, Cs	HNO ₃	125 mL	AACS1-125ML
		500 mL	AACS1-500ML
Chromium, Cr	HNO ₃	125 mL	AACR1-125ML
		500 mL	AACR1-500ML
Cobalt, Co	HNO ₃	125 mL	AAC01-125ML
		500 mL	AAC01-500ML
Copper, Cu	HNO ₃	125 mL	AACU1-125ML
		500 mL	AACU1-500ML
Dysprosium, Dy	HNO ₃	125 mL	AADY1-125ML
		500 mL	AADY1-500ML
Erbium, Er	HNO ₃	125 mL	AAER1-125ML
		500 mL	AAER1-500ML
Europium, Eu	HNO ₃	125 mL	AAEU1-125ML
		500 mL	AAEU1-500ML
Gadolinium, Gd	HNO ₃	125 mL	AAGD1-125ML
		500 mL	AAGD1-500ML
Gallium, Ga	HNO ₃	125 mL	AAGA1-125ML
		500 mL	AAGA1-500ML
Germanium, Ge	HNO ₃ / HF	125 mL	AAGE1-125ML
		500 mL	AAGE1-500ML
Gold, Au	HCl	125 mL	AAAU1-125ML
		500 mL	AAAU1-500ML
Hafnium, Hf	HNO ₃ / HF	125 mL	AAHF1-125ML
		500 mL	AAHF1-500ML
Holmium, Ho	HNO ₃	125 mL	AAHO1-125ML
		500 mL	AAHO1-500ML
Indium, In	HNO ₃	125 mL	AAIN1-125ML
		500 mL	AAIN1-500ML

1,000 µg/mL Standards

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Iridium, Ir	HCl	125 mL	AAIR1-125ML
		500 mL	AAIR1-500ML
Iron, Fe	HNO ₃	125 mL	AAFE1-125ML
		500 mL	AAFE1-500ML
Lanthanum, La	HNO ₃	125 mL	AALA1-125ML
		500 mL	AALA1-500ML
Lead, Pb	HNO ₃	125 mL	AAPB1-125ML
		500 mL	AAPB1-500ML
Lithium, Li	HNO ₃	125 mL	AALI1-125ML
		500 mL	AALI1-500ML
Lutetium, Lu	HNO ₃	125 mL	AALU1-125ML
		500 mL	AALU1-500ML
Magnesium, Mg	HNO ₃	125 mL	AAMG1-125ML
		500 mL	AAMG1-500ML
Manganese, Mn	HNO ₃	125 mL	AAMN1-125ML
		500 mL	AAMN1-500ML
Mercury, Hg	HNO ₃	125 mL	AAHG1-125ML
		500 mL	AAHG1-500ML
Molybdenum, Mo	NH ₄ OH	125 mL	AAMO1-125ML
		500 mL	AAMO1-525ML
Neodymium, Nd	HNO ₃	125 mL	AAND1-125ML
		500 mL	AAND1-500ML
Nickel, Ni	HNO ₃	125 mL	AANI1-125ML
		500 mL	AANI1-500ML
Niobium, Nb	HNO ₃ / HF	125 mL	AANB1-125ML
		500 mL	AANB1-500ML
Palladium, Pd	HCl	125 mL	AAPD1-125ML
		500 mL	AAPD1-500ML
Phosphorus, P	H ₂ O	125 mL	AAP1-125ML
		500 mL	AAP1-500ML
Platinum, Pt	HCl	125 mL	AAPT1-125ML
		500 mL	AAPT1-500ML
Potassium, K	HNO ₃	125 mL	AAK1-125ML
		500 mL	AAK1-500ML
Praseodymium, Pr	HNO ₃	125 mL	AAPR1-125ML
		500 mL	AAPR1-500ML
Rhenium, Re	HNO ₃	125 mL	AARE1-125ML
		500 mL	AARE1-500ML
Rhodium, Rh	HCl	125 mL	AARH1-125ML
		500 mL	AARH1-500ML
Rubidium, Rb	HNO ₃	125 mL	AARB1-125ML
		500 mL	AARB1-500ML
Ruthenium, Ru	HCl	125 mL	AARU1-125ML
		500 mL	AARU1-500ML
Samarium, Sm	HNO ₃	125 mL	AASM1-125ML
		500 mL	AASM1-500ML
Scandium, Sc	HNO ₃	125 mL	AASC1-125ML
		500 mL	AASC1-500ML

1,000 µg/mL

ANALYTE	MATRIX	VOLUME	CATALOG #
Selenium, Se	HNO ₃	125 mL	AASE1-125ML
		500 mL	AASE1-500ML
Silicon, Si	HNO ₃ / HF	125 mL	AASI1-125ML
		500 mL	AASI1-500ML
Silver, Ag	HNO ₃	125 mL	AAAG1-125ML
		500 mL	AAAG1-500ML
Sodium, Na	HNO ₃	125 mL	AANA1-125ML
		500 mL	AANA1-500ML
Strontium, Sr	HNO ₃	125 mL	AASR1-125ML
		500 mL	AASR1-500ML
Sulfur, S	H ₂ O	125 mL	AAS1-125ML
		500 mL	AAS1-500ML
Tantalum, Ta	HNO ₃ / HF	125 mL	AATA1-125ML
		500 mL	AATA1-500ML
Tellurium, Te	HCl	125 mL	AATE1-125ML
		500 mL	AATE1-500ML
Terbium, Tb	HNO ₃	125 mL	AATB1-125ML
		500 mL	AATB1-500ML
Thallium, Tl	HNO ₃	125 mL	AATL1-125ML
		500 mL	AATL1-500ML
Thorium, Th	HNO ₃	125 mL	AATH1-125ML
		500 mL	AATH1-500ML
Thulium, Tm	HNO ₃	125 mL	AATM1-125ML
		500 mL	AATM1-500ML
Tin, Sn	HNO ₃ / HF	125 mL	AASN1-125ML
		500 mL	AASN1-500ML
Titanium, Ti	HNO ₃ / HF	125 mL	AATI1-125ML
		500 mL	AATI1-500ML
Tungsten, W	HNO ₃ / HF	125 mL	AAW1-125ML
		500 mL	AAW1-500ML
Uranium, U	HNO ₃	125 mL	AAU1-125ML
		500 mL	AAU1-500ML
Vanadium, V	HNO ₃	125 mL	AAV1-125ML
		500 mL	AAV1-500ML
Ytterbium, Yb	HNO ₃	125 mL	AAYB1-125ML
		500 mL	AAYB1-500ML
Yttrium, Y	HNO ₃	125 mL	AAY1-125ML
		500 mL	AAY1-500ML
Zinc, Zn	HNO ₃	125 mL	AAZN1-125ML
		500 mL	AAZN1-500ML
Zirconium, Zr	HF	125 mL	AAZR1-125ML
		500 mL	AAZR1-500ML

MODIFIERS, BUFFERS & RELEASING AGENTS

Modifiers, Buffers & Releasing Agents

Custom modifiers, buffers and releasing agents are available upon request.

1% Lanthanum Releasing Agent*	
LACB1	Matrix: HCl
LACB1-500ML	Volume: 500 mL
Analyte	µg/mL
La	10,000

Used as a releasing agent (primarily for Ca in the presence of phosphate).

2% Lithium Ionization Buffer*	
LINB2	Matrix: HNO ₃
LINB2-125ML	Volume: 125 mL
Analyte	µg/mL
Li	20,000

Supplies an excess of electrons to plasma/ flame to minimize impact of ionization interferences.

1% Magnesium Nitrate Modifier*	
MM-MG-10 	Matrix: H ₂ O
MM-MG-10-125ML	Volume: 125 mL
Analyte	µg/mL
Mg(NO ₃) ₂	10,000

Used to change the volatility of the sample to prevent loss of analyte or to facilitate removal of interfering matrix components.

4% Phosphate Modifier*	
MM-P-40	Matrix: H ₂ O
MM-P-40-125ML	Volume: 125 mL
Analyte	µg/mL
PO ₄	40,000

Used to change the volatility of the sample to prevent loss of analyte or to facilitate removal of interfering matrix components.

0.5% Palladium Modifier*	
MM-PD-5	Matrix: HNO ₃
MM-PD-5-125ML	Volume: 125 mL
Analyte	µg/mL
Pd	5,000

Used to change the volatility of the sample to prevent loss of analyte or to facilitate removal of interfering matrix components.

1% Palladium Modifier*	
MM-PD-10	Matrix: HNO ₃
MM-PD-10-125ML	Volume: 125 mL
Analyte	µg/mL
Pd	10,000

Used to change the volatility of the sample to prevent loss of analyte or to facilitate removal of interfering matrix components.

0.3% Palladium / 0.2% Magnesium Nitrate Modifier*			
MM-PDMG-32		Matrix: HNO ₃	
MM-PDMG-32-125ML		Volume: 125 mL	
Analyte	µg/mL	Analyte	µg/mL
Mg(NO ₃) ₂	2,000	Pd	3,000

Used to change the volatility of the sample to prevent loss of analyte or to facilitate removal of interfering matrix components.

*Not to be used as a calibration standard, for analytical reagent use only. ISO 17034 Reference Material; Supplied with Product Information Sheet.

 Agilent/Varian

 Perkin Elmer

Our priority is your total satisfaction. Should you ever have a problem with any standard, Water QC or otherwise, let us know. We'll immediately investigate the problem by testing a retained sample of your solution. If the error is on our end, you'll be offered a full refund or a free replacement – your choice.



- ✓ **Five-year shelf life**
- ✓ **Traceable to NIST SRMs**
- ✓ **Produced under ISO 9001**
- ✓ **Produced under ISO 17025**
- ✓ **Produced under ISO 17034**
- ✓ **Assayed by optimal validated procedures**

Custom water standards for certain products are available upon request.

Bromate	
ICBR031	Matrix: H ₂ O
ICBR031-125ML ICBR031-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
BrO ₃ ⁻	1,000

Chlorite	
ICCL021	Matrix: H ₂ O
ICCL021-125ML ICCL021-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
ClO ₂ ⁻	1,000

NOTE: Contains less than 10ppm ClO₃⁻.

Bromide	
ICBR1	Matrix: H ₂ O
ICBR1-125ML ICBR1-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Br ⁻	1,000

Chromium ⁺⁶ Standard	
QCP-CR6-20ML	Volume: 20 mL Matrix: H ₂ O Dilution 1:100
Analyte	Range*
Cr ⁺⁶	90–900 ug/L

*Parts per billion

Cation Standard	
QCP-CAT-20ML	Volume: 20 mL Matrix: HNO ₃ Dilution 1:100
Analyte	Range
Ca ⁺²	10–100 mg/L
K ⁺	4.0–40 mg/L
Mg ⁺²	4.0–40 mg/L
Na ⁺	10–100 mg/L

Cyanide Standard	
QCP-CN-20ML	Volume: 20 mL Matrix: H ₂ O Dilution 1:200
Analyte	Range
Total Cyanide	0.1–1 mg/L
Free Cyanide	0.05–0.5 mg/L
Cyanide Amenable to Chlorination	0.05–0.5 mg/L

Chlorate	
ICCL031	Matrix: H ₂ O
ICCL031-125ML ICCL031-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
ClO ₃ ⁻	1,000

Hg Standard	
QCP-HG-20ML	Volume: 20 mL Matrix: HNO ₃ Dilution 1:200
Analyte	Range*
Hg	3–30 µg/L

Can be used in conjunction with QCP-TMS and QCP-MTL. *Parts per billion.

Mercury Standard	
MSHG-1PPM	Matrix: HCl
MSHG-1PPM-125ML MSHG-1PPM-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Hg	1

Metals Standard	
QCP-MTL-20ML	Volume: 20 mL Matrix: HNO ₃ Dilution 1:200
Analyte	Range*
Ag	100–1000 µg/L
Al	200–4000 µg/L
As	90–900 µg/L
Ba	100–2500 µg/L
Be	50–500 µg/L
Ca	10–100 mg/L
Cd	100–1000 µg/L
Cr3	100–1000 µg/L
Cu	100–1000 µg/L
Fe	200–4000 µg/L
Mn	200–2000 µg/L
Ni	200–2000 µg/L
Pb	100–1500 µg/L
Se	100–1000 µg/L
Tl	80–800 µg/L
Zn	300–2000 µg/L
Sb	90–900 µg/L

*Parts per billion

Minerals Standard	
QCP-MIN-500ML	Volume: 500 mL Matrix: H ₂ O Dilution: Ready to Use
Analyte	Range
Cl-	35–275 mg/L
F-	0.4–4 mg/L
K+	4–40 mg/L
Nitrate as Nitrogen	2–25 mg/L
Conductivity	200–1,200 µmhos
Alkalinity	25–400 mg/L
Na+	10–100 mg/L
SO₄⁻²	5–125 mg/L

Nitrite Standard	
QCP-NT-20ML	Volume: 20 mL Matrix: H ₂ O Dilution: 1:100
Analyte	Range
Nitrite as Nitrogen	0.4–4.0 mg/L

Oil & Grease Standard	
QCP-OG-A-20ML	Volume: 20 mL Matrix: Acetone Dilution 1:100
Analyte	Range
Oil & Grease	20–200 mg/L

Applicable to gravimetric methods only.

Oil & Grease Standard	
QCP-OG-W-250ML	Volume: 250 mL Matrix: H ₂ O Dilution: See Certificate of Analysis
Analyte	Range
Oil & Grease	20–200 mg/L

Simple Nutrients Standard	
QCP-NUT-1-20ML	Volume: 20 mL Matrix: H ₂ O Dilution 1:200
Analyte	Range
Phosphate as Phosphorus	0.5–5.5 mg/L
Nitrate as Nitrogen	2.0–25 mg/L
Ammonium as Nitrogen	1.0–20 mg/L

Simulated Rainwater Standard	
QCP-RAIN-125ML	Volume: 125 mL Matrix: H ₂ O Dilution: Ready to Use
Analyte	Range
Ca ⁺²	0.05–0.5 mg/L
Cl ⁻	0.1–5 mg/L
F ⁻	0.05–1 mg/L
K ⁺	0.05–1 mg/L
Mg ⁺²	0.05–0.5 mg/L
pH	3.5–4.5 units
Conductivity	20–120 µmhos
Na ⁺	0.2–2 mg/L
NH ₄ ⁺	0.1–1.5 mg/L
NO ₃ ⁻	0.1–10 mg/L
SO ₄ ⁻²	1–12 mg/L

Solids Standard	
QCP-SLD-450ML	Volume: 450 mL Matrix: H ₂ O Dilution: Ready to Use
Analyte	Range
Total Solids (total residue)	140–800 mg/L
Suspended Solids (nonfilterable residue)	20–100 mg/L
Dissolved Solids (filterable residue)	140–800 mg/L

Total Residual Chlorine Standard	
QCP-TRC-10ML	Volume: 10 mL Matrix: H ₂ O Dilution 1:200
Analyte	Range
Total Residual Chlorine	0.5–3.0 mg/L

Trace Metals Standard	
QCP-TMS-20ML	Volume: 20 mL Matrix: HNO ₃ Dilution 1:100
Analyte	Range*
Ag	100–1,000 µg/L
Al	200–4,000 µg/L
As	90–900 µg/L
B	800–2,000 µg/L
Ba	100–2,500 µg/L
Be	50–900 µg/L
Cd	100–1,000 µg/L
Co	100–1,000 µg/L
Cr	100–1,000 µg/L
Cu	100–1,000 µg/L
Fe	200–4,000 µg/L
Mn	200–2,000 µg/L
Mo	60–600 µg/L
Ni	200–2,000 µg/L
Pb	100–1,500 µg/L
Sb	90–900 µg/L
Se	100–1,000 µg/L
Sr	50–500 µg/L
Tl	80–800 µg/L
V	50–2,000 µg/L
Zn	300–2,000 µg/L

*Parts per billion

Turbidity Standard	
QCP-TURB-20ML	Volume: 20 mL Matrix: H ₂ O Dilution: 1:100
Analyte	Range
Turbidity	2–30 NTU

Water Hardness Standard	
QCP-WH-500ML	Volume: 500 mL Matrix: H ₂ O Dilution: Ready to Use
Analyte	Range
Ca	25–250 mg/L
Mg	4.0–40 mg/L
Total Hardness	40–415 mg/L

1,000 µg/mL Total Cyanide	
CN-1000-25-20ML	Volume: 20 mL Matrix: H ₂ O
Analyte	µg/mL
CN-	1,000

For the determination of cyanide in aqueous samples.

TOC Standards

ANALYTE	MATRIX	STARTING MATERIAL	µg/ML	VOLUME	CATALOG #
Total Organic Carbon, TOC	H ₂ O	KHP	1,000	125 mL	TOCKHP1-125ML
				500 mL	TOCKHP1-500ML

Each pH standard is compatible with your instrumentation and meets all requirements for calibration by a true Certified Reference Material. Each standard is traceable to a NIST SRM and is engineered for long-term stability. Manufactured under our ISO 17034 accreditation, each pH standard comes with a CoA and a temperature chart for your convenience. Each product is packaged in our TCT technology with a multi-year shelf life and a one-year* expiration date from opening. All product labels and SDS are GHS-compliant.



Contents

Wet Chemical Standards

Conductivity Standards	108
pH Standards and Colored pH Standards.....	109
Cyanide Standards	110
pH Buffers Specially Formulated for USP <791>	111

Sample Preparation

Dissolution Reagents & Stabilizers	112
Certified Titrants.....	113
Reagents	114
Need a Custom CRM?	16

At times, Wet Chemistry involves some difficult and unusual techniques. If you find yourself in a bind, give us a call. One of our experts will be happy to assist you. Plus, we offer analytical advice and in-depth technical guides on our website, inorganicventures.com.

- ✓ **Five-year shelf life**
- ✓ **Traceable to NIST SRMs**
- ✓ **Produced under ISO 9001**
- ✓ **Produced under ISO 17025**
- ✓ **Produced under ISO 17034**
- ✓ **Assayed by optimal validated procedures**

** For most products.*

For the calibration of analytical instruments and validation of analytical methods as appropriate. Preserved with antimicrobial agent. Custom conductivity standards are available upon request.

2 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON2-25-125ML	Volume: 125 mL
CON2-25-500ML	Volume: 500 mL

5 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON5-25-125ML	Volume: 125 mL
CON5-25-500ML	Volume: 500 mL

10 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON10-25-125ML	Volume: 125 mL
CON10-25-500ML	Volume: 500 mL

84 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON84-25-125ML	Volume: 125 mL
CON84-25-500ML	Volume: 500 mL
CON84-25-1L	Volume: 1 L

100 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON100-25-125ML	Volume: 125 mL
CON100-25-500ML	Volume: 500 mL
CON100-25-1L	Volume: 1 L

147 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON147-25-125ML	Volume: 125 mL
CON147-25-500ML	Volume: 500 mL
CON147-25-1L	Volume: 1 L

500 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON500-25-125ML	Volume: 125 mL
CON500-25-500ML	Volume: 500 mL
CON500-25-1L	Volume: 1 L

1,000 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON1000-25-125ML	Volume: 125 mL
CON1000-25-500ML	Volume: 500 mL
CON1000-25-1L	Volume: 1 L

1,200 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON1200-25-500ML	Volume: 500 mL
CON1200-25-1L	Volume: 1 L

1,400 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON1400-25-125ML	Volume: 125 mL
CON1400-25-500ML	Volume: 500 mL
CON1400-25-1L	Volume: 1 L

1,413 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON1413-25-125ML	Volume: 125 mL
CON1413-25-500ML	Volume: 500 mL
CON1413-25-1L	Volume: 1 L

1,430 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON1430-25-125ML	Volume: 125 mL
CON1430-25-500ML	Volume: 500 mL
CON1430-25-1L	Volume: 1 L

10,000 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON10000-25-125ML	Volume: 125 mL
CON10000-25-500ML	Volume: 500 mL
CON10000-25-1L	Volume: 1 L

100,000 $\mu\text{mhos/cm}$ Conductivity at 25°C	
Matrix: H ₂ O	
CON100000-25-125ML	Volume: 125 mL
CON100000-25-500ML	Volume: 500 mL
CON100000-25-1L	Volume: 1 L

For the calibration of analytical instruments and validation of analytical methods as appropriate. Preserved with antimicrobial agent. Custom pH standards are available upon request.

pH 1.68	
Potassium tetroxalate	
PH-1.68-250ML	Volume: 250 mL
PH-1.68-500ML	Volume: 500 mL
PH-1.68-1L	Volume: 1 L
PH-1.68-4L	Volume: 4 L

pH 2	
Potassium chloride and hydrochloric acid	
PH-2-250ML	Volume: 250 mL
PH-2-500ML	Volume: 500 mL
PH-2-1L	Volume: 1 L
PH-2-4L	Volume: 4 L
PH-2-10L	Volume: 10 L

pH 3	
Potassium acid phthalate and hydrochloric acid	
PH-3-250ML	Volume: 250 mL
PH-3-500ML	Volume: 500 mL
PH-3-1L	Volume: 1 L
PH-3-4L	Volume: 4 L
PH-3-10L	Volume: 10 L

pH 4	
Potassium hydrogen phthalate	
PH-4-250ML	Volume: 250 mL
PH-4-500ML	Volume: 500 mL
PH-4-1L	Volume: 1 L
PH-4-4L	Volume: 4 L
PH-4-10L	Volume: 10 L

pH 4 RED	
Potassium hydrogen phthalate	
PHRED-4-250ML	Volume: 250 mL
PHRED-4-500ML	Volume: 500 mL
PHRED-4-1L	Volume: 1 L
PHRED-4-4L	Volume: 4 L
PHRED-4-10L	Volume: 10 L

pH 4.01	
Potassium hydrogen phthalate	
PH-4.01-500ML	Volume: 500 mL

pH 5	
Potassium acid phthalate and sodium hydroxide	
PH-5-250ML	Volume: 250 mL
PH-5-500ML	Volume: 500 mL
PH-5-1L	Volume: 1 L
PH-5-4L	Volume: 4 L

pH 6	
Monobasic potassium phosphate and sodium hydroxide	
PH-6-250ML	Volume: 250 mL
PH-6-500ML	Volume: 500 mL
PH-6-1L	Volume: 1 L
PH-6-4L	Volume: 4 L
PH-6-10L	Volume: 10 L

pH 6.86	
Potassium phosphate and dibasic sodium phosphate	
PH-6.86-250ML	Volume: 250 mL
PH-6.86-500ML	Volume: 500 mL
PH-6.86-1L	Volume: 1 L
PH-6.86-4L	Volume: 4 L
PH-6.86-10L	Volume: 10 L

pH 7	
Monobasic potassium phosphate and sodium hydroxide	
PH-7-250ML	Volume: 250 mL
PH-7-500ML	Volume: 500 mL
PH-7-1L	Volume: 1 L
PH-7-4L	Volume: 4 L
PH-7-10L	Volume: 10 L

pH 7 YELLOW

Monobasic potassium phosphate and sodium hydroxide

PHYELLOW-7-250ML	Volume: 250 mL
PHYELLOW-7-500ML	Volume: 500 mL
PHYELLOW-7-1L	Volume: 1 L
PHYELLOW-7-4L	Volume: 4 L
PHYELLOW-7-10L	Volume: 10 L

pH 8

Monobasic potassium phosphate and sodium hydroxide

PH-8-250ML	Volume: 250 mL
PH-8-500ML	Volume: 500 mL
PH-8-1L	Volume: 1 L
PH-8-4L	Volume: 4 L
PH-8-10L	Volume: 10 L

pH 9

Boric acid, potassium chloride and sodium hydroxide

PH-9-250ML	Volume: 250 mL
PH-9-500ML	Volume: 500 mL
PH-9-1L	Volume: 1 L
PH-9-4L	Volume: 4 L
PH-9-10L	Volume: 10 L

pH 9.18

Sodium borate decahydrate

PH-9.18-250ML	Volume: 250 mL
PH-9.18-500ML	Volume: 500 mL
PH-9.18-1L	Volume: 1 L

pH 10

Sodium bicarbonate and sodium carbonate

PH-10-250ML	Volume: 250 mL
PH-10-500ML	Volume: 500 mL
PH-10-1L	Volume: 1 L
PH-10-4L	Volume: 4 L
PH-10-10L	Volume: 10 L

pH 10.01

Sodium bicarbonate and sodium carbonate

PH-10.01-500ML	Volume: 500 mL
-----------------------	----------------

pH 10 BLUE

Sodium bicarbonate and sodium carbonate

PHBLUE-10-250ML	Volume: 250 mL
PHBLUE-10-500ML	Volume: 500 mL
PHBLUE-10-1L	Volume: 1 L
PHBLUE-10-4L	Volume: 4 L
PHBLUE-10-10L	Volume: 10 L

pH 11

Dibasic sodium phosphate and sodium hydroxide

PH-11-250ML	Volume: 250 mL
PH-11-500ML	Volume: 500 mL
PH-11-1L	Volume: 1 L
PH-11-4L	Volume: 4 L
PH-11-10L	Volume: 10 L

pH 12

Potassium chloride and sodium hydroxide

PH-12-250ML	Volume: 250 mL
PH-12-500ML	Volume: 500 mL
PH-12-1L	Volume: 1 L
PH-12-4L	Volume: 4 L
PH-12-10L	Volume: 10 L

pH 12.45*

Sodium hydroxide and potassium chloride

PH-12.45-500ML	Volume: 500 mL
-----------------------	----------------

pH 12.47

Sodium hydroxide and potassium chloride

PH-12.47-250ML	Volume: 250 mL
PH-12.47-500ML	Volume: 500 mL
PH-12.47-1L	Volume: 1 L

Cyanide Standards

Custom cyanide standards are available upon request.

1,000 ug/mL Total Cyanide

CN-1000-25-20ML	Volume: 20 mL Matrix: H ₂ O
Analyte	µg/mL
CN-	1,000



Leave buffer preparation
to the experts.
**Your pH meter
will thank you.**

**pH BUFFERS SPECIALLY
FORMULATED FOR USP <791>**

For analysts in the pharmaceutical industry tasked with maintaining proper pH meter calibration in accordance with USP <791>, Inorganic Ventures' CRMs deliver confidence, control, and support.

Our pH buffers are NIST-traceable and manufactured and tested under ISO 17034 and ISO 17025 guidelines. Each standard is accompanied by a detailed Certificate of Analysis (CoA), displaying certified values for multiple temperatures. In addition, these solutions are formulated to meet USP <791> specifications* and are ready to use right out of the bottle.

pH Buffer Features:

- **NIST-traceable standards, certified within 5% of the nominal values and associated uncertainties of no more than 0.05 pH units.**
- **Packaged in Transpiration Control Technology (TCT) which guarantees scientific integrity for 5-years from the date of manufacture.**
- **Ready to use with no preparation required.**

**WE OFFER A COMPLETE LINE OF PH CALIBRATION
STANDARDS, SUITABLE FOR CALIBRATION AS
SPECIFIED IN USP <791>.**

PRODUCT OFFERINGS
pH-1.68
pH-4.01
pH-6.86
pH-9.18
pH-10.01
pH-12.45*

**To ensure solution stability, pH 12.45 is formulated using KCl/NaOH.*



**Do you need a specific pH value not listed
above or require a specific formulation?**

We have the scientific expertise to refine formulations and address your specific needs. Simply contact IV for a custom quote or to learn about our complete line of pH buffer stock standards.



CERTIFICATE# 0883.02

Dissolution Reagents & Stabilizers

Dissolution Reagents are designed for the preparation and measurement of samples containing silica mixed with fluoride insoluble elements, including zeolites, alumina and/or silica based catalysts, sand, limestone, coal fly ash and talc.

These products permit the simple dissolution of silicates without additional equipment, and are designed for ICP measurement of Si and other elements using traditional glass-based sample introduction systems.

The following products are intended to be used together; UA-1 for sample dissolution and UNS-1 for neutralization and stabilization. Please see the article titled Elemental Analysis of Zeolites on our website for more information. Custom dissolution reagents for specific applications are available upon request, and please contact us for more information.


Acid Dissolution Reagent	
UA-1-500ML*	Volume: 500 mL
Recommended for the dissolution of aluminosilicates, such as zeolites.	

Stabilizing Reagent	
UNS-1-500ML*	Volume: 500 mL
Designed for use with UA-1.	

*UA-1-500ML and UNS-1-500ML are custom solutions, made to order with 10-day turnaround time. Please request a custom quote.

Don't see exactly what you're looking for?

Give us a call. Custom reference materials are our specialty.



These Certified Titrants are traceable to NIST and accredited to ISO 17034. Custom certified titrants are available upon request.

0.05M EDTA	
0.05M-EDTA-500ML	Matrix: H ₂ O Volume: 500 mL
0.05M EDTA, 500mL	

0.5M EDTA	
0.5M-EDTA-500ML	Matrix: H ₂ O Volume: 500 mL
0.5M EDTA, 500mL	

0.1M Hydrochloric Acid	
0.1M-HCL-500ML	Matrix: H ₂ O Volume: 500 mL
0.1M Hydrochloric Acid, 500mL	

1.0M Hydrochloric Acid	
1.0M-HCL-500ML	Matrix: H ₂ O Volume: 500 mL
1.0M Hydrochloric Acid, 500mL	

0.1M Nitric Acid	
0.1M-HNO3-500ML	Matrix: H ₂ O Volume: 500 mL
0.1M Nitric Acid, 500mL	

1.0M Nitric Acid	
1.0M-HNO3-500ML	Matrix: H ₂ O Volume: 500 mL
1.0M Nitric Acid, 500mL	

0.1M Perchloric Acid	
0.1M-HCLO4-500ML	Matrix: H ₂ O/0.1M HClO ₄ in Glacial Acetic Acid Volume: 500 mL
0.1M Perchloric Acid, 500mL	

0.1N Silver Nitrate	
0.1N-AGNO3-500ML	Matrix: H ₂ O Volume: 500 mL
0.1N Silver Nitrate, 500mL	

0.1M Sodium Hydroxide	
0.1M-NAOH-500ML	Matrix: H ₂ O Volume: 500 mL
0.1M Sodium Hydroxide, 500mL	

0.1N Sodium Thiosulfate	
0.1N-NA2S2O3-500ML	Matrix: H ₂ O Volume: 500 mL
0.1N Sodium Thiosulfate 500 mL. Prepared and standardized according to USP specifications.	

1M Sodium Hydroxide	
1M-NAOH-500ML	Matrix: H ₂ O Volume: 500 mL
1M Sodium Hydroxide, 500mL	

Blank & Rinse Solutions

Blank & Rinse solutions are prepared using double-distilled reagents and 18 megohm (MΩ) deionized water. They come packaged in ultra-clean LDPE bottles and are ready to use. Custom solutions are available upon request.

2% (v/v) Nitric Acid Rinse	
CLP-MS-RINSE Ultra Pure	Matrix: HNO ₃ Dilution: Ready to Use
CLP-MS-RINSE-125ML	Volume: 125 mL
CLP-MS-RINSE-500ML	Volume: 500 mL

Deionized Blank	
IV-DI-BLANK 	Matrix: H ₂ O
IV-DI-BLANK-500ML	Volume: 500 mL
IV-DI-BLANK-1L	Volume: 1 L

For use with ICP-MS. Designed for ILM05.2 and ILM05.3.

5% (v/v) Nitric Acid Blank	
IV-ACID-BLANK Ultra Pure 	Matrix: HNO ₃
IV-ACID-BLANK-500ML	Volume: 500 mL
IV-ACID-BLANK-1L	Volume: 1 L

Platinum Cobalt 500 CU APHA/Hazen Color Standard	
PT-CO 	Matrix: HCl
PT-CO-500ML	Volume: 500 mL
Analyte	µg/mL
Pt	500
Co	250

Combined Color Reagent for Chloride Analysis by Discrete Analyzer	
DA-CL-COLOR	Matrix: 15% Methanol
DA-CL-COLOR-500ML	Volume: 500 mL

 Agilent/Varian

 HACH

 Perkin Elmer

AA Standards	
Modifiers, Buffers, & Releasing Agents	101
Single-Element Standards	98–100
Anion Standards	
Multi-Ion Standards	92–93
Single-Ion Standards	89–91
Buffers	
AA	101
ICP & ICP-MS	55
Calibration Blank & Rinse Solutions	69, 112
Cannabis Standards	58–60
Cation Standards	
Multi-Ion Standards	92–93
Single-Ion Standards	91
Conductivity Standards	108
Continuing and Initial Calibration	
Verification Standards (CICV)	62, 64, 67
Contract Required Detection Limit	
Standards (CRDL)	65
Contract Required Quantitation Limit	
Standards (CRQL)	67
Cyanide Standards	
Anion Standards	89
ICP & ICP-MS	40
Water Standards	103
Wet Chemistry	110
Dichloroacetate Standard	94, 96
Dissolution Reagents & Stabilizers	112
Eluents for Anions and Cations	93–94
EPA Methods	
200.7	69–78
200.8	79–81
300.0	94–96
300.1	94–96
314.0	96
1311	82
6020	83–87
ILM03.0	62–63
ILM04.0	64–65
ILM05.2	66–69
ILM05.3	66–69
Toxicity Characteristic Leachate Procedure (TCLP)	82
ICP, ICP-MS, and ICP-OES/AES Standards	
Multi-Element	41–55
Single-Element	26–40
Instrument Cross-Reference Table	18–24
Instrument Performance Check	74–75, 95
Interference Check Standards (ICS)	46, 63, 65, 67, 70, 83–86
Internal Standards	46–47, 51, 55, 68, 80, 83, 85–86
Ion Chromatography Standards	
Multi-Ion Standards	92–93
Single-Ion Standards	89–91
Ionization Buffers	55, 101
Laboratory Fortifying Stock Solutions	76
Matrix Modifiers	101
Mercury Preservation Solution	80
Mercury Standard	80, 105
Modifiers	101
pH Standards	
Calibration	109–111
Wet Chemistry	109–111
Quality Control Standards (QC)	
ICP & ICP-MS	71–72, 77–78, 81
Ion Chromatography	95
Reagents	114
Releasing Agents	101
Sample Preparation	
Dissolution Reagents	112
Blank & Rinse Solutions	114
Single-Ion Standards	89–91
Soil Spike Standards	63, 65, 68, 83, 85–86
Total Organic Carbon (TOC)	106
Toxicity Characteristic Leachate	
Procedure Standard (TCLP)	82
Tuning Solutions	44, 47, 52–55, 68, 81, 84–85, 87
USP <232> – Elemental Impurities	
Compliance Standards	56–57
Water Standards	103–106
Water Spike Standards	63, 65, 68, 84–85, 87
Wet Chemistry Standards	
Certified Titrants	113
Conductivity Standards	108
Cyanide Standards	110
Dissolution Reagents and Stabilizers	112
pH Standards and Colored pH Standards	109–111
Reagents	114

0.05M-EDTA-500ML.....	113	AAC01-500ML.....	98	AAPR1-500ML.....	99	CARB-100ML.....	94
0.1M-HCL-500ML.....	113	AACR1-125ML.....	98	AAPT1-125ML.....	99	CARB-500ML.....	94
0.1M-HCL04-500ML.....	113	AACR1-500ML.....	98	AAPT1-500ML.....	99	CCS-1-125ML.....	49
0.1M-HNO3-500ML.....	113	AACS1-125ML.....	98	AARB1-125ML.....	99	CCS-1-500ML.....	49
0.1M-NAOH-500ML.....	113	AACS1-500ML.....	98	AARB1-500ML.....	99	CCS-2-125ML.....	49
0.1N-AGNO3-500ML.....	113	AACU1-125ML.....	98	AARE1-125ML.....	99	CCS-4-125ML.....	49
0.1N-NA2S2O3-500ML.....	113	AACU1-500ML.....	98	AARE1-500ML.....	99	CCS-4-500ML.....	49
0.5M-EDTA-500ML.....	113	AACUCN-125ML.....	40	AARH1-125ML.....	99	CCS-5-125ML.....	49
1.0M-HCL-500ML.....	113	AACUCN-500ML.....	40	AARH1-500ML.....	99	CCS-5-500ML.....	49
1.0M-HNO3-500ML.....	113	AADY1-125ML.....	98	AARU1-125ML.....	99	CCS-6-125ML.....	50
1M-NAOH-500ML.....	113	AADY1-500ML.....	98	AARU1-500ML.....	99	CCS-6-500ML.....	50
2007ICS-1-125ML.....	70	AAER1-125ML.....	98	AAS1-125ML.....	100	CG6LI1-125ML.....	32
2007ICS-3-125ML.....	70	AAER1-500ML.....	98	AAS1-500ML.....	100	CG6LI1-30ML.....	32
2007ICS-4-125ML.....	70	AAEU1-125ML.....	98	AASB1-125ML.....	98	CGAG10-125ML.....	39
2008CAL-1-125ML.....	79	AAEU1-500ML.....	98	AASB1-500ML.....	98	CGAG10-500ML.....	39
2008CAL-2-125ML.....	79	AAFE1-125ML.....	99	AASC1-125ML.....	99	CGAG1-125ML.....	34
2008ISS-125ML.....	80	AAFE1-500ML.....	99	AASC1-500ML.....	99	CGAG1-30ML.....	34
2008ISS-500ML.....	80	AAGA1-125ML.....	98	AASE1-125ML.....	100	CGAG1-500ML.....	34
2008TS-125ML.....	68, 81	AAGA1-500ML.....	98	AASE1-500ML.....	100	CGAL10-125ML.....	36
300-CAL-A-125ML.....	94	AAGD1-125ML.....	98	AASI1-125ML.....	100	CGAL10-30ML.....	36
300-CAL-A-500ML.....	94	AAGD1-500ML.....	98	AASI1-500ML.....	100	CGAL10-500ML.....	36
300-LFS-A-125ML.....	95	AAGE1-125ML.....	98	AASM1-125ML.....	99	CGAL1-125ML.....	30
6020CAL-1-125ML.....	83, 84, 86	AAGE1-500ML.....	98	AASM1-500ML.....	99	CGAL1-30ML.....	30
6020ICS-0A-125ML.....	86	AAHF1-125ML.....	98	AASN1-125ML.....	100	CGAL1-500ML.....	30
6020ICS-0A-500ML.....	86	AAHF1-500ML.....	98	AASN1-500ML.....	100	CGALCL1-125ML.....	30
6020ICS-0B-125ML.....	86	AAHG1-125ML.....	99	AASR1-125ML.....	100	CGALCL1-30ML.....	30
6020ICS-8A-125ML.....	83	AAHG1-500ML.....	99	AASR1-500ML.....	100	CGALCL1-500ML.....	30
6020ICS-8A-500ML.....	83	AAHO1-125ML.....	98	AATA1-125ML.....	100	CGAS(3)1-125ML.....	30, 59
6020ICS-9A-125ML.....	84	AAHO1-500ML.....	98	AATA1-500ML.....	100	CGAS(3)1-30ML.....	30, 59
6020ICS-9A-500ML.....	84	AAIN1-125ML.....	98	AATB1-125ML.....	100	CGAS(3)1-500ML.....	30, 59
6020ICS-9B-125ML.....	85	AAIN1-500ML.....	98	AATB1-500ML.....	100	CGAS(5)1-125ML.....	30, 59
6020ISS-125ML.....	68, 83, 85, 86	AAIR1-125ML.....	99	AATE1-125ML.....	100	CGAS(5)1-30ML.....	30, 59
6020ISS-500ML.....	68, 83, 85, 86	AAIR1-500ML.....	99	AATE1-500ML.....	100	CGAS(5)1-500ML.....	30, 59
6020SPK-W-125ML.....	84, 85, 87	AAK1-125ML.....	99	AATH1-125ML.....	100	CGAS10-125ML.....	36, 59
AAAG1-125ML.....	100	AAK1-500ML.....	99	AATH1-500ML.....	100	CGAS10-30ML.....	36, 59
AAAG1-500ML.....	100	AALA1-125ML.....	99	AATI1-125ML.....	100	CGAS10-500ML.....	36, 59
AAAGCN-125ML.....	40	AALA1-500ML.....	99	AATI1-500ML.....	100	CGAS1-125ML.....	30, 59
AAAGCN-500ML.....	40	AALI1-125ML.....	99	AATL1-125ML.....	100	CGAS1-30ML.....	30, 59
AAAL1-125ML.....	98	AALI1-500ML.....	99	AATL1-500ML.....	100	CGAS1-500ML.....	30, 59
AAAL1-500ML.....	98	AALU1-125ML.....	99	AATM1-125ML.....	100	CGAU10-125ML.....	37
AAAS1-125ML.....	98	AALU1-500ML.....	99	AATM1-500ML.....	100	CGAU10-30ML.....	37
AAAS1-500ML.....	98	AAMG1-125ML.....	99	AAU1-125ML.....	100	CGAU10-500ML.....	37
AAAU1-125ML.....	98	AAMG1-500ML.....	99	AAU1-500ML.....	100	CGAU1-125ML.....	31
AAAU1-500ML.....	98	AAMN1-125ML.....	99	AAV1-125ML.....	100	CGAU1-30ML.....	31
AAAUCN-125ML.....	40	AAMN1-500ML.....	99	AAV1-500ML.....	100	CGAU1-500ML.....	31
AAAUCN-500ML.....	40	AAM01-125ML.....	99	AAW1-125ML.....	100	CGAUN1-125ML.....	31, 80
AAB1-125ML.....	98	AAM01-500ML.....	99	AAW1-500ML.....	100	CGAUN1-30ML.....	31, 80
AAB1-500ML.....	98	AANA1-125ML.....	100	AAY1-125ML.....	100	CGAUN1-500ML.....	31, 80
AABA1-125ML.....	98	AANA1-500ML.....	100	AAY1-500ML.....	100	CGB10-125ML.....	36
AABA1-500ML.....	98	AANB1-125ML.....	99	AAYB1-125ML.....	100	CGB10-500ML.....	36
AABE1-125ML.....	98	AANB1-500ML.....	99	AAYB1-500ML.....	100	CGB1-125ML.....	30
AABE1-500ML.....	98	AAND1-125ML.....	99	AAZN1-125ML.....	100	CGB1-30ML.....	30
AABI1-125ML.....	98	AAND1-500ML.....	99	AAZN1-500ML.....	100	CGB1-500ML.....	30
AABI1-500ML.....	98	AANI1-125ML.....	99	AAZNCN-125ML.....	40	CGBA10-125ML.....	36
AACA1-125ML.....	98	AANI1-500ML.....	99	AAZNCN-500ML.....	40	CGBA10-500ML.....	36
AACA1-500ML.....	98	AAP1-125ML.....	99	AAZR1-125ML.....	100	CGBA1-125ML.....	30
AACD1-125ML.....	98	AAP1-500ML.....	99	AAZR1-500ML.....	100	CGBA1-30ML.....	30
AACD1-500ML.....	98	AAPB1-125ML.....	99	AGI-TS-1-125ML.....	52	CGBA1-500ML.....	30
AACE1-125ML.....	98	AAPB1-500ML.....	99	AGI-TS-1-500ML.....	52	CGBE10-125ML.....	36
AACE1-500ML.....	98	AAPD1-125ML.....	99	BICARB-100ML.....	93	CGBE10-500ML.....	36
AAC01-125ML.....	98	AAPD1-500ML.....	99	BICARB-500ML.....	93	CGBE1-125ML.....	30
		AAPR1-125ML.....	99			CGBE1-30ML.....	30

BY CATALOG NUMBER

CGBE1-500ML	30	CGER1-30ML	31	CGIR1-125ML	32	CGND10-125ML	38
CGBI10-125ML	36	CGER1-500ML	31	CGIR1-30ML	32	CGND10-30ML	38
CGBI10-30ML	36	CGEU10-125ML	36	CGIR1-500ML	32	CGND10-500ML	38
CGBI10-500ML	36	CGEU10-30ML	36	CGK10-125ML	38	CGND1-125ML	32
CGBI1-125ML	30	CGEU10-500ML	36	CGK10-30ML	38	CGND1-30ML	32
CGBI1-30ML	30	CGEU1-125ML	31	CGK10-500ML	38	CGND1-500ML	32
CGBI1-500ML	30	CGEU1-30ML	31	CGK1-125ML	33	CGNI10-125ML	38
CGC10-125ML	36	CGEU1-500ML	31	CGK1-30ML	33	CGNI10-30ML	38
CGC1-125ML	30	CGFE10-125ML	37	CGK1-500ML	33	CGNI10-500ML	38
CGC1-500ML	30	CGFE10-30ML	37	CGLA10-125ML	37	CGNI1-125ML	32
CGCA10-125ML	36	CGFE10-500ML	37	CGLA10-30ML	37	CGNI1-30ML	32
CGCA10-30ML	36	CGFE1-125ML	32	CGLA10-500ML	37	CGNI1-500ML	32
CGCA10-500ML	36	CGFE1-30ML	32	CGLA1-125ML	32	CGOS1-125ML	33
CGCA1-125ML	30	CGFE1-500ML	32	CGLA1-30ML	32	CGOS1-30ML	33
CGCA1-30ML	30	CGGA10-125ML	37	CGLA1-500ML	32	CGOS1-500ML	33
CGCA1-500ML	30	CGGA10-500ML	37	CGLI10-125ML	37	CGP10-125ML	38
CGCD10-125ML	36, 59	CGGA1-125ML	31	CGLI10-30ML	37	CGP10-30ML	38
CGCD10-500ML	36, 59	CGGA1-30ML	31	CGLI10-500ML	37	CGP10-500ML	38
CGCD1-125ML	30, 59	CGGA1-500ML	31	CGLI1-125ML	32	CGP1-125ML	33
CGCD1-30ML	30, 59	CGGD10-125ML	37	CGLI1-30ML	32	CGP1-30ML	33
CGCD1-500ML	30, 59	CGGD10-30ML	37	CGLI1-500ML	32	CGP1-500ML	33
CGCE10-125ML	36	CGGD10-500ML	37	CGLU10-125ML	37	CGPB10-125ML	37, 59
CGCE10-30ML	36	CGGD1-125ML	31	CGLU10-30ML	37	CGPB10-30ML	37, 59
CGCE10-500ML	36	CGGD1-30ML	31	CGLU10-500ML	37	CGPB10-500ML	37, 59
CGCE1-125ML	30	CGGD1-500ML	31	CGLU1-125ML	32	CGPB1-125ML	32, 59
CGCE1-30ML	30	CGGE10-125ML	37	CGLU1-30ML	32	CGPB1-30ML	32, 59
CGCE1-500ML	30	CGGE10-500ML	37	CGLU1-500ML	32	CGPB1-500ML	32, 59
CGCO10-125ML	36	CGGE1-125ML	31	CGMG10-125ML	37	CGPD10-125ML	38
CGCO10-30ML	36	CGGE1-30ML	31	CGMG10-30ML	37	CGPD10-30ML	38
CGCO10-500ML	36	CGGE1-500ML	31	CGMG10-500ML	37	CGPD10-500ML	38
CGCO1-125ML	31	CGHF10-125ML	37	CGMG1-125ML	32	CGPD1-125ML	33
CGCO1-30ML	31	CGHF10-500ML	37	CGMG1-30ML	32	CGPD1-30ML	33
CGCO1-500ML	31	CGHF1-125ML	31	CGMG1-500ML	32	CGPD1-500ML	33
CGCR(3)10-125ML	36	CGHF1-30ML	31	CGMN10-125ML	37	CGPDN1-125ML	33
CGCR(3)10-30ML	36	CGHF1-500ML	31	CGMN10-30ML	37	CGPDN1-30ML	33
CGCR(3)10-500ML	36	CGHG10-125ML	37, 59	CGMN10-500ML	37	CGPDN1-500ML	33
CGCR(3)1-125ML	31	CGHG10-500ML	37, 59	CGMN1-125ML	32	CGPR10-125ML	38
CGCR(3)1-30ML	31	CGHG1-125ML	32, 59	CGMN1-30ML	32	CGPR10-30ML	38
CGCR(3)1-500ML	31	CGHG1-30ML	32, 59	CGMN1-500ML	32	CGPR10-500ML	38
CGCR(6)1-125ML	31	CGHG1-500ML	32, 59	CGMO10-125ML	37	CGPR1-125ML	33
CGCR(6)1-30ML	31	CGHO10-125ML	37	CGMO10-30ML	37	CGPR1-30ML	33
CGCR(6)1-500ML	31	CGHO10-30ML	37	CGMO10-500ML	37	CGPR1-500ML	33
CGCS10-125ML	36	CGHO10-500ML	37	CGMO1-125ML	32	CGPT10-125ML	38
CGCS10-500ML	36	CGHO1-125ML	31	CGMO1-30ML	32	CGPT10-30ML	38
CGCS1-125ML	31	CGHO1-30ML	31	CGMO1-500ML	32	CGPT10-500ML	38
CGCS1-30ML	31	CGH01-500ML	31	CGMSA10-125ML	39	CGPT1-125ML	33
CGCU10-125ML	36	CGICBR1-125ML	30	CGMSA10-500ML	39	CGPT1-30ML	33
CGCU10-30ML	36	CGICBR1-30ML	30	CGMSA1-125ML	34	CGPT1-500ML	33
CGCU10-500ML	36	CGICBR1-500ML	30	CGMSA1-500ML	34	CGPTN1-125ML	33
CGCU1-125ML	31	CGICCL1-125ML	31	CGNA10-125ML	39	CGPTN1-30ML	33
CGCU1-30ML	31	CGICCL1-500ML	31	CGNA10-30ML	39	CGPTN1-500ML	33
CGCU1-500ML	31	CGIC11-125ML	32	CGNA10-500ML	39	CGPTNO31-125ML	33
CGDY10-125ML	36	CGIC11-30ML	32	CGNA1-125ML	34	CGPTNO31-30ML	33
CGDY10-30ML	36	CGIC11-500ML	32	CGNA1-30ML	34	CGPTNO31-500ML	33
CGDY10-500ML	36	CGIN10-125ML	37	CGNA1-500ML	34	CGRB10-125ML	38
CGDY1-125ML	31	CGIN10-500ML	37	CGNB10-125ML	38	CGRB10-500ML	38
CGDY1-30ML	31	CGIN1-125ML	32	CGNB1-125ML	32	CGRB1-125ML	33
CGDY1-500ML	31	CGIN1-30ML	32	CGNB1-30ML	32	CGRB1-30ML	33
CGER10-125ML	36	CGIN1-500ML	32	CGNB1-500ML	32	CGRE10-125ML	38
CGER10-30ML	36	CGIR10-125ML	37	CGNB20510-125ML	38	CGRE1-125ML	33
CGER10-500ML	36	CGIR10-30ML	37	CGNB2051-125ML	33	CGRE1-30ML	33
CGER1-125ML	31	CGIR10-500ML	37	CGNB2051-500ML	33	CGRE1-500ML	33

CGRH10-125ML	38	CGSN1-125ML	35	CGW1-125ML	35	CON100-25-500ML	108
CGRH10-30ML	38	CGSN1-30ML	35	CGW1-30ML	35	CON10-25-125ML	108
CGRH10-500ML	38	CGSN1-500ML	35	CGW1-500ML	35	CON10-25-500ML	108
CGRH1-125ML	33	CGSNCL1-125ML	35	CGWH201-125ML	35	CON1200-25-125ML	108
CGRH1-30ML	33	CGSNCL1-30ML	35	CGY10-125ML	39	CON1200-25-1L	108
CGRH1-500ML	33	CGSR10-125ML	39	CGY10-30ML	39	CON1200-25-500ML	108
CGRHN1-125ML	33	CGSR10-500ML	39	CGY10-500ML	39	CON1400-25-125ML	108
CGRHN1-30ML	33	CGSR1-125ML	34	CGY1-125ML	35	CON1400-25-1L	108
CGRHN1-500ML	33	CGSR1-30ML	34	CGY1-30ML	35	CON1400-25-500ML	108
CGRU10-125ML	38	CGSR1-500ML	34	CGY1-500ML	35	CON1413-25-125ML	108
CGRU10-30ML	38	CGTA10-125ML	39	CGYB10-125ML	39	CON1413-25-1L	108
CGRU10-500ML	38	CGTA1-125ML	34	CGYB10-30ML	39	CON1413-25-500ML	108
CGRU1-125ML	33	CGTA1-30ML	34	CGYB10-500ML	39	CON1430-25-125ML	108
CGRU1-30ML	33	CGTA1-500ML	34	CGYB1-125ML	35	CON1430-25-1L	108
CGRU1-500ML	33	CGTB10-125ML	39	CGYB1-30ML	35	CON1430-25-500ML	108
CGS10-125ML	39	CGTB10-30ML	39	CGYB1-500ML	35	CON147-25-125ML	108
CGS10-30ML	39	CGTB10-500ML	39	CGZN10-125ML	40	CON147-25-1L	108
CGS10-500ML	39	CGTB1-125ML	34	CGZN10-30ML	40	CON147-25-500ML	108
CGS1-125ML	34	CGTB1-30ML	34	CGZN10-500ML	40	CON2-25-125ML	108
CGS1-30ML	34	CGTB1-500ML	34	CGZN1-125ML	35	CON2-25-500ML	108
CGS1-500ML	34	CGTE10-125ML	39	CGZN1-30ML	35	CON500-25-125ML	108
CGSB10-125ML	36	CGTE10-500ML	39	CGZN1-500ML	35	CON500-25-1L	108
CGSB10-500ML	36	CGTE1-125ML	34	CGZR10-125ML	40	CON500-25-500ML	108
CGSB1-125ML	30, 62, 64, 66, 70	CGTE1-30ML	34	CGZR10-30ML	40	CON5-25-125ML	108
CGSB1-30ML	30	CGTE1-500ML	34	CGZR10-500ML	40	CON5-25-500ML	108
CGSB1-500ML	30	CGTEN1-125ML	34	CGZR1-125ML	35	CON84-25-125ML	108
CGSBF1-125ML	30	CGTEN1-30ML	34	CGZR1-30ML	35	CON84-25-1L	108
CGSC10-125ML	38	CGTEN1-500ML	34	CGZR1-500ML	35	CON84-25-500ML	108
CGSC10-30ML	38	CGTH10-125ML	39	CGZRCL10-125ML	40	CSN-125-500ML	55
CGSC10-500ML	38	CGTH1-125ML	35	CGZRCL10-500ML	40	CSN-125-500ML	55
CGSC1-125ML	34	CGTH1-30ML	35	CIROS-OES-TS-125ML	52	DA-CL-COLOR-500ML	114
CGSC1-30ML	34	CGTH1-500ML	35	CLP-MS-RINSE-125ML	69, 114	ELUENT1817-500ML	93, 94
CGSC1-500ML	34	CGT110-125ML	39	CLP-MS-RINSE-500ML	69, 114	ELUENT3510-500ML	93
CGSE(4)1-125ML	34	CGT110-30ML	39	CLP-MS-SPK-125ML	68	ELUENT4514-500ML	93
CGSE(4)1-30ML	34	CGT110-500ML	39	CLPP-CAL-1-125ML	62, 64, 66	GENESIS-ICAL-125ML	52
CGSE(4)1-500ML	34	CGT11-125ML	35	CLPP-CAL-3-125ML	62, 64, 66	GENESIS-ICAL-500ML	52
CGSE(6)1-125ML	34	CGT11-30ML	35	CLPP-ICS-A-125ML	63, 65, 67	ICADP1-125ML	89
CGSE(6)1-30ML	34	CGT11-500ML	35	CLPP-ICS-A-500ML	63, 65, 67	ICBA1-125ML	91
CGSE10-125ML	38	CGTL10-125ML	39	CLPP-ICS-B-125ML	63	ICBEN1-125ML	89
CGSE10-30ML	38	CGTL10-500ML	39	CLPP-ICS-B4-125ML	65, 67	ICBR1-125ML	89, 96, 103
CGSE10-500ML	38	CGTL1-125ML	35	CLPP-SPK-1-125ML	63, 65, 68	ICBR1-500ML	89, 96, 103
CGSI10-125ML	38	CGTL1-30ML	35	CLPP-SPK-2-125ML	63, 69, 73	ICBRO31-125ML	89, 96, 103
CGSI10-30ML	38	CGTL1-500ML	35	CMS-1-125ML	48	ICBTR1-125ML	89
CGSI10-500ML	38	CGTM10-125ML	39	CMS-2-125ML	48	ICCA1-125ML	91
CGSI1-125ML	34	CGTM10-30ML	39	CMS-2-500ML	48	ICCA1-500ML	91
CGSI1-30ML	34	CGTM1-125ML	35	CMS-3-125ML	48	ICCL1-125ML	89
CGSI1-500ML	34	CGTM1-30ML	35	CMS-4-125ML	48	ICCL1-500ML	89
CGSINA1-125ML	34	CGTM1-500ML	35	CMS-4-500ML	48	ICCL021-125ML	89, 96, 103
CGSINA1-500ML	34	CGU10-125ML	39	CMS-5-125ML	48	ICCL021-500ML	89, 96, 103
CGSIO1-125ML	34	CGU10-30ML	39	CMS-5-500ML	48	ICCL031-125ML	89, 96, 103
CGSIO1-30ML	34	CGU10-500ML	39	CN-1000-25-20ML	89, 106, 110	ICCL031-500ML	89, 96, 103
CGSIO1-500ML	34	CGU1-125ML	35	CON10000-25-125ML	108	ICCL041-125ML	90, 96
CGSIONA1-125ML	34	CGU1-30ML	35	CON100000-25-125ML	108	ICCL041-500ML	90, 96
CGSIONA1-500ML	34	CGU1-500ML	35	CON100000-25-1L	108		
CGSM10-125ML	38	CGU10-125ML	39	CON100000-25-500ML	108		
CGSM10-30ML	38	CGV10-125ML	39	CON10000-25-125ML	108		
CGSM10-500ML	38	CGV10-30ML	39	CON10000-25-1L	108		
CGSM1-125ML	33	CGV10-500ML	39	CON10000-25-500ML	108		
CGSM1-30ML	33	CGV1-125ML	35	CON1000-25-125ML	108		
CGSM1-500ML	33	CGV1-30ML	35	CON1000-25-1L	108		
CGSN10-125ML	39	CGV1-500ML	35	CON1000-25-500ML	108		
CGSN10-30ML	39	CGW10-125ML	39	CON100-25-125ML	108		
CGSN10-500ML	39	CGW10-500ML	39	CON100-25-1L	108		

BY CATALOG NUMBER

ICCO31-125ML.....	89	ICSCN1-125ML.....	90	IV-STOCK-33-125ML.....	45	MSAS-10PPM-125ML.....	26, 58
ICCO31-500ML.....	89	IC-SCS1-125ML.....	92	IV-STOCK-33-500ML.....	45	MSAU-100PPM-125ML.....	28
ICCRO41-125ML.....	89	ICSO410-125ML.....	90	IV-STOCK-34-125ML.....	45	MSAU-100PPM-500ML.....	28
ICCS1-125ML.....	91	ICSO410-500ML.....	90	IV-STOCK-34-500ML.....	45	MSAU-10PPM-125ML.....	26
ICDCA-S-125ML.....	94, 96	ICSO41-125ML.....	90	IV-STOCK-35-125ML.....	45	MSAU-10PPM-500ML.....	26
ICDCA-S-500ML.....	94, 96	ICSO41-500ML.....	90	IV-STOCK-36-125ML.....	46	MSB-100PPM-125ML.....	28
ICDEA1-125ML.....	91	ICSR1-125ML.....	91	IV-STOCK-36-500ML.....	46	MSB-10PPM-125ML.....	26
ICDMA1-125ML.....	91	ICTEA1-125ML.....	91	IV-STOCK-38-125ML.....	56	MSBA-100PPM-125ML.....	28
ICF1-125ML.....	89	ICTMA1-125ML.....	91	IV-STOCK-40-125ML.....	56	MSBA-10PPM-125ML.....	26
ICF1-500ML.....	89	ICTMAH1-125ML.....	91	IV-STOCK-41-125ML.....	56	MSBE-100PPM-125ML.....	28
IC-FAS-1A-125ML.....	92	ICTRTR1-125ML.....	90	IV-STOCK-50-125ML.....	46	MSBE-10PPM-125ML.....	26
IC-FAS-1A-500ML.....	92	IV-19-125ML.....	72, 78	IV-STOCK-51-125ML.....	46	MSBI-100PPM-125ML.....	28
ICGLY1-125ML.....	89	IV-21-125ML.....	72, 78	IV-STOCK-52-125ML.....	46	MSBI-10PPM-125ML.....	26
ICGTR1-125ML.....	89	IV-26-125ML.....	72, 78	IV-STOCK-53-125ML.....	46	MSCA-100PPM-125ML.....	28
ICHCO1-125ML.....	89	IV-26-500ML.....	72, 78	IV-STOCK-54-125ML.....	46	MSCA-100PPM-500ML.....	28
ICHCO1-500ML.....	89	IV-28-125ML.....	72, 78	IV-STOCK-55-125ML.....	47	MSCA-10PPM-125ML.....	26
ICII-125ML.....	89	IV-28-500ML.....	72, 78	IV-STOCK-56-125ML.....	47	MSCA-10PPM-500ML.....	26
ICII-500ML.....	89	IV-7-125ML.....	71, 77	IV-STOCK-57-125ML.....	47	MSCD-100PPM-125ML.....	28, 58
ICK1-125ML.....	91	IV-ACID-BLANK-1L.....	114	IV-STOCK-58-125ML.....	47	MSCD-10PPM-125ML.....	26, 58
ICK1-500ML.....	91	IV-ACID-BLANK-500ML.....	114	IV-STOCK-59-125ML.....	92	MSCE-100PPM-125ML.....	28
ICKHP1-125ML.....	90	IV-DI-BLANK-1L.....	114	IV-STOCK-59-500ML.....	92	MSCE-10PPM-125ML.....	26
ICLCT1-125ML.....	89	IV-DI-BLANK-500ML.....	114	IV-STOCK-60-125ML.....	56	MSCO-100PPM-125ML.....	28
ICLI1-125ML.....	91	IV-ICPMS-71A-125ML.....	50	IV-STOCK-61-125ML.....	92	MSCO-10PPM-125ML.....	26
ICMEA1-125ML.....	91	IV-ICPMS-71A-500ML.....	50	IV-STOCK-61-500ML.....	92	MSCR(3)-100PPM-125ML.....	28
ICMEA1-500ML.....	91	IV-ICPMS-71B-125ML.....	50	IV-STOCK-62-125ML.....	92	MSCR(3)-10PPM-125ML.....	26
ICMG1-125ML.....	91	IV-ICPMS-71B-500ML.....	50	IV-STOCK-63-125ML.....	93	MSCR(6)-100PPM-125ML.....	28
ICMG1-500ML.....	91	IV-ICPMS-71C-125ML.....	51	IV-STOCK-64-125ML.....	93	MSCR(6)-10PPM-125ML.....	26
ICMLA1-125ML.....	89	IV-ICPMS-71C-500ML.....	51	IV-STOCK-64-500ML.....	93	MSCS-100PPM-125ML.....	28
ICMLE1-125ML.....	89	IV-ICPMS-71D-125ML.....	51	IV-STOCK-65-125ML.....	56, 60	MSCS-10PPM-125ML.....	26
ICML01-125ML.....	89	IV-ICPMS-71D-500ML.....	51	IV-STOCK-66-125ML.....	56	MSCU-100PPM-125ML.....	28
ICMMA1-125ML.....	91	IV-STOCK-2-125ML.....	41	IV-STOCK-67-125ML.....	57	MSCU-10PPM-125ML.....	26
ICMPA1-125ML.....	91	IV-STOCK-2-500ML.....	41	IV-STOCK-68-125ML.....	57	MSFE-100PPM-125ML.....	28
ICMSA1-125ML.....	89	IV-STOCK-3-125ML.....	41	IV-STOCK-69-125ML.....	57	MSFE-100PPM-500ML.....	28
ICNA1-125ML.....	91	IV-STOCK-3-500ML.....	41	IV-STOCK-70-125ML.....	57	MSFE-10PPM-125ML.....	26
ICNA1-500ML.....	91	IV-STOCK-4-125ML.....	41	IV-STOCK-72-125ML.....	80	MSGE-100PPM-125ML.....	28
ICNH41-125ML.....	91	IV-STOCK-4-500ML.....	41	IV-STOCK-73-125ML.....	80	MSGE-10PPM-125ML.....	26
ICNH41-500ML.....	91	IV-STOCK-6-125ML.....	41	IV-STOCK-74-500ML.....	47	MSHF-100PPM-125ML.....	28
ICNNH41-125ML.....	91	IV-STOCK-7-125ML.....	42, 92	IV-STOCK-75-125ML.....	47	MSHF-10PPM-125ML.....	26
ICNNH41-500ML.....	91	IV-STOCK-7-500ML.....	42, 92	IV-STOCK-77-500ML.....	47	MSHG-100PPM-125ML.....	28, 58
ICNNO21-125ML.....	89	IV-STOCK-8-125ML.....	42	IV-STOCK-78-125ML.....	57	MSHG-10PPM-125ML.....	27, 51, 58
ICNNO21-500ML.....	89	IV-STOCK-8-500ML.....	42	IV-STOCK-79-125ML.....	57	MSHG-10PPM-500ML.....	27, 51, 58
ICNNO31-125ML.....	89	IV-STOCK-9-125ML.....	42	IV-STOCK-80-500ML.....	47	MSHG-1PPM-125ML.....	80, 104
ICNNO31-500ML.....	89	IV-STOCK-10-125ML.....	42	LACB1-500ML.....	101	MSHG-1PPM-500ML.....	80, 104
ICNO2-100PPM-125ML.....	90	IV-STOCK-12-125ML.....	42	MM-MG-10-125ML.....	101	MSHGN-100PPM-125ML.....	29, 58
ICNO21-125ML.....	89	IV-STOCK-13-125ML.....	42	MM-P-40-125ML.....	101	MSHGN-10PPM-125ML.....	27, 58
ICNO21-500ML.....	89	IV-STOCK-14-500ML.....	43	MM-PD-10-125ML.....	101	MSHGN-10PPM-500ML.....	27, 58
ICNO31-125ML.....	89	IV-STOCK-1643-125ML.....	52	MM-PD-5-125ML.....	101	MSHO-100PPM-125ML.....	28
ICNO31-500ML.....	89	IV-STOCK-1643-500ML.....	52	MM-PDMG-32-125ML.....	101	MSHO-10PPM-125ML.....	26
ICNTA1-125ML.....	89	IV-STOCK-17-125ML.....	43	MS10B-10PPM-100ML.....	26	MSIN-100PPM-125ML.....	28
ICOAC1-125ML.....	89	IV-STOCK-18-125ML.....	43	MS11B-10PPM-100ML.....	26	MSIN-10PPM-125ML.....	26
ICOAC1-500ML.....	89	IV-STOCK-21-125ML.....	43	MS6LI-100PPM-100ML.....	26	MSK-100PPM-125ML.....	29
ICOPR1-125ML.....	90	IV-STOCK-21-500ML.....	43	MS6LI-100PPM-125ML.....	28	MSK-100PPM-500ML.....	29
ICOXA1-125ML.....	90	IV-STOCK-22-125ML.....	44	MS6LI-10PPM-125ML.....	26	MSK-10PPM-125ML.....	27
ICOXA1-500ML.....	90	IV-STOCK-23-500ML.....	44	MSAELUENT-100ML.....	94	MSLI-100PPM-125ML.....	28
ICPO41-125ML.....	90	IV-STOCK-24-125ML.....	44	MSAELUENT-500ML.....	94	MSLI-100PPM-500ML.....	28
ICPO41-500ML.....	90	IV-STOCK-24-500ML.....	44	MSAG-100PPM-125ML.....	29	MSLI-10PPM-125ML.....	26, 51
ICPPO41-125ML.....	90	IV-STOCK-26-125ML.....	44	MSAG-100PPM-500ML.....	29	MMSG-100PPM-125ML.....	28
ICPPO41-500ML.....	90	IV-STOCK-27-125ML.....	44	MSAG-10PPM-125ML.....	27	MMSG-10PPM-125ML.....	26
ICRB1-125ML.....	91	IV-STOCK-28-125ML.....	45	MSAL-100PPM-125ML.....	28	MMSG-10PPM-500ML.....	26
ICS2031-125ML.....	90	IV-STOCK-29-125ML.....	45	MSAL-10PPM-125ML.....	26	MSMN-100PPM-125ML.....	28
ICS2031-500ML.....	90	IV-STOCK-30-125ML.....	45	MSAS-100PPM-125ML.....	28, 58	MSMN-10PPM-125ML.....	26
ICSCC1-125ML.....	90	IV-STOCK-31-125ML.....	45			MSMN-10PPM-500ML.....	26

MSMO-100PPM-125ML.....	29	MSV-10PPM-125ML.....	27	PH-4-500ML.....	109	QCP-CN-20ML.....	103
MSMO-10PPM-125ML.....	27	MSW-100PPM-125ML.....	29	PH-5-1L.....	109	QCP-CR6-20ML.....	103
MSNA-100PPM-125ML.....	29	MSW-10PPM-125ML.....	27	PH-5-250ML.....	109	QCP-HG-20ML.....	103
MSNA-100PPM-500ML.....	29	MSY-100PPM-125ML.....	29	PH-5-4L.....	109	QCP-MIN-500ML.....	104
MSNA-10PPM-125ML.....	27	MSY-100PPM-500ML.....	29	PH-5-500ML.....	109	QCP-MTL-20ML.....	104
MSNA-10PPM-500ML.....	27	MSY-10PPM-125ML.....	27	PH-6.86-10L.....	109	QCP-NT-20ML.....	104
MSNI-100PPM-125ML.....	29	MSZN-100PPM-125ML.....	29	PH-6.86-1L.....	109	QCP-NUT-1-20ML.....	105
MSNI-10PPM-125ML.....	27	MSZN-100PPM-500ML.....	29	PH-6.86-250ML.....	109	QCP-OG-A-20ML.....	104
MSOS-100PPM-125ML.....	29	MSZN-10PPM-125ML.....	27	PH-6.86-4L.....	109	QCP-OG-W-250ML.....	104
MSOS-10PPM-125ML.....	27	MSZN-10PPM-500ML.....	27	PH-6.86-500ML.....	109	QCP-QCS-1-125ML.....	71, 77
MSP-100PPM-125ML.....	29			PH-6-10L.....	109	QCP-QCS-2-125ML.....	71, 77
MSP-100PPM-500ML.....	29	PE-CHK-1-125ML.....	52	PH-6-1L.....	109	QCP-QCS-3-125ML.....	81
MSP-10PPM-125ML.....	27	PE-TS-1-125ML.....	53	PH-6-250ML.....	109	QCP-QCS-3-500ML.....	81
MSPB-100PPM-125ML.....	28, 58	PE-TS-1-500ML.....	53	PH-6-4L.....	109	QCP-QCS-4-125ML.....	81
MSPB-100PPM-500ML.....	28, 58	PH-1.68-1L.....	109	PH-6-500ML.....	109	QCP-QCS-5-125ML.....	95
MSPB-10PPM-125ML.....	26, 58	PH-1.68-250ML.....	109	PH-7-10L.....	109	QCP-RAIN-125ML.....	105
MSPB-10PPM-125ML.....	26, 58	PH-1.68-4L.....	109	PH-7-1L.....	109	QCP-SLD-450ML.....	105
MSPT-100PPM-125ML.....	29	PH-1.68-500ML.....	109	PH-7-250ML.....	109	QCP-TMS-20ML.....	105
MSPT-10PPM-125ML.....	27	PH-10.01-500ML.....	110	PH-7-4L.....	109	QCP-TRC-10ML.....	105
MSRH-100PPM-125ML.....	29	PH-10-10L.....	110	PH-7-500ML.....	109	QCP-TURB-20ML.....	105
MSRH-10PPM-125ML.....	27	PH-10-1L.....	110	PH-8-10L.....	110	QCP-WH-500ML.....	106
MSRHN-100PPM-125ML.....	29	PH-10-250ML.....	110	PH-8-1L.....	110		
MSRHN-10PPM-125ML.....	27	PH-10-4L.....	110	PH-8-250ML.....	110	TCLP-1REV-125ML.....	82
MSS-100PPM-125ML.....	29	PH-10-500ML.....	110	PH-8-4L.....	110	TCLP-AA-HG-125ML.....	82
MSS-10PPM-125ML.....	27	PH-11-10L.....	110	PH-8-500ML.....	110	THERMO-4AREV-1L.....	53
MSSB-100PPM-125ML.....	28	PH-11-1L.....	110	PH-9.18-1L.....	110	THERMO-4AREV-500ML.....	53
MSSB-10PPM-125ML.....	26	PH-11-250ML.....	110	PH-9.18-250ML.....	110	THERMO-5A-125ML.....	53
MSSC-100PPM-125ML.....	29	PH-11-4L.....	110	PH-9.18-500ML.....	110	THERMO-5A-250ML.....	53
MSSC-100PPM-500ML.....	29	PH-11-500ML.....	110	PH-9-10L.....	110	TOCKHP1-125ML.....	30, 106
MSSC-10PPM-125ML.....	27	PH-12.45-500ML.....	110	PH-9-1L.....	110	TOCKHP1-500ML.....	30, 106
MSSC-10PPM-125ML.....	27	PH-12.47-1L.....	110	PH-9-250ML.....	110	TUNE F-X-SERIES-125ML.....	54
MSSC-10PPM-125ML.....	27	PH-12.47-250ML.....	110	PH-9-4L.....	110		
MSSC-10PPM-500ML.....	29	PH-12.47-500ML.....	110	PH-9-500ML.....	110	VAR-CAL-1-125ML.....	54
MSSC-10PPM-500ML.....	29	PH-12-10L.....	110	PHBLUE-10-10L.....	110	VAR-CAL-1-500ML.....	54
MSSC-10PPM-500ML.....	29	PH-12-1L.....	110	PHBLUE-10-1L.....	110	VAR-CAL-2-125ML.....	54
MSSC-10PPM-500ML.....	29	PH-12-250ML.....	110	PHBLUE-10-250ML.....	110	VAR-CAL-7-125ML.....	54
MSSC-10PPM-500ML.....	29	PH-12-4L.....	110	PHBLUE-10-4L.....	110	VAR-CAL-7-500ML.....	54
MSSC-10PPM-500ML.....	29	PH-12-500ML.....	110	PHBLUE-10-500ML.....	110	VAR-IS-1-125ML.....	55
MSSC-10PPM-500ML.....	29	PH-2-10L.....	109	PHRED-4-10L.....	109	VAR-TS-MS-125ML.....	55
MSSC-10PPM-500ML.....	29	PH-2-1L.....	109	PHRED-4-1L.....	109	WW-CAL-1A-125ML.....	69, 73
MSSC-10PPM-500ML.....	29	PH-2-250ML.....	109	PHRED-4-250ML.....	109	WW-CAL-2-125ML.....	69, 73
MSSC-10PPM-500ML.....	29	PH-2-4L.....	109	PHRED-4-4L.....	109	WW-CAL-3-125ML.....	69, 73
MSSC-10PPM-500ML.....	29	PH-2-500ML.....	109	PHRED-4-500ML.....	109	WW-CAL-4A-125ML.....	70, 73
MSSC-10PPM-500ML.....	29	PH-3-10L.....	109	PHYELLOW-7-10L.....	110	WW-CAL-4B-125ML.....	70, 73
MSSC-10PPM-500ML.....	29	PH-3-1L.....	109	PHYELLOW-7-1L.....	110	WW-CAL-5-125ML.....	70, 74
MSSC-10PPM-500ML.....	29	PH-3-250ML.....	109	PHYELLOW-7-250ML.....	110	WW-IPC-1-125ML.....	74
MSSC-10PPM-500ML.....	29	PH-3-4L.....	109	PHYELLOW-7-4L.....	110	WW-IPC-2-125ML.....	75
MSSC-10PPM-500ML.....	29	PH-3-500ML.....	109	PHYELLOW-7-500ML.....	110	WW-IPC-3-125ML.....	75
MSSC-10PPM-500ML.....	29	PH-4.01-500ML.....	109	PT-CO-500ML.....	114	WW-LFS-1-125ML.....	76
MSSC-10PPM-500ML.....	29	PH-4-10L.....	109			WW-LFS-2-125ML.....	76
MSSC-10PPM-500ML.....	29	PH-4-1L.....	109	QCP-CAT-20ML.....	103	WW-MSCAL-1-125ML.....	79
MSSC-10PPM-500ML.....	29	PH-4-250ML.....	109	QCP-CICV-1-125ML.....	62, 64, 66	WW-MSCAL-2-125ML.....	79
MSSC-10PPM-500ML.....	29	PH-4-4L.....	109	QCP-CICV-3-125ML.....	62, 64, 66		

Pricing

Visit our website for all current pricing information:

inorganicventures.com

Shipping

Orders shipped within the U.S. are sent FOB shipping point via FedEx. Other carriers are also available. If you prefer a specific carrier, let us know.

Terms

Our terms are net 30 days. To establish your account, contact us during regular business hours. 8:00 a.m. to 5:00 p.m. EST.

Credit Orders

We gladly accept MasterCard, Visa or American Express.



Cautionary Notice

OUR PRODUCTS ARE NOT FOR DRUG, FOOD OR HOUSEHOLD APPLICATIONS. They are intended for laboratory use only by qualified individuals trained in the proper handling of such materials. The consumer assumes all responsibility for the safe storage, handling, disposal and application of any Inorganic Ventures product.

Copyright

All content and graphics appearing herein are © 1985–2025 by Inorganic Ventures, Inc. and may not be reproduced for any purpose without Inorganic Ventures' express permission.

* Refunds, replacements and exchanges are considered at management's discretion.

Customer Service

Representatives are available Monday through Friday, between 8:00 a.m. and 5:00 p.m. EST.

Phone: 800.669.6799 (US & Canada)
+1.540.585.3030 (International)
Fax: 540.585.3012
Email: info@inorganicventures.com
Online: inorganicventures.com

Export Orders

We have representation in countries across the globe. Contact your distributor for regional pricing, orders and additional catalogs. Visit our website to find the distributor in your area.

Our Guarantee

As stated in our Declaration of Integrity (see pg. 7), if you're dissatisfied with your order for any reason, we'll resolve the situation in whatever way works best for you:

- A full refund;
- Complimentary technical services; or
- A replacement item rushed to you at no cost.*

Returns, Exchanges and Damaged Goods

Situations not covered under the terms of our guarantee (above) may require a return and/or exchange of the product. Refunds, replacements and exchanges are considered at management's discretion. Items will only be accepted for return with a valid return authorization number. Contact us to receive this number.

All transportation and packaging expenses on returned or exchanged items are assumed by the customer. Credit for returned items is subject to a modest 15% restocking charge. Credit will not be issued until it is determined by our inspection that returned items are unused and undamaged. All returned items must be received within 60 days of the original invoice date.

If carton is received in a conspicuously damaged condition, alert carrier and refuse acceptance.



INORGANIC CUSTOM &
STOCK CERTIFIED
REFERENCE MATERIALS



300 Technology Drive | Christiansburg, VA 24073, USA
Phone: 800.669.6799 or 540.585.3030 | Fax: 540.585.3012
info@inorganicventures.com | inorganicventures.com

 Please recycle this catalog.